



Red DCI[®] and DCIP[™]

Adult and Pediatric Reusable SpO₂ finger sensors

Images	2
en English	3-5
fr French	6-8
de German	9-11
it Italian	12-14
es Spanish	15-17
sv Swedish	18-20
nl Dutch	21-23
da Danish	24-26
pt Portuguese	27-29
zh Chinese	30-32
ja Japanese	33-35
fi Finnish	36-38
no Norwegian	39-41
cs Czech	42-44
hu Hungarian	45-47
pl Polish	48-50
ro Romanian	51-53
sk Slovak	54-56
tr Turkish	57-59
el Greek	60-62
ru Russian	63-65
ar Arabic	68-66
fa Farsi	71-69

Red DCI[®] and DCIP[™]

Adult and Pediatric Reusable SpO₂ Sensors

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

Body weight

Do not discard

INDICATIONS

Red reusable DCI[®] and DCIP[™] sensors are indicated for either spot “check” or the continuous non-invasive monitoring of arterial oxygen saturation (SpO₂) and pulse rate.

Red DCI and DCIP sensors are for use only with instruments containing Masimo Rainbow[®] SET[™] Technology or licensed to use Rainbow compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular instruments and sensor models. Each instrument manufacturer is responsible for determining whether their instruments are compatible with each sensor model.

NOTE: Do not attempt to use Red reusable sensors on technology that does not have Masimo Rainbow SET Technology.

The sensors are indicated as follows:

	Red DCI	Red DCIP
	Masimo SET Technology	Masimo SET Technology
	> 30 kg	10 - 50 kg
Application site	Finger or Toe	Finger
Saturation Accuracy, No Motion	± 2%	± 2%
Saturation Accuracy, Motion	± 3%	± 3%
Pulse Rate Accuracy, No Motion	± 3 bpm	± 3 bpm
Pulse Rate Accuracy, Motion	± 5 bpm	± 5 bpm
Low Perfusion Accuracy	SpO ₂ ± 2%	SpO ₂ ± 2%
	Pulse ± 3 bpm	Pulse ± 3 bpm

The Red DCI and DCIP have been validated on the Radical-7[™] Pulse CO-Oximeter containing Masimo Rainbow SET Technology.

CONTRAINDICATIONS

Red DCI and DCIP sensors are contraindicated for use on mobile patients or for prolonged periods of use. They must be removed and repositioned to a different monitoring site at least every four (4) hours.

INSTRUCTIONS

A) Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring or middle finger of the non-dominant hand is preferred.
- Alternatively, the other digits on the non-dominant hand may be used. Always choose a site that will completely cover the sensor's detector window. The great toe or long toe (next to the great toe) may be used on restrained patients or patients whose hands are unavailable.
- Site should be cleaned of debris and dry prior to sensor placement.

B) Attaching the sensor to the patient

- Open the sensor by pressing on hinge tabs. Place the selected digit over the sensor window of the sensor. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. The tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop (Fig 1).
- The hinged tabs of the sensor should open to evenly distribute the grip of the sensor along the length of the finger. Check position of sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data (Fig 2).
- Orient the sensor so that the cable will be running towards the top of the patient's hand (Fig 3). Connect the sensor connector to a patient cable or directly to the device.

C) Disconnecting the Sensor from the Patient Cable

- Place thumb and index finger on gray buttons on either side of the patient cable connector.
- Press firmly on the gray buttons and pull to remove the sensor (Fig. 4).

CAUTION

To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not sterilize by irradiation, steam, autoclave or ethylene oxide.

CLEANING: To clean the sensor, first remove it from the patient and disconnect it from the patient cable. You may then clean the sensor by wiping it with a 70% isopropyl alcohol pad. Allow the sensor to dry prior to placement on a patient.

Red DCI[®] and DCIP[™]

Adult and Pediatric Reusable SpO₂ Sensors

en

WARNINGS

- The site must be checked at least every four (4) hours to ensure proper adhesion, skin integrity and proper alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every two (2) hours with poorly perfused patients.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the readings may read lower than core arterial oxygen saturation.
- Misapplied sensors or sensors that become partially dislodged may cause either over or under reading of actual arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage or damage the sensor.
- The sensor should be free of visible defects. Never use a damaged sensor or one with exposed electrical circuitry.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Intravascular dyes or externally applied coloring (such as nail polish) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MethHb) will lead to inaccurate SpO₂ measurements.
- Failure to apply the sensor properly may cause incorrect measurements.
- Do not use the sensor during MRI scanning.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should be not below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor)
- Venous pulsations may cause erroneous low readings (e.g. tricuspid valve regurgitation).
- Do not modify or alter the sensor in any way. Alterations or modification may affect performance and/or accuracy.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal[™] technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

SPECIFICATIONS

When used with Pulse CO-Oximetry monitors containing Masimo Rainbow SET Technology Pulse CO-Oximetry monitors or with licensed Masimo Rainbow SET pulse co-oximetry modules, during no motion, the accuracy of the Red reusable DCI and DCIP sensors is from 70% to 100% SpO₂ is ± 2 digits (± 1 Std. Dev.) for adults. Pulse rate accuracy from 25-240 bpm is ± 3 digits, (± 1 Std. Dev.).

The Saturation accuracy of Red reusable DCI and DCIP sensors during motion from 70% to 100% SpO₂ is ± 3 digits (± 1 Std. Dev.) for adults. The pulse rate accuracy of Red reusable DCI and DCIP sensors during motion from 25-240 bpm is ± 5 digits (± 1 Std. Dev.). Red reusable DCI and DCIP sensors have been validated on the Radical-7 Pulse CO-Oximeter.

INSTRUMENT CAPABILITY

 These sensors are intended for use only with instruments containing Masimo Rainbow SET Technology or pulse oximetry monitors licensed to use Rainbow compatible sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original instrument manufacturer. Use of this sensor with other instruments may result in no or improper performance.

WARRANTY

Masimo warrants to the initial buyer only that each product it manufactures, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

Red DCI[®] and DCIP[™]

Adult and Pediatric Reusable SpO₂ Sensors


en

NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RED DCI AND DCIP SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Caution		Storage temperature range
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		Storage humidity limitation
	Mark of conformity to European medical device directive 93/42/EEC		Do not use if package is damaged
	Manufacturer		Separate collection for electrical and electronic equipment (WEEE).
	Date of manufacture YYYY-MM-DD		Keep dry
	Authorized representative in the European community		Not made with natural rubber latex
	Catalogue number (model number)		Atmospheric pressure limitation
	Masimo reference number		Greater than
	Lot code		Less than
	Body weight		Fragile, handle with care
	Non-Sterile		Do not discard
	Follow Instructions for Use		Pediatric patient

Covered by one or more of the following U.S. Patents: 6,580,086, 6,397,091, 6,011,986, 5,758,644, international equivalents, or one or more of the patents referenced at www.masimo.com/patents. Other patents pending.

Masimo SET technology and Rainbow technology under license from Masimo Corporation.

Masimo, SET, , Rainbow and DCI are federally registered trademarks of Masimo Corporation.

Radical-7 and DCIP are trademarks of Masimo Corporation.

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© 2016 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH