
RD SedLine® Patient Cable

Images	2
en English	3-5
fr French	6-8
de German	9-11
it Italian	12-14
es Spanish	15-17
sv Swedish	18-20
nl Dutch	21-23
da Danish	24-26
pt Portuguese	27-29
zh Chinese	30-32
ja Japanese	33-35
fi Finnish	36-38
no Norwegian	39-41
cs Czech	42-44
hu Hungarian	45-47
pl Polish	48-50
ro Romanian	51-53
sk Slovak	54-56
tr Turkish	57-59
el Greek	60-62
ru Russian	63-65
ar Arabic	68-66
fa Farsi	71-69

RD SedLine® Patient Cable

RD SedLine® Patient Cable

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

DESCRIPTION

The RD SedLine Patient Cable transfers analog EEG signals collected from the Sensor to the Module for processing. The Patient Cable is a reusable cable and may be used from patient to patient.

Always refer to the operator's manual for complete instructions or additional instructions.

The RD SedLine Patient Cable is compatible with the RD SedLine EEG Sensor.

INSTRUCTIONS

A) Connecting

To connect the Patient Cable to the Module

1. Align the ridged Patient Cable connector end with the available Module connector end.
2. Push to insert.

To connect the Patient Cable to the Sensor

1. Align the yellow indicator connector on the end of the RD SedLine Patient Cable with the yellow indicator on the Sensor Connector.
2. Refer to Fig. 1. Push together until the connections fit securely.

B) Disconnecting the Sensor

To disconnect the Patient Cable from the Module

- Gently pull the Patient Cable out of the Module by grasping the patient cable connector end.

To disconnect the Patient Cable from the Sensor

- Gently pull the Sensor out of the Patient Cable connector from the Sensor's connector end.

CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad, quaternary ammonium chloride solution, or a 10% (1:10) chlorine bleach to water solution.
3. Wipe all surfaces of the cable.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.
5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

The surface of the cable has been tested to be chemically resistant to the following solutions: Quaternary Ammonium Chloride solution wipe, Accelerated Hydrogen Peroxide, Chlorhexidine, Cidex, 70% and 90% Isopropyl Alcohol, and 10% (1:10) chlorine bleach to water solution.

CAUTION

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

SPECIFICATIONS

Cable length	3 feet (0.9 meters)
Patient module connector	9 pin proprietary
Sensor connector	Custom connector

ENVIRONMENTAL

Storage/Transport Temperature	-40° C to +70° C, ambient humidity
Operating Temperature	+5° C to +40° C, ambient humidity
Storage/Transport Humidity	10% to 95%, non-condensing
Operating Humidity	10% to 95%, non-condensing

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

This single-patient use sensor is licensed to you under the patents owned by masimo for single-patient use only. By acceptance or use of this product, you acknowledge and agree that no license is granted for use of this product with more than a single patient.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Non-Sterile
	Authorized representative in the European community		Not made with natural rubber latex
	Storage temperature range		Separate collection for electrical and electronic equipment (WEEE).
	Storage humidity limitation		Keep dry
	Atmospheric pressure limitation		Use by date YYYY-MM-DD
	Consult instructions for use		Body weight
	Date of manufacture YYYY-MM-DD		Do not re-use/Single patient use only
	Manufacturer		Do not use if package is damaged
	Fragile, handle with care		Lot code
	Mark of conformity to European medical device directive 93/42/EEC		Catalogue number (model number)
	Masimo reference number		Serial number
Rx ONLY	Caution: Federal law restricts this instrument to sale by or on the order of a physician (USA audiences only)		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.masimo.com/TechDocs Note: eFU is not available for CE mark countries.

Patents: <http://www.masimo.com/patents.htm>

Masimo and SedLine are federally registered trademarks of the Masimo Corporation.

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© 2016 Masimo Corporation
MDSS GmbH