M-LNCS™, LNCS® Series, DCI® & DCI-P








Adult & Pediatric Reusable Finger Sensor

Images	2
en English	3-7
fr French	8-12
de German	13-17
it Italian	18-22
es Spanish	23-27
sv Swedish	28-32
nl Dutch	33-37
da Danish	38-42
pt Portuguese	43-47
zh Chinese	48-52
ja Japanese	53-57
fi Finnish	58-62
no Norwegian	63-67
cs Czech	68-72
hu Hungarian	73-77
pl Polish	78-82
ro Romanian	83-87
sk Slovak	88-92
tr Turkish	93-97
el Greek	98-102
ru Russian	103-107
ar Arabic	112-108
fa Farsi	117-113
Performance Specification	118

M-LNCS™, LNCS® Series, DCI® & DCI-P Adult & Pediatric Reusable Finger Sensor

M-LNCS™, LNCS® Series, DCI® & DCI-P

en

Adult & Pediatric Reusable Finger Sensor

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

INDICATIONS

When used with Masimo SET®:

The M-LNCS™, LNCS® DCI® and DCI-P reusable sensors are indicated for either “spot check” or continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

When used with Nellcor® and Nellcor Compatible Pulse Oximeters:

The M-LNCS, LNCS DCI and DCI-P reusable sensors are indicated for either “spot check” or continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

The DCI and DCI-P sensors are contraindicated for use on active patients or for prolonged periods of use. They are not intended for long-term monitoring. They must be removed and repositioned to a different monitoring site at least every four (4) hours. If extended monitoring is required, use of a M-LNCS or LNCS disposable sensor is recommended.

DESCRIPTION

The DCI and DCI-P are for use with devices containing Masimo SET oximetry or licensed to use M-LNCS, LNCS sensors and also with Nellcor and Nellcor compatible pulse oximeters, except Nellcor OxiMax® enabled instruments. Consult individual instrument manufacturer for compatibility of particular instrument and sensor models. Each instrument manufacturer is responsible for determining whether its devices are compatible with each sensor model.

The DCI and DCI-P sensors have been verified with Masimo SET Oximetry Technology and on Nellcor's N-200 Pulse Oximeter.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be checked frequently or per clinical protocol to ensure adequate circulation, skin integrity and correct optical alignment.
- Exercise extreme caution; skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess site as frequently as every (1) hour and move the sensor if there are signs of loss of skin integrity and/or loss of circulation or perfusion.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Carefully route sensor and patient cable to reduce the possibility of patient entanglement or strangulation.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect readings.
- Misapplications due to wrong sensor types can cause inaccurate or no readings.
- Inaccurate SpO₂ readings may be caused by abnormal venous pulsation or venous congestion.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can affect the pulse rate displayed on the oximeter. Verify patient's pulse rate against the ECG heart rate.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.

- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment as it may result in physical harm.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- Inaccurate readings may be caused by EMI radiation interference.
- Abnormal fingers, Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate or no readings.
- High levels of COHb or Methb may occur with a seemingly normal SpO₂. When elevated levels of COHb or Methb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ readings.
- Elevated levels of Methemoglobin (Methb) may lead to inaccurate SpO₂ readings.
- Elevated levels of total bilirubin may lead to inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by severe anemia, very low arterial perfusion, or extreme motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by vasospastic disease such as Raynaud's, and peripheral vascular disease.
- Inaccurate SpO₂ readings may be caused by elevated levels of dyshemoglobin, hypocapnic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- SpO₂ readings may be affected under very low perfusion conditions at the monitored site.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the sensor in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide as it will damage the sensor.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal™ technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A) Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring finger of the non-dominant hand is preferred. Alternatively, the other digits on the non-dominant hand may be used.
- Always choose a site that will completely cover the sensor's detector window. The great toe or long toe (next to the great toe) may be used on restrained patients or patients whose hands are unavailable.
- Site should be cleaned of debris prior to sensor placement.

CAUTION: Before using the sensor, ensure that the sensor is physically intact, with no broken or frayed wires or damaged parts.

B) Attaching the sensor to the patient

1. Refer to **Fig. 1**. Open the sensor by pressing on hinge tabs. Place the selected digit over the sensor window of the DCI or DCI-P. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. On a finger site, the tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop.
2. Refer to **Fig. 2**. The hinged tabs of the sensor should open to evenly distribute the grip of the sensor along the length of the finger. Check position of sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to **Fig. 3**. Orient the sensor so that the cable will be running towards the top of the patient's hand.

NOTE: With smaller digits, in order to completely cover the detector window, it may be necessary to use a lower weight range sensor. The sensor is not intended for use on the thumb or across a child's hand or foot.

C) Attaching the Sensor to the Patient Cable

1. Refer to **Fig. 4 M-LNCS, Fig. 5 LNCS**. Properly orient the sensor connector (1) and insert the sensor connector completely into the patient cable connector (2).
2. Close the protective cover (3).

D) Disconnecting the Sensor from the Patient Cable

1. Refer to **Fig. 6 M-LNCS, Fig. 7 LNCS**. Lift the protective cover (1) to gain access to the sensor connector (2).
2. Pull firmly on the sensor connector (2) to remove it from the patient cable (3).

NOTE: To avoid damage, pull on the sensor connector, not the cable.

CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the sensor by wiping with a 70% isopropyl alcohol pad or mild detergent.
3. Allow the sensor to dry thoroughly prior to placement on a patient.

or

1. If low-level disinfection is required, use a 1:10 bleach / water solution.
2. Saturate a cloth or gauze pad with the cleaning solution and wipe all surfaces of the sensor and cable.
3. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the sensor and cable.
4. Dry the sensor and cable with a clean cloth or dry gauze pad.

CAUTIONS:

- Do not use undiluted bleach (5% - 5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor may occur.
- Do not immerse the sensor or connector in any liquid solution.
- Do not sterilize by irradiation, steam autoclave or ethylene oxide.

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, with licensed Masimo SET pulse oximetry modules, or with Nellcor technology, M-LNCS/LNCS DCI and DCI-P reusable sensors are intended for the following performance specifications:

	Masimo SET Technology		Nellcor Technology ⁵	
	DCI	DCI-P	DCI	DCI-P
 Body Weight	> 30 kg	10 - 50 kg	> 30 kg	10 - 50 kg
Application site	Finger or toe	Finger or toe	Finger or toe	Finger or toe
Saturation Accuracy, No Motion ¹ (70-100%)	2%	2%	2%	2%
Saturation Accuracy, Motion ² (70-100%)	3%	3%	N/A	N/A
Pulse Rate Accuracy, No Motion ³ (70-100%)	3 bpm	3 bpm	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ³ (25 - 240 bpm)	5 bpm	5 bpm	N/A	N/A
Low Perfusion Accuracy ⁴ (70-100%)	SpO ₂ ± 2%	SpO ₂ ± 2%	N/A	N/A
	Pulse ± 3 bpm	Pulse ± 3 bpm	N/A	N/A

NOTE: The table indicates ARMS Accuracy which is calculated based upon measurement values that are statistically distributed; approximately 68% of the measured values fell within +/- the ARMS value when compared to the reference device in a controlled study.

¹ The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

² The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

³ The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁴ The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁵ When used with Nellcor and Nellcor compatible pulse oximeters during no motion, the accuracy of the M-LNCS, LNCS DCI and DCI-P sensors from 70-100% SpO₂ is ± 2 digits (1 Std. Dev.) for adults/pediatrics/infants. Pulse rate accuracy from 25-240 bpm is ± 3 bpm (1 Std. Dev.). The M-LNCS, LNCS DCI and DCI-P sensors have been validated on Nellcor's N-200 pulse oximeter.

COMPATIBILITY

This sensor is intended for use only with devices containing Masimo SET Technology or pulse oximetry monitors licensed to use SET compatible sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE M-LNCS, LNCS SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
 For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult Instructions for Use		Do not discard
	Follow Instructions for Use		Lot code
	Manufacturer		Catalogue number (model number)
	Date of Manufacture YYYY-MM-DD		Masimo reference number
	Non-sterile		Body weight
	Not made with natural rubber latex		Greater than
Rx ONLY	Federal law (U.S.A) restricts this device to sale by or on the order of a physician		Less than
	Mark of Conformity to European Medical Device Directive 93/42/EEC		Storage Humidity Limitation
	Authorized representative in the European community		Keep Dry
	Separate collection for electrical and electronic equipment (WEEE).		Storage Temperature Range
	Do not use if package is damaged		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.

Patents: <http://www.masimo.com/patents.htm>

M-LNCS and X-Cal are trademarks of Masimo Corporation

Masimo, SET, LNCS, DCI and  are federally registered trademarks of Masimo Corporation

Nellcor and Oximax are registered trademarks of Nellcor Puritan Bennett Incorporated

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© 2016 Masimo Corporation
MDSS GmbH