






# RD SedLine® EEG Sensor


Images	2
<b>en</b> English	3-5
<b>fr</b> French	6-8
<b>de</b> German	9-11
<b>it</b> Italian	12-14
<b>es</b> Spanish	15-17
<b>sv</b> Swedish	18-20
<b>nl</b> Dutch	21-23
<b>da</b> Danish	24-26
<b>pt</b> Portuguese	27-29
<b>zh</b> Chinese	30-32
<b>ja</b> Japanese	33-35
<b>ko</b> Korean	36-38
<b>fi</b> Finnish	39-41
<b>no</b> Norwegian	42-44
<b>cs</b> Czech	45-47
<b>hu</b> Hungarian	48-50
<b>pl</b> Polish	51-53
<b>ro</b> Romanian	54-56
<b>sk</b> Slovak	57-59
<b>tr</b> Turkish	60-62
<b>el</b> Greek	63-65
<b>ru</b> Russian	66-68
<b>ar</b> Arabic	71-69
<b>fa</b> Farsi	74-72

# RD SedLine® EEG Sensor

## DIRECTIONS FOR USE

Do not re-use

Not made with natural rubber latex

Non-sterile

Leave on for no more than 24 hours

Prior to using this Sensor the user should read and understand the Operator's Manual for the Device or Monitor and this Directions for Use.

### INDICATIONS

The RD SedLine EEG Sensor electrodes are applied directly to the patient's skin to enable recording of electrophysiologic signals (such as EEG).

### CONTRAINDICATIONS

The RD SedLine EEG Sensor is not intended for use in children less than 18 years of age.

### DESCRIPTION

The RD SedLine EEG Sensor (1) is comprised of 6 gelled electrodes, 4 active channels (R1, R2, L1, L2), one reference channel (CT), and one ground channel (CB). The Sensor is a single-use, non-sterile product that does not contain natural rubber latex.

Refer to the SedLine Operator's Manual for additional information including Warnings and Cautions.

**WARNING:** Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

### WARNINGS, CAUTIONS, AND NOTES

- **Warning:** If skin irritation is noticed, discontinue use of the Sensor.
- **Warning:** To avoid cross contamination, only use Masimo single-use sensor on the same patient.
- **Warning:** Do not place electrodes between surgical site and the electro-surgical return electrode. Doing so increases the risk of burns in case of a defect in the electro-surgical return electrode.
- **Caution:** Federal law restricts this device to sale by or on the order of a physician.
- **Caution:** For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.
- **Caution:** Use the Sensor immediately after opening to prevent gel from drying out.
- **Caution:** Do not use the Sensor if pouch or product is defective.
- **Caution:** Do not allow electrode gel to come into contact with the SedLine Sensor connector as this may cause impedance problems and inaccurate PSI readings.
- **Caution:** Each Sensor is approved for use up to 24 hours. If monitoring in excess of this period is desired, replace the Sensor every 24 hours.
- **Caution:** SedLine is designed to operate with the sensor. Simultaneous use with other electrodes on the forehead is not recommended.
- **Caution:** Do not use scissors to open pouch as this may cause damage to the Sensor.
- **Caution:** Do not cut the Sensor as this may result in improper operation.
- **Caution:** Applying excessive pressure to electrodes may lead to decreased signal quality, decreased PSI reliability, and poor adhesion.
- **Caution:** The conductive parts of electrodes, sensors, or connectors should not contact other conductive parts including earth ground.
- Do not attempt to reprocess, recondition or recycle any Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal™ technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 24 hours of patient monitoring. After single-patient use, discard sensor.

---

## INSTRUCTIONS

### A) Connecting the Sensor

1. Prepare the forehead skin with alcohol and dry before applying the sensor.
2. Detach the sensor connector from the back of the tray (1).
3. Snap the outer edges of the tray from the center tray (2).
4. Peel the center of the sensor from the tray and discard.
5. Center CT/CB electrodes on the forehead above the nose. Press the black lining around the electrodes firmly on the forehead to secure the sensor. Do not press directly on the electrode otherwise the gel may leak out (3).
6. Apply R1/ L1 to either side of CB.
7. Peel back the R2/ L2 from the outer tray edges and apply R2 /L2 in hairless area just above the temple (4).
8. Align the yellow indicator on the end of the patient cable to the yellow indicator on the sensor connection (5).
9. Push together until the sensor connector fits into the patient cable.
10. Ensure the SedLine Module is connected to the SedLine patient cable and the Root monitor.
11. If any of the electrodes are displayed as red, yellow, or blue push and wiggle the white lining around the electrode. Do not press directly on the electrode, otherwise gel may leak out.
12. For additional details, see the SedLine Operators Manual.

### B) Disconnecting the Sensor

To disconnect the Sensor, gently pull the Sensor connector out of the Patient Cable connector.

### C) Sensor Removal

To remove the sensor from the patient's forehead, pull each electrode tab gently away from the skin.

If the sensor is difficult to remove, use alcohol to assist in removal.

**Note:** The Sensor is not reusable. Dispose of the sensor according to local laws.

### D) Applying with the O3 Sensor series when using O3 Regional Oximetry

Following the application of the RD SedLine EEG Sensor, align the rounded section of each O3 sensor with the curves of the O3 sensor marked on the t-section of the RD SedLine EEG Sensor (6a) and position the above L1/R1 and L2/R2 (6b).

---

## WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

---

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

---

## NO IMPLIED LICENSE


This single-patient use sensor is licensed to you under the patents owned by Masimo for single-patient use only. By acceptance or use of this product, you acknowledge and agree that no license is granted for use of this product with more than a single patient.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

For professional use. See instructions for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

**The following symbols may appear on the product or product labeling:**


SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Non-Sterile
	Authorized representative in the European community		Not made with natural rubber latex
	Storage temperature range		Separate collection for electrical and electronic equipment (WEEE).
	Storage humidity limitation		Leave on for no more than 24 hours
	Atmospheric pressure limitation		Keep dry
	Consult instructions for use		Use by date YYYY-MM-DD
	Date of manufacture YYYY-MM-DD		Body weight
	Manufacturer		Do not re-use/Single patient use only
	Fragile, handle with care		Do not use if package is damaged
	Mark of conformity to European medical device directive 93/42/EEC		Lot code
	Masimo reference number		Catalogue number (model number)
<b>Rx ONLY</b>	<b>Caution:</b> Federal law restricts this instrument to sale by or on the order of a physician (USA audiences only)		Serial number
	Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.masimo.com/TechDocs">http://www.masimo.com/TechDocs</a> <b>Note: eIFU is not available for CE mark countries.</b>		

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET and SedLine are federally registered trademarks of the Masimo Corporation.

X-Cal is a trademark of Masimo Corporation.

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© 2017 Masimo Corporation
MDSS GmbH