









# Radius PPG™

## Chip and Receiver

Images	2
<b>en</b> English	3-6
<b>fr</b> French	7-10
<b>de</b> German	11-14
<b>it</b> Italian	15-18
<b>es</b> Spanish	19-22
<b>sv</b> Swedish	23-26
<b>nl</b> Dutch	27-30
<b>da</b> Danish	31-34
<b>pt</b> Portuguese	35-38
<b>zh</b> Chinese	39-42
<b>ja</b> Japanese	43-46
<b>fi</b> Finnish	47-50
<b>no</b> Norwegian	51-54
<b>cs</b> Czech	55-58
<b>hu</b> Hungarian	59-62
<b>pl</b> Polish	63-66
<b>ro</b> Romanian	67-70
<b>sk</b> Slovak	71-74
<b>tr</b> Turkish	75-78
<b>el</b> Greek	79-82
<b>ru</b> Russian	83-86
<b>ko</b> Korean	87-90
<b>ar</b> Arabic	94-91

## Chip and Receiver

## DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this product the user should read and understand the Operator's Manual for the device/monitor, The Directions for Use for the sensor, and this Directions for Use.

**INDICATIONS**

Masimo Radius PPG™ is intended for the non-invasive continuous monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>), pulse rate (PR). Masimo Radius PPG is indicated for the continuous monitoring of functional arterial oxygen saturation of hemoglobin (SpO<sub>2</sub>) and pulse rate (PR) for use with adult, pediatric and neonatal patients during both no motion and motion conditions and for patients who are well or poorly perfused in hospital, hospital-type facilities and home environments.

Devices with Masimo technology are only to be used with Masimo sensors and cables.

**CONTRAINDICATIONS**

The Radius PPG is contraindicated for patients who exhibit allergic reactions to foam rubber products and/or adhesive tape.

**DESCRIPTION**

Radius PPG consists of three parts:

- Radius PPG wireless receiver
- Radius PPG reusable chip
- Radius PPG adhesive sensor (*shipped separately, see Directions for Use for the Radius PPG Adhesive Sensor*)

Radius PPG is a wireless sensor for use with devices containing Masimo technology MX Version 7.14.8.x. or higher. Consult individual device manufacturers for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether their device is compatible with each sensor model.

**WARNINGS, CAUTIONS AND NOTES**

- Do not leave the sensor components unattended around children. Small items may become choking hazards.
- To ensure continued monitoring, routinely verify the wireless connection.
- When using multiple Radius PPG sensors, repeat pairing before monitoring to ensure proper wireless connection.
- When using Radius PPG, keep it within the recommended range from the connected host (see Wireless Technology Information for details); moving outside of this range may cause a loss in connection with the host device.
- When using Radius PPG, relocate the devices away from sources that may interfere with the Bluetooth connection. The presence of other devices that may create radio frequency interference (RFI) may result in loss of Quality of Service (see Specifications for details) of the Bluetooth connection. Devices that may cause RFI include but are not limited to the following: electrocautery equipment, diathermy equipment, other cellular telephones, wireless PC and tablets, pagers, RFID devices, MRI, and electromagnetic security systems.

**INSTRUCTIONS****a) Initial setup**

1. Verify you have all components. Refer to **Fig. 1**.

**A.** Radius PPG adhesive sensor    **B.** Radius PPG wireless bluetooth receiver    **C.** Radius PPG reusable chip    **D.** Radius PPG reusable chip holder for device

2. Turn on the patient monitor.
3. Plug the cable into the patient monitor. Refer to **Fig. 2**. The light on wireless receiver will be white.
4. Attach the wireless receiver to the side of the patient monitor using the adhesive provided. Refer to **Fig. 3**.
5. Avoid covering speakers or holes used for mounting when attaching the module. Refer to **Fig. 4**.
6. Attach the chip holder near the receiver on the patient monitor. Refer to **Fig. 5**.
7. Verify that the wireless receiver is attached and plugged in to patient monitor. Verify that the chip holder is attached to patient monitor. Refer to **Fig. 6**.

**b) Pairing the reusable transmitter chip with the wireless receiver**

1. Ensure the device is powered on. Refer to **Fig. 7**.
2. Hold the reusable chip to the indent on the wireless receiver until the Bluetooth symbol on the wireless receiver turns green. Refer to **Fig. 8**.
3. Insert the reusable chip into the sensor attachment strap until there is a tactile or audible click of connection. Refer to **Fig. 9**.
4. Verify the light on the wireless receiver turns blue. (See **LIGHT INDICATOR GUIDE** section.)

**c) Disconnecting**

1. Push down on the tab to release the reusable chip from the sensor. Refer to **Fig. 10**.
2. After cleaning, store the reusable chip in the chip holder attached to the patient monitor. Refer back to **Fig. 6**.
3. Discard the adhesive sensor and strap.

**CLEANING**

**WARNING:** Before cleaning, make sure the device is off and is not applied to a patient.

Thoroughly clean the Radius PPG reusable chip and wireless receiver before applying to a new patient.

**To surface clean the reusable chip and wireless receiver:**

1. Remove the sensor from the patient and disconnect the reusable chip.
2. Wipe all surfaces of the reusable chip and wireless receiver with one of the following:
  - a. 70% Isopropyl alcohol
  - b. 10% (1:10) chlorine bleach to water solution
  - c. Quaternary ammonium chloride solution
3. Allow the reusable chip and wireless receiver to dry thoroughly before using again.

**CAUTIONS:**

- To avoid permanent damage to the reusable chip and wireless receiver, do not use undiluted bleach (5–5.25% sodium hypochlorite) or any other cleaning solution not recommended.
- Do not immerse the reusable chip and wireless receiver in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave or ethylene oxide.

## LIGHT INDICATOR GUIDE

Color	Wireless receiver	Transmitter chip	Description	Next steps
No light	--	--	<ul style="list-style-type: none"> <li>Wireless receiver cable is not connected to host device with power</li> <li>Chip not connected to sensor with battery</li> </ul>	<ul style="list-style-type: none"> <li>Turn on patient monitor and plug cable into patient monitor</li> <li>See <i>Instructions</i>, section a) for set up</li> </ul>
White	solid	--	<ul style="list-style-type: none"> <li>Wireless receiver is connected to host device with power ready to initiate pairing with transmitter chip</li> <li>Pairing search period has expired</li> </ul>	<ul style="list-style-type: none"> <li>Hold reusable chip to the indent on the wireless receiver to initiate pairing</li> <li>See <i>Instructions</i>, section b) for pairing</li> </ul>
Green	solid (2 seconds)	--	<ul style="list-style-type: none"> <li>Chip and receiver are linked</li> </ul>	<ul style="list-style-type: none"> <li>Insert reusable chip into sensor attachment strap to complete pairing</li> <li>See <i>Instructions</i>, section b) for pairing</li> </ul>
	flashing (30 seconds)	flashing (30 seconds)	<ul style="list-style-type: none"> <li>Pairing search period</li> </ul>	
Blue	solid	flashing	<ul style="list-style-type: none"> <li>Successful pairing of receiver and chip</li> </ul>	<ul style="list-style-type: none"> <li>Verify sensor attachment so host device can receive data</li> </ul>
Purple	flashing	flashing	<ul style="list-style-type: none"> <li>Battery seal tab has not been removed to activate battery</li> <li>Battery is obstructed</li> </ul>	<ul style="list-style-type: none"> <li>Remove tab to activate battery</li> <li>Refer to Directions for Use for the Radius PPG Adhesive Sensor</li> <li>Disconnect reusable chip from sensor, wait 30 seconds, insert chip into sensor (Refer to <b>Figs. 9 and 10</b>.)</li> </ul>
Orange	flashing	flashing	<ul style="list-style-type: none"> <li>Low sensor battery</li> </ul>	<ul style="list-style-type: none"> <li>Consider replacing sensor, do not discard reusable chip</li> <li>See <i>Instructions</i>, section c) for disconnecting</li> </ul>
Red	flashing	flashing	<ul style="list-style-type: none"> <li>Depleted sensor battery</li> <li>Hardware or sensor failure, chip blinking board failure code</li> </ul>	<ul style="list-style-type: none"> <li>Replace sensor, do not discard reusable chip. If issue persists, replace reusable chip</li> <li>See <i>Instructions</i>, section c) for disconnecting</li> <li>Contact Masimo Technical Support, or replace sensor and chip</li> </ul>

## ACCURACY SPECIFICATIONS (ARMS\*)

Oxygen Saturation (SpO <sub>2</sub> )		
SpO <sub>2</sub> Accuracy, No Motion <sup>1</sup> (70–100%)	Adults, Pediatrics, Infants	2%
	Neonates <sup>5</sup>	3%
SpO <sub>2</sub> Accuracy, Motion <sup>2</sup>	All patient populations	3%
SpO <sub>2</sub> Accuracy, Low Perfusion <sup>3</sup>	All patient populations	2%

Pulse Rate <sup>4</sup> (PR) 25–240 bpm	
All patient populations	
No Motion	3 bpm
Motion	5 bpm
Low Perfusion	3 bpm

\* **NOTE:** ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within  $\pm$  ARMS of the reference measurements in a controlled study.

<sup>1</sup> The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO<sub>2</sub> against a laboratory co-oximeter.

<sup>2</sup> The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70–100% SpO<sub>2</sub> against a laboratory co-oximeter.

<sup>3</sup> The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

<sup>4</sup> The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25–240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

<sup>5</sup> The saturation accuracy of the Neonate and Preterm sensors were validated on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

## ENVIRONMENTAL SPECIFICATIONS

Storage/Transport Temperature	-40 – +70°C @ ambient humidity
Operating Temperature	0–40°C @ ambient humidity
Storage/Transport Humidity	5–95% non-condensing
Operating Humidity	5–95% non-condensing
Atmospheric Pressure	540 to 1060 mBar @ ambient temperature and humidity

## BATTERY LIFE

Battery Life	96 hours in typical continuous usage
--------------	--------------------------------------

## WIRELESS TECHNOLOGY INFORMATION

Bluetooth LE Wireless Technology Information	
Modulation Type	GFSK
Max. Output Power	+8 dBm
Frequency Range	2402–2480 MHz
Antenna Peak Gain	1 dBi
Recommended Range	100 ft (~30 meters) line-of-sight
Quality of Service (QoS)	Delay < 30 seconds
Security	Proprietary binary protocol

FCC ID are as follows: Chip: VKF-AIRT01 Receiver: VKF-AIRDG01

IC IDs are as follows: Chip: 7362A-AIRT01 Receiver: 7362A-AIRDG01

EU Radio Equipment Directive (RED 2014/53/EU)		
Radius PPG Receiver:		Radius PPG Chip:
EN 301 489-17 V3.1.1	EN 300 328 V2.2.1	EN 301 489-17 V3.1.1
EN 301 489-3 V2.1.1	EN 300 330 V2.1.0	EN 300 328 V2.2.1

**WARNING:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**WARNING:** The frequency bands of this device (2.4 GHz) are only for indoor use, in accordance with international telecommunication requirements.

**CAUTION:** Keep the Radius PPG away from electrical equipment that emits radio frequencies to minimize radio interference. Radio interference may result in no or inaccurate readings.

**Note:** This device complies with part 15 of FCC Rules and Industry Canada's license-exempt RSS's. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Note:** When using Radius PPG consideration should be taken to local government frequency allocations and technical parameters to minimize the possibility of interference to/from other wireless devices.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

**RF Radiation Exposure Statement:** This equipment has been exempted from FCC RF radiation exposure testing set forth for an uncontrolled environment.

**Note:** Users are advised that high-power radars are allocated as primary users (i.e. priority users) of the bands 5250–5350 MHz and 5650–5850 MHz and that these radars could cause interference and/or damage to LE-LAN devices

**RF Radiation Exposure Statement:** Radius PPG Receiver has been exempted from IC RSS 102 RF radiation exposure limits set forth for an uncontrolled environment.

**WARRANTY**

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

**WARRANTY EXCLUSIONS**

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to products that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

**NO IMPLIED LICENSE**

PURCHASE OR POSSESSION OF THIS PRODUCT DOES NOT CARRY ANY EXPRESS OR IMPLIED LICENSE TO USE WITH ANY DEVICE THAT IS NOT AN AUTHORIZED DEVICE OR SEPARATELY AUTHORIZED TO USE MASIMO PRODUCTS.

**RECOMMENDED SEPARATION DISTANCES**

**RECOMMENDED SEPARATION DISTANCE BETWEEN PORTABLE AND MOBILE RF COMMUNICATION EQUIPMENT AND THE ME EQUIPMENT**

The ME Equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ME Equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ME Equipment as recommended below, according to the maximum output power of the communication equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (M)		
	150 kHz to 80 MHz $d = 1.17 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \cdot \sqrt{P}$	800 MHz a 2.5 GHz $d = 2.33 \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.37
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE)	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use	<b>LOT</b>	Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer	<b>REF</b>	Catalogue number (model number)	<b>EC REP</b>	Authorized representative in the European community
	Date of manufacture YYYY-MM-DD	<b>###</b>	Masimo reference number		Body weight
<b>SN</b>	Serial Number		Light Emitting Diode (LED) LED emits light when current flows through		Storage temperature range
	Do not discard	<b>&gt;</b>	Greater than		Keep dry
	Non-Sterile	<b>&lt;</b>	Less than		Do not use if package is damaged
	Not made with natural rubber latex		Storage humidity limitation		Atmospheric pressure limitation
	Caution		Bluetooth		Warning
	Fragile, handle with care	<b>FC</b>	Federal Communications Commission (FCC) Licensing	<b>FCC ID:</b>	Identifies unit has been registered as a radio device
<b>IP47</b>	Protected against solid foreign objects of 1.0 mm diameter and greater and protected against the effects of temporary immersion in water	<b>IP22</b>	Protected against solid foreign objects of 12.5 mm diameter and greater and protection against vertically falling water drops when enclosure is tilted at 15 degrees		Instructions/Directions for Use/ Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eIFU is not available in all countries.</b>

Patents: <http://www.masimo.com/patents.htm>

Radius PPG is a trademark of Masimo Corporation.

Masimo, SET, X-Cal, and  are federally registered trademarks of Masimo Corporation.

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© 2019 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH