LNCS-II™ rainbow® DCI & DCIP 8λ SpHb Series









Reusable SpHb®, SpMet® & SpO2 Sensors

Images	2
en English	3-7
fr French	8-12
de German	13-17
it Italian	18-22
es Spanish	23-27
sv Swedish	28-32
nl Dutch	33-37
da Danish	38-42
pt Portuguese	43-47
zh Chinese	48-51
ja Japanese	52-56
fi Finnish	57-61
no Norwegian	62-66
cs Czech	67-71
hu Hungarian	72-76
pl Polish	77-81
ro Romanian	82-86
sk Slovak	87-91
tr Turkish	92-96
el Greek	97-101
ru Russian	102-106
ko Korean	107-110
ar Arabic	113-111

LNCS-II™ rainbow® DCI & DCIP 8λ SpHb Series

Reusable SpHb®, SpMet® & SpO2 Sensors

LNCS-II™ rainbow® DCI & DCIP 8λ SpHb Series

en

Reusable SpHb®, SpMet® & SpO2 Sensors

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this sensor, the user should read and understand the Operator's Manual for the device and this Directions for Use.

INDICATIONS

The LNCS-II™ rainbow® DCI® and DCIP 8λ SpHb reusable sensors are indicated for spot check or continuous noninvasive monitoring of arterial saturation (SpO₂), arterial methemoglobin (SpMet®), total hemoglobin (SpHb®), and pulse rate for use with adult and pediatric patients during both motion and no motion conditions, and for patients who are well or poorly perfused in hospital, hospital-type facilities, mobile, and home environments. SpMet and SpHb accuracies were not validated under motion or low perfusion conditions.

WARNING: Masimo® sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

CONTRAINDICATIONS

LNCS-II rainbow® DCI and DCIP 8λ SpHb reusable sensors are contraindicated for use on active patients or for prolonged periods of use. The sensor site must be inspected at least every four (4) hours or sooner; and if the circulatory condition or skin integrity is compromised, the sensor should be applied to a different site.

DESCRIPTION

LNCS-II rainbow® DCI and DCIP 8λ SpHb reusable sensors are for use only with devices containing Masimo rainbow SET™ technology Version 7.4 or higher, or licensed to use rainbow compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular instruments and sensor models. Each instrument manufacturer is responsible for determining whether their instruments are compatible with each sensor model.

NOTE: Though this sensor is capable of reading all parameters, it is limited by the parameters on the device. SpHb parameter must be on the device.

NOTE: The sensor will not function beyond its specified duration of use.

The LNCS-II rainbow® DCI and DCIP 8λ SpHb reusable sensors have been verified using Masimo rainbow SET technology.

WARNINGS, CAUTIONS AND NOTES

- Laboratory diagnostic tests using blood samples should be conducted prior to clinical decision making to completely understand the patient's condition. Comparisons between SpHb measurements and laboratory diagnostic hemoglobin measurements may be affected by sample type, collection technique, physiological, and other factors.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- Circulation distal to the sensor site should be checked routinely.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor, Trendelenburg position).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid value regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.

- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- High intensity extreme lights (such as pulsating strobe lights) directed on the sensor, may not allow the Pulse CO-Oximeter to obtain vital sign readings.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Inaccurate SpCO and SpMet readings can be caused by abnormal hemoglobin levels, low arterial perfusion, low arterial oxygen saturation levels including altitude induced hypoxemia, motion artifact.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ and SpHb measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ and SpHb measurements.
- Elevated levels of total bilirubin may lead to inaccurate SpO₂, SpHb and SpMet measurements.
- Inaccurate SpO₂ and SpHb readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause erroneous SpHb readings.
- Inaccurate SpHb readings may be caused by elevated PaO₂ levels and low arterial oxygen saturation levels.
- Inaccurate SpHb readings may be caused by vasospastic disease such as Raynaud's, peripheral vascular disease, liver disease, EMI radiation interference.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate SpO₂, SpHb, and SpMet measurements.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize the sensor.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. Reusable spot-check sensors are available with different quantities of spot-check measurements. The sensor will stop functioning after the specified number of spot-check measurements have been completed.

INSTRUCTIONS

A) Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring or middle finger of the non-dominant hand is preferred.
- Alternatively, the other digits on the non-dominant hand may be used. Always choose a site that will completely cover the sensor's detector window.
- Site should be cleaned of debris and dry prior to sensor placement.

Slender Digit Gauge Instructions

- The Slender Digit Gauge on the sensor cable is an aid for selecting an appropriate digit for sensor application.
 - Use of this gauge is recommended only for patients weighing > 30 kg that have slender digits. Remove the gauge from the digit BEFORE sensor application.
1. It is preferable to start with the ring finger of the non-dominant hand. If this is too slender, select a larger digit. Slide the gauge circle on the digit. If the gauge circle stops at any point of the nail bed before the cuticle, the DCI sensor should be used on that digit (Fig A).
 2. If the gauge slides past the cuticle, the digit is too slender for this sensor. Select a different digit, or use a pediatric/slender digit sensor (DCIP) on this patient (Fig A).

B) Attaching the sensor to the patient

1. Refer to **Fig. 1**. Open the sensor by pressing on hinge tabs. Position the selected digit so that it corresponds to that of the finger shown on top of the sensor. Place the selected digit over the sensor window of the sensor. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. The tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop.
2. Refer to **Fig. 2**. The hinged tabs of the sensor should open to evenly distribute the grip pressure of the sensor along the length of the finger. Check the arrangement of the sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to **Fig. 3**. Orient the sensor so that the cable will be running towards the top of the patient's hand.

NOTE: The sensor is not intended for use across a child's hand or foot.

C) Attaching the sensor to patient cable

1. Refer to **Fig. 4a**. Properly orient the sensor connector and insert the sensor completely into the patient cable connector.
2. Refer to **Fig. 4b**. Close the protective latch completely.

NOTE: Check the remaining duration of the sensor on the device prior to the start of monitoring. Consult device user manual for additional information.

Disconnecting the sensor from patient cable

1. Refer to **Fig. 5a**. Lift up the protective latch.
2. Refer to **Fig. 5b**. Pull firmly on the sensor connector to remove it from the patient cable.

NOTE: To avoid damage, pull on the sensor connector, not the cable.

Removing the sensor from the patient

1. Open the sensor by pressing on hinge tabs. Remove the sensor from the digit and follow the cleaning instructions, if necessary.
2. Store the sensor away for its next use.

CLEANING: To clean the sensor, first remove it from the patient and disconnect it from the instrument. You may then clean the sensor by wiping it with a 70% isopropyl alcohol pad. Allow the sensor to dry prior to placement on a patient.

SPECIFICATIONS

When used with Masimo rainbow SET technology monitors or with licensed Masimo rainbow SET technology modules using RD/LNC-II rainbow patient cables, the LNCS-II rainbow® DCI and DCIP 8λ SpHb sensors have the following performance specifications:

DC Series Sensor	DCI	DCIP
Application Site	finger	finger or thumb
  Body Weight	> 30 kg	10-50 kg
SpO ₂ Accuracy, No Motion ¹	60-80% 3% 70-100% 2%	60-80% 3% 70-100% 2%
SpO ₂ Accuracy, Motion ²	3%	3%
Pulse Rate Accuracy, No Motion ³	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ³	5 bpm	5 bpm
SpO ₂ Accuracy, Low Perfusion ⁴	2%	2%
Pulse Rate Accuracy, Low Perfusion ³	3 bpm	3 bpm
SpMet Accuracy, No Motion ⁵	1%	1%
SpHb Accuracy, No Motion ⁶	1 g/dL	1 g/dL

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within +/- ARMS of the reference measurements in a controlled study.

¹The Masimo rainbow SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 60%-100% SpO₂ against a laboratory co-oximeter.

²The Masimo rainbow SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

³The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 60% to 100%.

⁴The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁵SpMet accuracy was determined by testing on healthy adult volunteers with light to dark skin pigmentation in the range of 1%-15% MetHb against a laboratory CO-Oximeter.

⁶SpHb accuracy, against that of a Coulter Counter, was determined by testing on healthy adult volunteers with light to dark skin pigmentation in the range 8 to 17 g/dL.

COMPATIBILITY

These sensors are intended for use only with devices containing Masimo rainbow SET technology or pulse oximetry monitors licensed to use rainbow compatible sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. In no event shall Masimo be liable for any damages associated a product that has been reprocessed, reconditioned or recycled. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE



PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE LNCS-II rainbow SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
 (blue background)	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code	 0123	Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Caution		Storage temperature range
	Do not discard		Greater than		Keep dry
	Non-Sterile		Less than		Do not use if package is damaged and consult instructions for use
	Not made with natural rubber latex		Storage humidity limitation		Atmospheric pressure limitation
	Fragile, handle with care		Medical device		Unique device identifier
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available in all countries.				

Patents: <http://www.masimo.com/patents.htm>

Masimo, rainbow, SET, , SpMet, rainbow SET, SpHb, X-Cal, and DCI are federally registered trademarks of Masimo Corporation.

LNCS-II is a trademark of Masimo Corporation.

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© 2020 Masimo Corporation
MDSS GmbH