



LiDCO® Module Cables

Images	2
en English	3-5
fr French	6-8
de German	9-11
it Italian	12-14
es Spanish	15-17
sv Swedish	18-20
nl Dutch	21-23
da Danish	24-26
pt Portuguese	27-29
zh Chinese	30-32
ja Japanese	33-35
fi Finnish	36-38
no Norwegian	39-41
cs Czech	42-44
hu Hungarian	45-47
pl Polish	48-50
ro Romanian	51-53
sk Slovak	54-56
tr Turkish	57-59
el Greek	60-63
ru Russian	64-67
ko Korean	68-70
ar Arabic	73-71

LiDCO®

Module Cables

Module Cables

DIRECTIONS FOR USE

Reusable

 Not made with natural rubber latex Non-Sterile

INDICATIONS

LiDCO® module cables have the same indications for use as the LiDCO Module. Please see the compatible LiDCO Module's Directions for Use for indications for use and prescribed information.

DESCRIPTION

The LiDCO Module and cables support minimally invasive hemodynamic measurements that derive nominal stroke volume and heart rate from the patient's existing arterial pressure waveform using the LiDCO algorithm. The module connects to a host device (such as Masimo® Root®) via MOC9-® connector and receives its power from the host device. The LiDCO Module cables provide connection from the LiDCO Module to a 3rd party patient monitor and/or a patient transducer.

WARNING: Masimo LiDCO cables are designed for use with modules containing or licensed to use Masimo or LiDCO technology.

WARNINGS, CAUTIONS, AND NOTES

- Always refer to the LiDCO module operator's manual for complete instructions or additional instructions.
- Check the cable before use to ensure the cable is physically intact, with no broken or frayed wires or damaged parts and discard if cracks or discoloration are found.
- Do not use the cable in the event of any damage.
- LiDCO Module cables are designed for use with specific monitors and/or transducers. Verify the compatibility of the monitor, cable and transducer before use, otherwise degraded performance and/or patient injury can result.
- Carefully route patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the monitor or the LiDCO module to the cable may result in intermittent readings, inaccurate results, or no reading.
- Always insert and remove cables using the connector. When connecting do not try to force, instead use the mark to align the connector.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- To prevent damage, do not soak or immerse the cable in any liquid solution. Do not attempt to sterilize the cable.
- Do not attempt to reprocess, recondition or recycle Masimo cables as these processes may damage the electrical components, potentially leading to patient harm.
- Always use the correct LiDCO Module cable for connection to the 3rd party patient monitor.
- Make sure that the cable is appropriately arranged to minimize the risk of injury to patients, users and the equipment.

INSTRUCTIONS

A. Connect the LiDCO cable to the LiDCO Module

1. Refer to **Fig. 1**. Connect the module connector of the LiDCO cable to the mating connector of the LiDCO Module.
2. Ensure the cable connector is fully engaged with the module. A tactile or audible click of connection should be heard.

B. Connect the LiDCO cable to a 3rd party patient monitor

1. Identify whether the cable is a LiDCO BP Module cable or an analog cable.
 - a. If the cable has 3 connectors, it is a BP Module cable (see **Fig 2**).
 - b. Otherwise, it is an analog cable (see **Fig 3**).
2.
 - a. If it is a BP Module cable, refer to **Fig. 2**. Connect one end to a cable connected to the 3rd party patient monitor as shown.
 - b. If it is an analog cable, refer to **Fig 3**. Connect the cable to a 3rd party patient monitor as shown.

3. Ensure the LiDCO cable connector is fully engaged with the connector of the 3rd party patient monitor. A tactile or audible click of connection should be heard.

C. Connect the LiDCO cable to a transducer (BP Module Cable only)

1. Orient and connect the transducer end of the LiDCO cable to the transducer.

Refer to the manufacturer's user manual for detailed instructions on transducers.

D. Disconnect the LiDCO cable

1. Pull firmly on the LiDCO cable connector to remove it from the 3rd party patient monitor or LiDCO Module/transducer connector.

Note: If using a LiDCO cable specific to a 3rd party patient monitor, use an analog cable. Different 3rd party patient monitors require a specific analog cable.

Note: If using a LiDCO cable not specific to a 3rd party patient monitor, use a BP Module Cable. Different transducers require a specific BP Module cable.

Please contact your Masimo rep for details.

CLEANING

1. Remove all connections from LiDCO cable before cleaning.

2. Clean the surface of the cable by wiping it with an approved cleaning solution:

70% isopropyl alcohol, Glutaraldehyde, 10% bleach (sodium hypochlorite) to 90% water solution, Quaternary ammonium chloride, Accelerated Hydrogen Peroxide (e.g. Oxivir Tb)

3. Wipe all surfaces of the cable.

4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.

5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

CAUTION

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

ENVIRONMENTAL

Operating Temperature 32°F to 104°F (0°C to 40°C) @ ambient humidity

Storage Temperature -40°F to 158°F (-40°C to 70°C) @ ambient humidity

Relative Humidity 10% to 95% non-condensing @ ambient temperature

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

Purchase or possession of these cables does not carry any express or implied license to use these cables with any device that is not an authorized device or separately authorized to use LiDCO cables.

COMPATIBILITY

LiDCO cables are intended for use only with compatible devices and/or transducers. Consult individual system manufacturers for compatibility of particular models. Each device manufacturer is responsible for determining whether their devices are compatible with LiDCO cables. Use of this cable with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use	LOT	Lot code	CE	European Union Conformity Mark
	Manufacturer	REF	Catalogue number (model number)	EC/REP	Authorized representative in the European community
	Date of manufacture YYYY-MM-DD	####	Masimo reference number		Body weight
	Do not discard		Caution		Storage temperature range
	Non-sterile	>	Greater than		Keep dry
	Not made with natural rubber latex	<	Less than		Do not use if package is damaged and consult instructions for use
	Fragile, handle with care		Storage humidity limitation		Atmospheric pressure limitation
MD	Medical device	UDI	Unique device identifier		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available in all countries.
	Importer		Distributor		

Patents: <http://www.masimo.com/patents.htm>

LiDCO is a trademark of Masimo Corporation.

Masimo, Root, MOC9-, and  are federally registered trademarks of Masimo Corporation.

All other products, logos, or company names mentioned herein may be trademarks and/or registered trademarks of their respective companies.

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Authorized representative for Masimo Corporation:
MDSS GmbH