




Red DCI® and DCI-P

Adult and Pediatric Reusable SpO2 Sensors

Images	2
en English	3-6
fr French	7-10
de German	11-14
it Italian	15-18
es Spanish	19-22
sv Swedish	23-26
nl Dutch	27-30
da Danish	31-34
pt Portuguese	35-38
zh Chinese	39-42
ja Japanese	43-46
fi Finnish	47-50
no Norwegian	51-54
cs Czech	55-58
hu Hungarian	59-62
pl Polish	63-66
ro Romanian	67-70
sk Slovak	71-74
tr Turkish	75-78
el Greek	79-82
ru Russian	83-86
ko Korean	87-90
ar Arabic	94-91

Red DCI® and DCI-P

Adult and Pediatric Reusable SpO₂ Sensors

en

DIRECTIONS FOR USE

Reusable

Do not discard

Not made with natural rubber latex

Non-sterile

Prior to using this sensor, the user should read and understand the Operator's Manual for the device or monitor and this Directions for Use.

INDICATIONS

Red reusable DCI® and DCI-P sensors are indicated for either spot "check" or the continuous non-invasive monitoring of arterial oxygen saturation (SpO₂) and pulse rate.

CONTRAINDICATIONS

Red DCI and DCI-P sensors are contraindicated for use on mobile patients or for prolonged periods of use. They must be removed and repositioned to a different monitoring site at least every four (4) hours.

DESCRIPTION

Red DCI and DCI-P sensors are for use only with devices containing Masimo® rainbow SET® Technology or licensed to use rainbow® compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor model.

The DCI and DCI-P sensors have been verified with Masimo rainbow SET Oximetry Technology.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- The site must be checked at least every four (4) hours to ensure proper adhesion, skin integrity and proper alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every two (2) hours with poorly perfused patients.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the readings may read lower than core arterial oxygen saturation.
- Misapplied sensors or sensors that become partially dislodged may cause either over or under reading of actual arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage or damage the sensor.
- The sensor should be free of visible defects. Never use a damaged sensor or one with exposed electrical circuitry.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Intravascular dyes or externally applied coloring (such as nail polish) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Failure to apply the sensor properly may cause incorrect measurements.
- Do not use the sensor during MRI scanning.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should be not below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor)
- Venous pulsations may cause erroneous low readings (e.g. tricuspid valve regurgitation).
- Do not modify or alter the sensor in any way. Alterations or modification may affect performance and/or accuracy.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.

- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A. Site selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring or middle finger of the non-dominant hand is preferred.
- Alternatively, the other digits on the non-dominant hand may be used. Always choose a site that will completely cover the sensor's detector window. The great toe or long toe (next to the great toe) may be used on restrained patients or patients whose hands are unavailable.
- Site should be cleaned of debris and dry prior to sensor placement.

B. Attaching the sensor to the patient

1. Refer to **Fig. 1**. Open the sensor by pressing on hinge tabs. Place the selected digit over the sensor window of the sensor. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. The tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop.
2. Refer to **Fig. 2**. The hinged tabs of the sensor should open to evenly distribute the grip of the sensor along the length of the finger. Check position of sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to **Fig. 3**. Orient the sensor so that the cable will be running towards the top of the patient's hand. Connect the sensor connector to a patient cable or directly to the device.

C. Disconnecting the Sensor from the Patient Cable

1. Place thumb and index finger on gray buttons on either side of the patient cable connector.
2. Refer to **Fig. 4**. Press firmly on the gray buttons and pull to remove the sensor.

CLEANING

To clean the sensor, first remove it from the patient and disconnect it from the patient cable. You may then clean the sensor by wiping it with a 70% isopropyl alcohol pad. Allow the sensor to dry prior to placement on a patient.

CAUTION: To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not sterilize by irradiation, steam, autoclave or ethylene oxide.

SPECIFICATIONS

When used with Masimo rainbow SET pulse oximetry monitors or with licensed Masimo rainbow SET pulse oximetry modules Red DCI and DCI-P reusable sensors have the following specifications:

Sensor	Red DCI	Red DCI-P
  Body Weight	> 30 kg	10 - 50 kg
Application site	Finger	Finger
SpO ₂ Accuracy, No Motion (70 - 100%) ¹	2%	2%
SpO ₂ Accuracy, Motion ²	3%	3%
SpO ₂ Accuracy, Low Perfusion ³	2%	2%
Pulse Rate Accuracy, No Motion (25 - 240 bpm) ⁴	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ⁴	5 bpm	5 bpm
Pulse Rate Accuracy, Low Perfusion ⁴	3 bpm	3 bpm

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within \pm ARMS of the reference measurements in a controlled study.

¹The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory co-oximeter. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

²The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory co-oximeter. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

³The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

⁴The Masimo SET Technology has been validated for pulse rate accuracy for the range of 30-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

COMPATIBILITY

 These sensors are intended for use only with devices containing Masimo rainbow SET technology or pulse oximetry monitors licensed to use rainbow compatible sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

ENVIRONMENTAL

Operating Temperature	5°C to 40°C (41°F to 104°F)
Storage Temperature	-40°C to 70°C (-40°F to 158°F)
Relative Humidity	5% to 95% noncondensing

WARRANTY

Masimo warrants to the initial buyer only that each product it manufactures, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE


PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RED DCI AND DCI-P SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Caution		Storage temperature range
	Do not discard		Greater than		Keep dry
	Non-Sterile		Less than		Do not use if package is damaged and consult instructions for use
	Not made with natural rubber latex		Storage humidity limitation		Atmospheric pressure limitation
	Fragile, handle with care		Medical device		Unique device identifier
	Importer		Distributor		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available in all countries.

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , rainbow SET, X-Cal, and DCI are federally registered trademarks of Masimo Corporation.

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© 2021 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH