






# LNOP® MP Series Patient Cables

Images	2
<b>en</b> English	3-6
<b>fr</b> French	7-10
<b>de</b> German	11-14
<b>it</b> Italian	15-18
<b>es</b> Spanish	19-22
<b>sv</b> Swedish	23-26
<b>nl</b> Dutch	27-30
<b>da</b> Danish	31-34
<b>pt</b> Portuguese	35-38
<b>zh</b> Chinese	39-42
<b>ja</b> Japanese	43-46
<b>fi</b> Finnish	47-50
<b>no</b> Norwegian	51-54
<b>cs</b> Czech	55-58
<b>hu</b> Hungarian	59-62
<b>pl</b> Polish	63-66
<b>ro</b> Romanian	67-70
<b>sk</b> Slovak	71-74
<b>tr</b> Turkish	75-78
<b>el</b> Greek	79-82
<b>ru</b> Russian	83-86
<b>ko</b> Korean	87-90
<b>ar</b> Arabic	94-91

# LNOP<sup>®</sup> MP Series Patient Cables

## DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

### INDICATIONS

The LNOP<sup>®</sup> MP Series Patient Cables and Masimo Sensors are indicated for the continuous noninvasive monitoring of arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate for adult, pediatric, infant and neonatal patients.

### WARNINGS

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Failure to properly connect the sensor or the oximeter module to the cable will result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Refer to the specific sensor's *Directions for Use* for detailed sensor instructions and warnings.
- Probes and cables are designed for use with specific monitors. The user and/or operator needs to verify the compatibility of the monitor, probe, and cable before use, otherwise patient injury can result.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.

### INSTRUCTIONS

#### Attaching LNOP MP Patient Cables to the Instrument

Refer to Figures 1a and 1b. Orient the connector of the cable to mate with the patient cable connector on the SpO<sub>2</sub> Module (**Fig. 1a**) or the MMS (**Fig. 1b**).

**CAUTION: ENSURE YOU ARE CONNECTING AN LNOP MP PATIENT CABLE TO EITHER A MASIMO SET OR PHILIPS FAST SpO<sub>2</sub> SOCKET. THERE ARE A NUMBER OF SIMILAR CONNECTORS WITH DIFFERENT COLORS AND DIFFERENT MECHANICAL KEYING. NEVER FORCE A PATIENT CABLE CONNECTOR INTO THE MODULE OR MMS. FAILURE TO USE AN LNOP MP PATIENT CABLE MAY RESULT IN DAMAGE TO THE MODULE OR MMS, INACCURATE READINGS, OR NO READINGS.**

**Note:** There is a difference in color shade between the LNOP MP Patient Cable and the Philips MMS/SpO<sub>2</sub> Module receptacle. However, this is an acceptable configuration. Masimo SET performance is not available when connected to an instrument with Philips FAST SpO<sub>2</sub> Technology.

#### Attaching the LNOP MP Patient Cable to the Masimo Sensor

1. Refer to **Fig. 2**. Orient the connecting tab of the sensor so that the "shiny" contacts are pointed up. Mate the logo on the sensor to the logo on the patient cable.
2. Refer to **Fig. 3**. Insert the tab of the sensor into the patient cable connector until there is a tactile or audible click of connection. Gently tug on the patient cable connector to ensure a positive contact. Tape may be used to secure the cable to the patient for ease of movement.

#### Disconnecting the LNOP MP Patient Cable from the Masimo Sensor

1. Refer to **Fig. 4**. Place thumb and index finger on the buttons on either side of the patient cable connector.
2. Press firmly on the buttons and pull to remove the sensor.

### CLEANING

Clean the LNOP MP Patient Cable by wiping it with a 70% isopropyl alcohol pad and allow it to dry.

### CAUTION

Do not soak or immerse the cable in any liquid solution. Do not attempt to sterilize.

## SPECIFICATIONS

LNOP MP Patient Cables are for use with the following devices and sensors:

When used with Masimo SET® Technology:

- Masimo SET Intellivue Module
- Masimo SET Intellivue Multi-Measurement Server (MMS)
- Masimo SET Intellivue Multi-Measurement Server X2 (MMS X2)

Sensors	Weight Range	Saturation Accuracy (70-100% SpO <sub>2</sub> )		Pulse Rate Accuracy (25-240 bpm)		Low Perfusion Accuracy (70-100% SpO <sub>2</sub> )	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNOP Adt/Adtx	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP Pdt/Pdtx	10-50 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP Inf-L	3-20 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP Neo-L	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
	> 40 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP NeoPt-L	< 1 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
LNOP Trauma	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP Newborn Neonatal	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
LNOP Newborn Infant/Pediatric	3-10 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
	10-30 kg						
LNOP DCI	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP DCIP	10-50 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP DC-195	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP YI	1-3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
	> 3 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
LNOP TC-I	> 30 kg	3.5%	N/A	3 bpm	N/A	3.5%	3 bpm
LNOP TF-I®	> 30 kg	2%	N/A	3 bpm	N/A	N/A	N/A

Sensors	Weight Range	Saturation Accuracy (SpO <sub>2</sub> )			Pulse Rate Accuracy (25-240 bpm)	Low Perfusion Accuracy (70-100% SpO <sub>2</sub> )	
		No Motion				Saturation	Pulse Rate
		80-100%	70-100%	60-80%			
LNOP Blue	2.5-30 kg	3%	3.3%	4%	3 bpm	3%	3 bpm

The LNOP TF-I and TC-I sensors were not validated under motion conditions.

The LNOP YI was not validated under low perfusion conditions.

**When used with Philips FAST-SpO<sub>2</sub> Technology:**

- FAST-SpO<sub>2</sub> Intellivue Module M1020B #A01
- FAST-SpO<sub>2</sub> M3/M4 Multi-Measurement Server (MMS) M3000A #A01
- FAST-SpO<sub>2</sub> Intellivue Multi-Measurement Server (MMS) M3001A #A01
- FAST-SpO<sub>2</sub> Intellivue Multi-Measurement Server X2 (MMS X2) M3002A #A01
- FAST-SpO<sub>2</sub> Intellivue MP2, M8102A
- FAST-SpO<sub>2</sub> Intellivue MP5/MP5T/MP5SC, M8105A/AT/AS

Sensors	Weight Range	Saturation Accuracy (70-100% SpO <sub>2</sub> )	Pulse Rate Accuracy (25-240 bpm)
		No Motion	No Motion
LNOP Adt/Adtx	> 30 kg	2%	3 bpm
LNOP Pdt/Pdtx	10-50 kg	2%	3 bpm
LNOP Inf-L	3-20 kg	2%	3 bpm
LNOP Neo-L	< 3 kg	3%	3 bpm
	> 40 kg	2%	3 bpm
LNOP NeoPt-L	< 1 kg	3%	3 bpm
LNOP DCI	> 30 kg	2%	3 bpm
LNOP DCIP	10-50 kg	2%	3 bpm
LNOP DC-195	> 30 kg	2%	3 bpm
LNOP YI	1-3 kg	3%	3 bpm
	> 3 kg	2%	3 bpm
LNOP TC-I	> 30 kg	3.5%	3 bpm
LNOP Trauma <sup>1</sup>	> 30 kg	2%	3 bpm
LNOP Newborn Neonatal <sup>1</sup>	< 3 kg	3%	3 bpm
LNOP Newborn Infant/Pediatric <sup>1</sup>	3-10 kg	2%	3 bpm
	10-30 kg		

<sup>1</sup>Provides saturation and pulse rate values when used with Philips technology but does not operate in the specialty mode unless used with Masimo technology.

ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within  $\pm$  ARMS of the reference measurements in a controlled study. LNOP MP Series Patient Cables and the above Masimo Sensors have been validated with the Masimo SET and Philips FAST-SpO<sub>2</sub> devices listed above. The saturation accuracy of the neonatal sensors were validated on adult volunteers and 1% was added to account for the properties of fetal hemoglobin. The LNOP YI multisite sensor was validated on adult digits.

**ENVIRONMENTAL**

Operating Temperature	41°F to 104°F (5°C to 40°C)
Storage Temperature	-40°F to 158°F (-40°C to 70°C)
Relative Humidity	5% to 95% noncondensing

**WARRANTY**

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

## NO IMPLIED LICENSE

Purchase or possession of this patient cable does not carry any express or implied license to use this cable with any device that is not an authorized device or separately authorized to use LNOP MP Patient Cables.

## CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
 (blue background)	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use	<b>LOT</b>	Lot code	<b>CE</b>	Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer	<b>REF</b>	Catalogue number (model number)	<b>EC REP</b>	Authorized representative in the European community
	Date of manufacture YYYY-MM-DD	<b>####</b>	Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Storage humidity limitation		Storage temperature range
	Do not discard		Do not use if package is damaged		Keep dry
	Pediatric patient	<b>&gt;</b>	Greater than	<b>&lt;</b>	Less than
	Non-Sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eIFU is not available for CE mark countries.</b>				

Patents: <http://www.masimo.com/patents.htm>

Masimo,  LNOP, TF-I, DCI and SET are federally registered trademarks of Masimo Corporation.

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© 2017 Masimo Corporation
PHILIPS is a trademark of Koninklijke Philips Electronics, N.V.
MDSS GmbH