






LNOP® MP Series Patient Cables

Images	2
en English	3-6
fr French	7-10
de German	11-14
it Italian	15-18
es Spanish	19-22
sv Swedish	23-26
nl Dutch	27-30
da Danish	31-34
pt Portuguese	35-38
zh Chinese	39-42
ja Japanese	43-46
fi Finnish	47-50
no Norwegian	51-54
cs Czech	55-58
hu Hungarian	59-62
pl Polish	63-66
ro Romanian	67-70
sk Slovak	71-74
tr Turkish	75-78
el Greek	79-82
ru Russian	83-86
ko Korean	87-90
ar Arabic	94-91

LNOP[®] MP Series Patient Cables

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this Cable, the user should read and understand the Operator's Manual for the Device, this Directions for Use and the Directions for Use for the Sensor.

INDICATIONS

The LNOP[®] MP Series Patient Cables and Masimo Sensors are indicated for continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by a SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients in hospitals, hospital-type facilities, mobile, and home environments.

DESCRIPTION

The LNOP MP Series patient cables are for use on Philips modules with Masimo SET[®] Technology and devices with Philips FAST-SpO₂ Technology.

On Philips modules with Masimo SET Technology, LNOP MP Patient Cables are for use with LNOP sensors.

On devices with Philips FAST-SpO₂ Technology, LNOP MP Patient Cables are for use with LNOP sensors. Masimo SET Technology is not available on devices with Philips FAST-SpO₂ Technology.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET[®] oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully route patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the sensor or the oximeter module to the cable will result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- To prevent damage, do not soak or immerse the sensor or cable in any liquid solution. Do not attempt to sterilize the cable or sensor.
- Do not attempt to reprocess, recondition or recycle Masimo patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Refer to the specific sensor's *Directions for Use* for detailed sensor instructions and warnings.
- Probes and cables are designed for use with specific monitors. The user and/or operator needs to verify the compatibility of the monitor, probe, and cable before use, otherwise patient injury can result.
- **Caution:** Replace the cable when a replace cable or equivalent message is displayed consistently. Consult monitoring device operator's manual for more information.
- **Note:** The cable is provided with X-Cal[®] technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The cable will provide up to 17,520 hours of patient monitoring time when connected to a monitoring device that supports X-Cal technology. Replace the cable when the patient monitoring time is exhausted. Some legacy monitoring devices do not distinguish between the cable and the sensor. In these devices, a replace sensor or similar message could pertain either to the sensor or the cable.

INSTRUCTIONS

Attaching LNOP MP Patient Cables to the Instrument

Refer to Figures 1a and 1b. Orient the connector of the cable to mate with the patient cable connector on the SpO₂ Module (**Fig. 1a**) or the MMS (**Fig. 1b**).

CAUTION: ENSURE YOU ARE CONNECTING AN LNOP MP PATIENT CABLE TO EITHER A MASIMO SET OR PHILIPS FAST SPO₂ SOCKET. THERE ARE A NUMBER OF SIMILAR CONNECTORS WITH DIFFERENT COLORS AND DIFFERENT MECHANICAL KEYING. NEVER FORCE A PATIENT CABLE CONNECTOR INTO THE MODULE OR MMS. FAILURE TO USE AN LNOP MP PATIENT CABLE MAY RESULT IN DAMAGE TO THE MODULE OR MMS, INACCURATE READINGS, OR NO READINGS.

Note: There is a difference in color shade between the LNOP MP Patient Cable and the Philips MMS/SpO₂ Module receptacle. However, this is an acceptable configuration. Masimo SET performance is not available when connected to an instrument with Philips FAST SpO₂ Technology.

Attaching the LNOP MP Patient Cable to the Masimo Sensor

1. Refer to **Fig. 2**. Orient the connecting tab of the sensor so that the «shiny» contacts are pointed up. Mate the logo on the sensor to the logo on the patient cable.
2. Refer to **Fig. 3**. Insert the tab of the sensor into the patient cable connector until there is a tactile or audible click of connection. Gently tug on the patient cable connector to ensure a positive contact. Tape may be used to secure the cable to the patient for ease of movement.

Disconnecting the LNOP MP Patient Cable from the Masimo Sensor

1. Refer to **Fig. 4**. Place thumb and index finger on the buttons on either side of the patient cable connector.
2. Press firmly on the buttons and pull to remove the sensor.

SPECIFICATIONS

LNOP MP Patient Cables are for use with the following devices and sensors:

When used with Masimo SET® Technology:

Sensors	 Body Weight	Saturation Accuracy (70-100% SpO ₂)			Pulse Rate Accuracy (25-240 bpm)		
		No Motion	Motion	Low Perfusion	No Motion	Motion	Low Perfusion
LNOP Adt/Adtx	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP Pdt/Pdtx	10-50 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP Inf-L	3-20 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP Neo-L	< 3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
	> 40 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP NeoPt-L	< 1 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
LNOP Trauma	> 30 kg	2%	3%	N/A	3 bpm	5 bpm	N/A
LNOP Newborn Neonatal	< 3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
LNOP Newborn Infant/ Pediatric	3-10 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
	10-30 kg						
LNOP Blue	2.5-30 kg	3% (for 80-100)	N/A	3%	3 bpm	N/A	3 bpm
		4% (for 60-80)					
		3.3% (for 70-100)					
LNOP DC-195	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP YI	1-3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
	> 3 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNOP TC-I	> 30 kg	3.5%	N/A	3.5%	3 bpm	N/A	3 bpm
LNOP TF-1®	> 30 kg	2%	N/A	N/A	3 bpm	N/A	N/A

When used with Philips FAST-SpO₂ Technology:

Sensors	 Body Weight	Saturation Accuracy (70-100% SpO ₂)	Pulse Rate Accuracy (25-240 bpm)
		No Motion	No Motion
LNOP Adt/Adtx	> 30 kg	2%	3 bpm
LNOP Pdt/Pdtx	10-50 kg	2%	3 bpm
LNOP Inf-L	3-20 kg	2%	3 bpm
LNOP Neo-L	< 3 kg	3%	3 bpm
	> 40 kg	2%	3 bpm
LNOP NeoPt-L	< 1 kg	3%	3 bpm
LNOP DCI	> 30 kg	2%	3 bpm
LNOP DCI-P	10-50 kg	2%	3 bpm
LNOP DC-195	> 30 kg	2%	3 bpm
LNOP YI	1-3 kg	3%	3 bpm
	> 3 kg	2%	3 bpm
LNOP TC-I	> 30 kg	3.5%	3 bpm
LNOP Trauma ¹	> 30 kg	2%	3 bpm
LNOP Newborn Neonatal ¹	< 3 kg	3%	3 bpm
LNOP Newborn Infant/ Pediatric ¹	3-10 kg	2%	3 bpm
	10-30 kg		

¹Provides saturation and pulse rate values when used with Philips technology but does not operate in the specialty mode unless used with Masimo technology.

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within \pm ARMS of the reference measurements in a controlled study. The saturation accuracy of the neonatal sensors were verified on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad.
3. Wipe all surfaces of the cable.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.
5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

CAUTION

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

ENVIRONMENTAL

Operating Temperature	41°F to 104°F (5°C to 40°C)
Storage Temperature	-40°F to 158°F (-40°C to 70°C)
Relative Humidity	5% to 95% noncondensing

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PPRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.



NO IMPLIED LICENSE

Purchase or possession of this patient cable does not carry any express or implied license to use this cable with any device that is not an authorized device or separately authorized to use LNOP MP Patient Cables.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
 (blue background)	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use	LOT	Lot code	CE	Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer	REF	Catalogue number (model number)	EC REP	Authorized representative in the European community
	Date of manufacture YYYY-MM-DD	####	Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Storage humidity limitation		Storage temperature range
	Do not discard		Do not use if package is damaged		Keep dry
	Pediatric patient	>	Greater than	<	Less than
	Non-Sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available in all countries.				

Patents: <http://www.masimo.com/patents.htm>

Masimo, , LNOP, TF-I, DCI and SET are federally registered trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
PHILIPS is a registered trademark of Koninklijke Philips N.V.
MDSS GmbH