





# LNC MP Series Patient Cables

Images	2
<b>en</b> English	3-6
<b>fr</b> French	7-10
<b>de</b> German	11-14
<b>it</b> Italian	15-18
<b>es</b> Spanish	19-22
<b>sv</b> Swedish	23-26
<b>nl</b> Dutch	27-30
<b>da</b> Danish	31-34
<b>pt</b> Portuguese	35-38
<b>zh</b> Chinese	39-42
<b>ja</b> Japanese	43-46
<b>fi</b> Finnish	47-50
<b>no</b> Norwegian	51-54
<b>cs</b> Czech	55-58
<b>hu</b> Hungarian	59-62
<b>pl</b> Polish	63-66
<b>ro</b> Romanian	67-70
<b>sk</b> Slovak	71-74
<b>tr</b> Turkish	75-78
<b>el</b> Greek	79-82
<b>ru</b> Russian	83-86
<b>ko</b> Korean	87-90
<b>ar</b> Arabic	94-91

# LNC MP Series Patient Cables

## DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

**Prior to using this Cable, the user should read and understand the Operator's Manual for the Device, this Directions for Use and the Directions for Use for the Sensor.**

### INDICATIONS

The LNC MP Series Patient Cables and Masimo Sensors are indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, infant, and neonatal patients in hospitals, hospital-type facilities, mobile, and home environments.

### DESCRIPTION

The LNC MP Series patient cables are for use on Philips modules with Masimo SET® Technology and devices with Philips FAST-SpO<sub>2</sub> Technology.

On Philips modules with Masimo SET Technology, LNC MP Patient Cables are for use with LNCS sensors. On devices with Philips FAST-SpO<sub>2</sub> Technology, LNC MP Patient Cables are for use with LNCS sensors. Masimo SET Technology is not available on devices with Philips FAST-SpO<sub>2</sub> Technology.

**WARNING:** Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

### WARNINGS, CAUTIONS, AND NOTES

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully route patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the sensor or the oximeter module to the cable will result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- To prevent damage, do not soak or immerse the sensor or cable in any liquid solution. Do not attempt to sterilize the cable or sensor.
- Do not attempt to reprocess, recondition or recycle Masimo patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Refer to the specific sensor's *Directions for Use* for detailed sensor instructions and warnings.
- Probes and cables are designed for use with specific monitors. The user and/or operator needs to verify the compatibility of the monitor, probe, and cable before use, otherwise patient injury can result.
- **Caution:** Replace the cable when a replace cable or equivalent message is displayed consistently. Consult monitoring device operator's manual for more information.
- **Note:** The cable is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The cable will provide up to 17,520 hours of patient monitoring time when connected to a monitoring device that supports X-Cal technology. Replace the cable when the patient monitoring time is exhausted. Some legacy monitoring devices do not distinguish between the cable and the sensor. In these devices, a replace sensor or similar message could pertain either to the sensor or the cable.

### INSTRUCTIONS

#### A. Attaching LNC MP Patient Cables to the Device

Refer to Fig. 1a and 1b. Orient the connector of the cable to mate with the patient cable connector on the SpO<sub>2</sub> Module (Fig. 1a) or the MMS (Fig. 1b).

**CAUTION:** ENSURE YOU ARE CONNECTING A LNC MP PATIENT CABLE TO EITHER A MASIMO SET OR PHILIPS FAST SPO<sub>2</sub> SOCKET. THERE ARE A NUMBER OF SIMILAR CONNECTORS WITH DIFFERENT COLORS AND DIFFERENT MECHANICAL KEYING. NEVER FORCE A PATIENT CABLE CONNECTOR INTO THE MODULE OR MMS. FAILURE TO USE A LNC MP PATIENT CABLE MAY RESULT IN DAMAGE TO THE MODULE OR MMS, INACCURATE READINGS, OR NO READINGS.

**Note:** There is a difference in color shade between the LNC MP Patient Cable and the Philips MMS/SpO<sub>2</sub> Module receptacle. However, this is an acceptable configuration. Masimo SET performance is not available when connected to an instrument with Philips FAST SpO<sub>2</sub> Technology.

## B. Attaching the LNC MP Patient Cable to a LNCS Sensor

1. Refer to Fig. 2. Insert the Sensor cable connector completely into the LNC MP Patient Cable.
2. Refer to Fig. 3. Completely close the protective cover.

## C. Disconnecting the LNC MP Patient Cable from a LNCS Sensor

1. Refer to Fig. 4. Lift the protective cover to gain access to the sensor connector.
2. Refer to Fig. 5. Pull firmly on the sensor connector to remove from the patient cable.

## SPECIFICATIONS

LNC MP Patient Cables are indicated for use with the following devices and sensors:

When used with Masimo SET® Technology:

Sensors	 Body Weight	Saturation Accuracy (70–100% SpO <sub>2</sub> )			Pulse Rate Accuracy (25–240 bpm)		
		No Motion	Motion	Low Perfusion	No Motion	Motion	Low Perfusion
LNCS Aidx/Aidx-3	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS Pidx/Pidx-3	10–50 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS Inf/Inf-3/Inf-L	3–20 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS Neo/Neo-3/Neo-L	< 3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
	> 40 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS NeoPt/NeoPt-3/NeoPt-L/ NeoPt-500	< 1 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
LNCS Trauma	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS Newborn Neonatal	< 3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
LNCS Newborn Infant/Pediatric	3–10 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
	10–30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS DCI	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS DCIP	10–50 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS TC-1	> 30 kg	3.5%	N/A	3.5%	3 bpm	N/A	3 bpm
LNCS TF-1*	> 30 kg	2%	N/A	N/A	3 bpm	N/A	N/A
LNCS YI	1–3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
	> 3 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
LNCS DBI	> 30 kg	2%	N/A	2%	3 bpm	N/A	3 bpm
LNCS E1	> 30 kg	2.5%	N/A	2%	3 bpm	N/A	3 bpm
LNCS TFA-1	> 10 kg	2%	N/A	2%	3 bpm	N/A	3 bpm
LNCS S-ROS 3U	universal	See applicable S-DOS sensor below for accuracy specifications when used in combination with S-ROS optical sensor.					
ReSposable S-DOS 25	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 20	10–50 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 25L	< 3 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
	> 30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 20L	3–10 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
	10–30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 30L	3–10 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
	> 40 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 25Pt-L	< 1 kg	3%	3%	3%	3 bpm	5 bpm	3 bpm
ReSposable S-DOS 20L-500	3–10 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm
	10–30 kg	2%	3%	2%	3 bpm	5 bpm	3 bpm

When used with Philips FAST SpO<sub>2</sub> Technology:

Sensors	 Body Weight	Saturation Accuracy (70–100% SpO <sub>2</sub> )	Pulse Rate Accuracy (25–240 bpm)
		No Motion	No Motion
LNCS Adtx/Adtx-3	> 30 kg	2%	3 bpm
LNCS Pdtx/Pdtx-3	10–50 kg	2%	3 bpm
LNCS Inf/Inf-3/Inf-L	3–20 kg	2%	3 bpm
LNCS Neo/Neo-3/Neo-L	< 3 kg	3%	3 bpm
	> 40 kg	2%	3 bpm
LNCS NeoPt/NeoPt-3/NeoPt-L	< 1 kg	3%	3 bpm
LNCS Trauma <sup>1</sup>	> 30 kg	2%	3 bpm
LNCS Newborn Neonatal <sup>1</sup>	< 3 kg	3%	3 bpm
LNCS Newborn Infant/Pediatric <sup>1</sup>	3–10 kg	2%	3 bpm
	10–30 kg	2%	3 bpm
LNCS DCI	> 30 kg	2%	3 bpm
LNCS DCIP	10–50 kg	2%	3 bpm
LNCS TC-I	> 30 kg	3.5%	3 bpm
LNCS YI	1–3 kg	3%	3 bpm
	> 3 kg	2%	3 bpm
LNCS DBI <sup>2</sup>	> 30 kg	2%	3 bpm

<sup>1</sup>Provides saturation and pulse rate values when used with Philips technology but does not operate in the specialty mode unless used with Masimo technology.

<sup>2</sup>DBI sensors have not been tested to be compatible with obsolete Philips products (Philips 24CT, CMS, M3).

**NOTE:** ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within  $\pm$  ARMS of the reference measurements in a controlled study. The saturation accuracy of the neonatal sensors were verified on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

## CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad.
3. Wipe all surfaces of the cable.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.
5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

## CAUTION

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

## ENVIRONMENTAL

Operating Temperature	41°F to 104°F (5°C to 40°C)
Storage Temperature	-40°F to 158°F (-40°C to 70°C)
Relative Humidity	5% to 95% noncondensing

## WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

## NO IMPLIED LICENSE




Purchase or possession of this patient cable does not carry any express or implied license to use this cable with any device that is not an authorized device or separately authorized to use LNC MP Patient Cables.

## CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use	<b>LOT</b>	Lot code	<b>CE</b>	Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer	<b>REF</b>	Catalogue number (model number)	<b>EC REP</b>	Authorized representative in the European community
	Date of manufacture YYYY-MM-DD	<b>####</b>	Masimo reference number		Body weight
	Caution		Storage humidity limitation		Storage temperature range
	Do not discard		Do not use if package is damaged and consult instructions for use		Keep dry
	Fragile, handle with care	<b>&gt;</b>	Greater than	<b>&lt;</b>	Less than
	Non-sterile		Not made with natural rubber latex		Atmospheric pressure limitation
<b>MD</b>	Medical device	<b>UDI</b>	Unique device identifier		
	Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eIFU is not available in all countries.</b>				

Patents: <http://www.masimo.com/patents.htm>

Masimo, , SET, X-Cal, TF-I and LNCs are federally registered trademarks of Masimo Corporation.

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© 2020 Masimo Corporation
PHILIPS is a registered trademark of Koninklijke Philips N.V.
EU Authorized Representative for Masimo Corporation:
MDSS GmbH