rainbow® DC Series

Reusable SpCO®, SpMet® & SpO2 Sensors






Images	2
en English	3-7
fr French	8-12
de German	13-17
it Italian	18-22
es Spanish	23-27
sv Swedish	28-32
nl Dutch	33-37
da Danish	38-42
pt Portuguese	43-47
zh Chinese	48-51
ja Japanese	52-56
fi Finnish	57-61
no Norwegian	62-66
cs Czech	67-71
hu Hungarian	72-76
pl Polish	77-81
ro Romanian	82-86
sk Slovak	87-91
tr Turkish	92-96
el Greek	97-101
ru Russian	102-106
ko Korean	107-110
ar Arabic	114-111

rainbow® DC Series

Reusable SpCO®, SpMet® & SpO2 Sensors

rainbow[®] DC Series

Reusable SpCO[®], SpMet[®] & SpO₂ Sensors

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

> 10 kg

Prior to using this sensor, the user should read and understand the Operator's Manual for the device and this Directions for Use.

INDICATIONS

rainbow[®] DCI[®] and DCI-P reusable sensors are indicated for spot check or continuous non-invasive monitoring of arterial saturation (SpO₂), arterial carboxyhemoglobin saturation (SpCO[®]), arterial methemoglobin (SpMet[®]), and pulse rate. SpCO and SpMet accuracies were not validated under motion or low perfusion conditions.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET[®] oximetry or licensed to use Masimo sensors.

CONTRAINDICATIONS

rainbow DCI and DCI-P reusable sensors are contraindicated for use on mobile patients or for prolonged periods of use. The sensor site must be inspected at least every four (4) hours or sooner; and if the circulatory condition or skin integrity is compromised, the sensor should be applied to a different site.

DESCRIPTION

rainbow DCI and DCI-P reusable sensors are for use only with devices containing Masimo rainbow SET[®] technology Version 7.1 or higher, devices containing Masimo SET technology or licensed to use rainbow compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular devices and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor model.

NOTE: Though this sensor is capable of reading all parameters, it is limited by the parameters on the device.

The rainbow DCI and DCI-P reusable sensors have been validated on the Radical-7[®] Pulse CO-Oximeter[®] containing Masimo rainbow SET technology.

WARNINGS, CAUTIONS AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- Circulation distal to the sensor site should be checked routinely.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor, Trendelenburg position).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid value regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.

- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- High intensity extreme lights (such as pulsating strobe lights) directed on the sensor, may not allow the Pulse CO-Oximeter to obtain vital sign readings.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High levels of COHb or MethHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MethHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Inaccurate SpCO and SpMet readings can be caused by abnormal hemoglobin levels, low arterial perfusion, low arterial oxygen saturation levels including altitude induced hypoxemia, motion artifact.
- SpCO readings may not be provided if there are low arterial oxygen saturation levels or elevated methemoglobin levels
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated Total Bilirubin levels may lead to inaccurate SpO₂ measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate SpO₂, SpCO and SpMet measurements.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize the sensor.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal[®] technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

CABLE AND SENSOR SELECTION

NOTE: The maximum length for the sensor, cable and extension cable combination is 15 feet. Refer to the table below for correct examples:

Cable	Sensor	
	3 feet	> 3 feet
Length	3 feet	> 3 feet
RC-1	Yes	Yes
RC-4	Yes	Yes
RC-12	No	No

INSTRUCTIONS

A. Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring or middle finger of the non-dominant hand is preferred.
- Alternatively, the other digits on the non-dominant hand may be used. Always choose a site that will completely cover the sensor's detector window.
- Site should be cleaned of debris and dry prior to sensor placement.

Slender Digit Gauge Instructions

- The Slender Digit Gauge on the sensor cable is an aid for selecting an appropriate digit for sensor application.
- Use of this gauge is recommended only for patients weighing > 30 kg that have slender digits. Remove the gauge from the digit BEFORE sensor application.
 1. It is preferable to start with the ring finger of the non-dominant hand. If this is too slender, select a larger digit. Slide the gauge circle on the digit. If the gauge circle stops at any point of the nail bed before the cuticle, the DCI sensor should be used on that digit (Fig A).
 2. If the gauge slides past the cuticle, the digit is too slender for this sensor. Select a different digit, or use a pediatric/slender digit sensor (DCI-P) on this patient (Fig A).

B. Attaching the sensor to the patient

1. Refer to Fig 1. Open the sensor by pressing on hinge tabs. Position the selected digit so that it corresponds to that of the finger shown on top of the sensor. Place the selected digit over the sensor window of the sensor. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. The tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop.
2. Refer to Fig 2. The hinged tabs of the sensor should open to evenly distribute the grip pressure of the sensor along the length of the finger. Check the arrangement of the sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to Fig 3. Orient the sensor so that the cable will be running towards the top of the patient's hand.

NOTE: With smaller digits, in order to completely cover the detector window, it may be necessary to use a lower weight range sensor. The sensor is not intended for use on the thumb or across a child's hand or foot.

C. Attaching the sensor to patient cable

1. Refer to Fig 4. Orient sensor connector (1) to patient cable connector (2) as shown.
2. Insert sensor connector (1) completely into patient cable connector (2). Close sensor connector cover (3) over patient cable connector (2) and lock in place as shown in Fig 5.

Disconnecting the sensor from the patient cable

1. Refer to Fig 5. While holding the sides of patient cable connector (2), lift sensor connector cover (3) from patient cable connector (2).
2. Disconnect sensor connector (1) from patient cable connector (2).

Removing the sensor from the patient

1. Open the sensor by pressing on hinge tabs. Remove the sensor from the digit and follow the cleaning instructions, if necessary.
2. Store the sensor away for its next use.

CLEANING

To clean the sensor, first remove it from the patient and disconnect it from the instrument. You may then clean the sensor by wiping it with a 70% isopropyl alcohol pad. Allow the sensor to dry prior to placement on a patient.

SPECIFICATIONS

When used with Masimo rainbow SET technology monitors or with licensed Masimo rainbow SET technology modules using rainbow Series patient cables, the DCI-mini sensors have the following performance specifications:

Sensor	DCI	DCI-P
Application Site	finger	finger or thumb
 Weight Range	> 30 kg	10-50 kg
SpO ₂ Accuracy, No Motion ¹	60-80% 3% 70-100% 2%	60-80% 3% 70-100% 2%
SpO ₂ Accuracy, Motion ²	3%	3%
Pulse Rate Accuracy, No Motion ³	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ³	5 bpm	5 bpm
SpO ₂ Accuracy, Low Perfusion ⁴	2%	2%
Pulse Rate Accuracy, Low Perfusion ³	3 bpm	3 bpm
SpCO Accuracy, No Motion ⁵	3%	3%
SpMet Accuracy, No Motion ⁶	1%	1%

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within +/- ARMS of the reference measurements in a controlled study.

¹ The Masimo rainbow SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 60%-100% SpO₂ against a laboratory co-oximeter.

² The Masimo rainbow SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

³ The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 60% to 100%.

⁴ The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁵ SpCO accuracy was determined by testing on healthy adult volunteers in the range of 1%-40% SpCO against a laboratory CO-Oximeter.

⁶ SpMet accuracy was determined by testing on healthy adult volunteers with light to dark skin pigmentation in the range of 1%-15% Methb against a laboratory CO-Oximeter.

COMPATIBILITY

These sensors are intended for use only with devices containing Masimo rainbow SET technology or pulse oximetry monitors licensed to use rainbow compatible sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free from defects in materials and workmanship for a period of one (1) year.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

Purchase or possession of this sensor confers no express or implied license to use the sensor with any device which is not separately authorized to use rainbow sensors.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
 For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult instructions for use		Separate collection for electrical and electronic equipment (WEEE).
	Follow instructions for use		Lot code
	Manufacturer		Catalogue number (model number)
	Date of manufacture YYYY-MM-DD		Masimo reference number
	Use-by YYYY-MM-DD		Pediatric patient
	Do not discard		Body weight
	Non-Sterile		Greater than
	Not made with natural rubber latex		Less than
Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Storage humidity limitation
	Mark of conformity to European Medical Device Directive 93/42/EEC		Storage temperature range
	Authorized representative in the European community		Keep dry
	Atmospheric pressure limitation		Fragile, handle with care
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.		

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , Radical-7, DCI, rainbow, SpMet, X-Cal, and SpCO are federally registered trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
MDSS GmbH