
M-LNCS™/LNCS®/LNOP® TC-I

Reusable Tip Clip Sensor








Images	2-3
en English	4-8
fr French	9-13
de German	14-18
it Italian	19-23
es Spanish	24-28
sv Swedish	29-33
nl Dutch	34-38
da Danish	39-43
pt Portuguese	44-48
zh Chinese	49-52
ja Japanese	53-57
fi Finnish	58-62
no Norwegian	63-67
cs Czech	68-72
hu Hungarian	73-77
pl Polish	78-82
ro Romanian	83-87
sk Slovak	88-92
tr Turkish	93-97
el Greek	98-102
ru Russian	103-107
ar Arabic	112-108
fa Farsi	116-113
Performance Specifications	117





M-LNCS™/LNCS®/LNOP® TC-I

Reusable Tip Clip Sensor

M-LNCS™/LNCS®/LNOP® TC-I

Reusable Tip Clip Sensor

Reusable Tip Clip Sensor

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

Prior to using this sensor, the user should read and understand the Operator's Manual for the device and this Directions for Use.

INDICATIONS

When used with Masimo SET:

The M-LNCS™/LNCS®/LNOP® TC-I, Reusable Sensor is indicated for either "spotcheck" or continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients, (weighing > 30 kg), who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

When used with NELLCOR AND NELLCOR COMPATIBLE PULSE OXIMETERS:

The M-LNCS™/LNCS®/LNOP® TC-I reusable sensor is indicated for either "spot check" or continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients, (weighing > 30 kg), in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

The TC-I reusable sensor is contraindicated for patients with pierced ears at the measuring site.

DESCRIPTION

The TC-I is for use only with devices containing Masimo SET® oximetry or licensed to use Masimo sensors and also with Nellcor and Nellcor compatible pulse oximeters. Consult individual device manufacturer for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether its devices are compatible with each sensor model.

The TC-I sensor has been verified using Masimo SET Oximetry Technology and on Nellcor's N-200 Pulse Oximeter.

If extended monitoring is required, use of a Masimo disposable sensor is recommended.

The TC-I reusable sensor must be removed and repositioned to a different monitoring site at least every four (4) hours. Because individual skin conditions and perfusion levels affect the ability of the site to tolerate sensor placement, it may be necessary to move the sensor more frequently.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- The site must be changed at least every four (4) hours.
- Circulation distal to the sensor site should be checked routinely.
- Do not use the TC-I on any site other than the ear lobe or pinna. This may result in inaccurate readings due to tissue thickness.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every two (2) hours with poorly perfused patients.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage and/or pressure necrosis or damage the sensor.
- Misapplied sensors or sensors that become partially dislodged may cause either over or under reading of actual arterial oxygen saturation.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Failure to apply the sensor properly may cause incorrect measurements.
- Elevated levels of Carboxyhemoglobin (COHb) will lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated levels of Total Bilirubin may lead to inaccurate SpO₂ measurements.
- Intravascular dyes or externally applied coloring may lead to inaccurate SpO₂ measurements.
- Sensors applied too tightly may cause erroneously low readings.

- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level.
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- Venous pulsations may cause erroneous low readings (e.g. tricuspid valve regurgitation, Trendelenberg position).
- The sensor should be free of visible defects. Never use a damaged sensor or one with exposed electrical circuitry.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- To prevent damage, do not soak or immerse the connector in any liquid solution. Do not attempt to sterilize the sensor.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reprocess or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal™ technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A) Site selection

- The preferred measuring sites are either the ear lobe or pinna.
- The site should be free of debris prior to sensor placement.

CAUTION: Before using the sensor, ensure that the sensor is physically intact, with no broken or frayed wires or damaged parts.

B) Attaching the Sensor to the Patient

1. To improve perfusion to the ear, rub the earlobe vigorously for 25-30 seconds. The ear lobe can also be rubbed with rubefacient cream (10-30% methylsalicylate and 2-10% menthol).

CAUTION: Do not use strong vasodilator creams such as nitroglycerine paste.

2. Refer to **Fig. 1**. Clip the sensor onto the ear lobe or pinna. Orient the cable so that it runs down the neck toward the body. If the TC-I sensor does not fit properly on the ear, consider using a Masimo disposable sensor or reusable finger clip sensor on another measuring site.

CAUTION: Discontinue use of the TC-I sensor if there is loss of spring tension sufficient enough to allow slippage or movement of the sensor from its proper position on the ear lobe or pinna.

CAUTION: Avoid positioning the patient so that external pressure is applied to the sensor at the measuring site.

C) Attaching the sensor to the patient cable

M-LNCS connector

1. Refer to **Fig. 2**. Orient sensor connector to patient cable connector as shown.
2. Refer to **Fig. 3**. Insert sensor connector completely into patient cable connector. Close sensor connector cover over patient cable connector until it locks in place.

LNCS connector

1. Refer to **Fig. 4**. Insert the sensor connector completely into the patient cable connector.
2. Refer to **Fig. 5**. Completely close the protective cover.

LNOP connector

1. Refer to **Fig. 6**. Orient the connecting tab so that the "shiny" contacts are facing up and mate the logo to the logo on the patient cable.
2. Insert the LNOP connector into the patient cable connector until there is a tactile or audible click of connection. Gently tug on the connectors to ensure a positive connection. Tape may be used to secure the cable to the patient for ease of movement.

D) Stress and Exercise Testing

1. Attach the TC-I sensor to the patient following the steps in B (Attaching the Sensor to the Patient).
2. Minimize unnecessary sensor motion during excessive patient movement. Using a headband, loop the ear sensor cable under the chin, and secure the cable under the headband on the side of the head opposite the ear sensor.
3. Clip the clothing clip to the patient gown.

E) Disconnecting the sensor

From the patient cable

M-LNCS connector

1. Refer to **Fig. 7**. While holding the sides of patient cable connector, lift sensor connector cover from patient cable connector.
2. Refer to **Fig. 8**. Pull firmly on the sensor connector to remove it from the patient cable.

LNCS connector

1. Refer to **Fig. 9**. Lift the protective cover to gain access to the sensor connector.
2. Refer to **Fig. 10**. Pull firmly on the sensor connector to remove from the patient cable.

LNOP connector

1. Place thumb and index finger on buttons on either side of the patient cable connector.
2. Refer to **Fig. 11**. Press firmly on the buttons and pull to remove the sensor.

From the patient

To remove the TC-I sensor from the patient's ear, gently press open the ear clip and remove.

CAUTION: To avoid discomfort to the patient or damage to the sensor, do not remove the sensor from the patient's ear by pulling on the cable.

CLEANING

CAUTION:

- Do not use undiluted bleach (5%-5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.
- Do not sterilize by irradiation, steam, autoclave, Gluteraldehyde (Cidex) or ethylene oxide

To surface clean the sensor

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the TC-I sensor by wiping it with a 70% isopropyl alcohol pad or mild detergent.
3. Allow the sensor to dry prior to placement on a patient.

or

1. If low-level disinfection is required, wipe all surfaces of the TC-I sensor and cable with a cloth or gauze pad saturated with a 1:10 bleach/water solution.
2. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the TC-I sensor and cable.
3. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

To clean or disinfect the sensor using a soaking method

1. Place the sensor in the cleaning solution (1:10 bleach/water solution), so that the sensor and desired length of cable are completely immersed.

WARNING: Do not immerse the connector end of the sensor cable as this may damage the sensor.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours. Do not immerse the connector.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable with a clean cloth or dry gauze pad.

SPECIFICATIONS

The TC-I Sensor's have the following performance specifications:

When used with:	Masimo SET technology	Nellcor technology
Sensor	TC-I	TC-I
 Body Weight	> 30 kg	> 30 kg
Application Site	Ear Lobe or Pinna	Ear Lobe or Pinna
SpO ₂ Accuracy, No Motion (70-100%) ¹	3.5%	3.5%
SpO ₂ Accuracy, Low Perfusion ²	3.5%	3.5%
Pulse Rate Accuracy, No Motion (25-240 bpm) ³	3 bpm	3 bpm
Pulse Rate Accuracy ³ , Low Perfusion	3 bpm	3 bpm

NOTE: The Arms Accuracy is calculated based upon measurement values that are statistically distributed; approximately 68% of the measured values fell within +/- the Arms value when compared to the reference device in a controlled study.

¹The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

²The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

³The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use Masimo sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE M-LNCS SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
 For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult instructions for use		Do not use if package is damaged
	Follow instructions for use		Lot code
	Manufacturer		Catalogue number (model number)
	Date of Manufacture YYYY-MM-DD		Masimo reference number
	Use By YYYY-MM-DD		Pediatric patient
	Do not discard		Body weight
	Non-sterile		Greater than
	Not made with natural rubber latex		Less than
Rx ONLY	Federal law (U.S.A) restricts this device to sale by or on the order of a physician		Storage humidity limitation
	Mark of conformity to European Medical Device Directive 93/42/EEC		Storage temperature range
	Authorized representative in the European community		Keep dry
	Separate collection for electrical and electronic equipment (WEEE).		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.

Patents: <http://www.masimo.com/patents.htm>

M-LNCS is a trademark of Masimo Corporation.

Masimo, SET, , LNCS and LNOP are federally registered trademarks of Masimo Corporation.

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© 2016 Masimo Corporation
MDSS GmbH