
DBI Series

LNCS[®] and Direct Connect

Reusable Soft Sensors






Images	2
en English	3-6
fr French	7-10
de German	11-14
it Italian	15-18
es Spanish	19-22
sv Swedish	23-26
nl Dutch	27-30
da Danish	31-34
pt Portuguese	35-38
zh Chinese	39-42
ja Japanese	43-46
fi Finnish	47-50
no Norwegian	51-54
cs Czech	55-58
hu Hungarian	59-62
pl Polish	63-66
ro Romanian	67-70
sk Slovak	71-74
tr Turkish	75-78
el Greek	79-82
ru Russian	83-86
ko Korean	87-90
ar Arabic	94-91
Performance Specifications	95

DBI Series

LNCS® and Direct Connect

Reusable Soft Sensors

DBI Series LNCS® and Direct Connect

Reusable Soft Sensors

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this sensor, the user should read and understand the Operator's Manual for the device or monitor and this Directions for Use.

INDICATIONS

When used with Masimo SET® or Masimo rainbow® SET:

DBI Series Reusable Soft Sensors are indicated for the continuous noninvasive monitoring and spot-checking of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients during no motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

When used with Nellcor and Nellcor Compatible Pulse Oximeters or Philips FAST SpO₂ Technology:

LNCS DBI Reusable Soft Sensors are indicated for the continuous noninvasive monitoring and spot-checking of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with patients during no motion conditions, in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

DBI Series Reusable Soft Sensors are contraindicated for use on active patients or for prolonged periods of use. The sensor site must be inspected at least every four (4) hours or sooner, and if the circulatory condition or skin integrity is compromised, the sensor should be applied to a different site.

DESCRIPTION

The DBI Reusable Sensor is a soft-bodied reusable finger sensor for use with systems incorporating Masimo-compatible or Masimo-approved pulse oximetry technology. Masimo SET performance specifications are met when the DBI sensor is used with Masimo patient cables and systems incorporating Masimo SET Pulse Oximetry technology. The LNCS DBI sensor may also be used on Nellcor 395 systems and below.

DBI Series sensors have been verified using Masimo rainbow SET™ technology.

The LNCS® DBI sensors have been verified on the Nellcor N-200 and N-395 Pulse Oximeter.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be checked frequently or per clinical protocol to ensure adequate circulation, skin integrity and correct optical alignment.
- Circulation distal to the sensor site should be checked routinely.
- Exercise extreme caution; skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess site as frequently as every (1) hour and move the sensor if there are signs of loss of skin integrity and/or loss of circulation or perfusion.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Carefully route sensor and patient cable to reduce the possibility of patient entanglement or strangulation.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect readings.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).
- Venous pulsations may cause inaccurate readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can affect the pulse rate displayed on the oximeter. Verify patient's pulse rate against the ECG heart rate.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment as it may result in physical harm.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- Inaccurate readings may be caused by EMI radiation interference.
- Abnormal fingers, intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate or no readings.
- High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ readings.
- Elevated levels of Methemoglobin (MetHb) may lead to inaccurate SpO₂ readings.
- Elevated levels of total bilirubin may lead to inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by severe anemia, very low arterial perfusion, or extreme motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by vasospastic disease such as Raynaud's, and peripheral vascular disease.
- Inaccurate SpO₂ readings may be caused by elevated levels of dysmethemoglobin, hypocapnic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- SpO₂ readings may be affected under very low perfusion conditions at the monitored site.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the sensor in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide as it will damage the sensor.
- Do not attempt to remanufacture, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.

- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A. Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring finger of the non-dominant hand is preferred.
- Site should be cleaned of debris prior to sensor placement.

B. Attaching the Sensor to the Patient

1. Refer to Fig. 1. Hold the sensor with star and cable on top. Open the sensor by pressing on the sides of the sensor. Insert the selected digit into the sensor.
2. Refer to Fig. 2. Orient the sensor so that the cable runs toward the top of the patient's hand. The tip of the finger should touch the raised digit stop inside the bottom of the sensor. If the fingernail is long, it may extend through the opening at the end of the sensor. Check the arrangement of the sensor to verify correct positioning. Release the sides of the sensor. Complete coverage of the detector is needed to ensure accurate data.

C. Attaching the Sensor to the Patient Cable or Device

LNCS DBI

1. Refer to Fig. 3. Properly orient sensor connector (1) and insert it completely into the patient cable connector (2). Completely close the protective cover (3).

Direct Connect to Device

1. Connect the sensor directly to a device that contains (or licensed to use) Masimo SET technology or Masimo rainbow SET technology.

NOTE: Consult device user manual for additional information.

D. Disconnecting the Sensor from the Patient Cable or Device

NOTE: To avoid damage, pull on the sensor connector, not the cable.

LNCS DBI

1. Refer to Fig. 4. Lift the protective cover (1) to gain access to the sensor connector (2). Pull firmly on the sensor connector (2) to remove from the patient cable (3).

Direct Connect to Device

1. Refer to Fig. 5. Place thumb and index finger on the buttons (1) on either side of the sensor device connector (2). Press firmly on the buttons and pull to remove the connector from the patient.

E. Removing the Sensor from the Patient:

LNCS DBI and Direct Connect to Device

1. Open the sensor by pressing on the sides and pull gently to remove the sensor from the patient.
2. Follow the cleaning instructions.

CLEANING

To surface clean the sensor:

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the sensor surface by wiping it with a 70% isopropyl alcohol pad.
3. Allow the sensor to dry prior to placement on a patient. If low-level disinfection is required, proceed to steps 4 through 7.
4. Saturate a cloth or gauze pad with one of the following solutions:
 - Glutaraldehyde (such as Metrex, Metricide 28 or Cidex 2250)
 - Ammonium Chlorides (such as ECOLAB, Huntington Brand, Aspeti-Wipe)
 - 1:10 bleach/water solution
5. Wipe all surfaces of the sensor and cable.
6. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the sensor and cable.
7. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

To clean or disinfect the sensor using a soaking method:

1. Place the sensor in one of the following cleaning solutions, so that the sensor and desired length of cable are completely immersed.
 - Glutaraldehyde (such as Metrex, Metricide 28 or Cidex 2250)
 - Ammonium Chlorides (such as ECOLAB, Huntington Brand, Aspeti-Wipe)
 - 1:10 bleach/water solution
2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours. Do not immerse the connector.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

CAUTION:

- Do not use undiluted bleach (5–5.25% sodium hypochlorite) or any cleaning solution other than those recommended herein because permanent damage to the sensor could occur.
- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.

SPECIFICATIONS

DBI Series Reusable Soft Sensors have the following specifications:

When Used with:	Masimo SET Technology		Nellcor Technology
DBI Sensor	LNCS	Direct Connect Red Direct Connect	LNCS
 Body Weight	> 30 kg	> 30 kg	> 30 kg
Application Site	Middle or ring finger	Middle or ring finger	Middle or ring finger
SpO ₂ Accuracy, No Motion	2% ¹	2% ¹	3% ⁴
SpO ₂ Accuracy, Low Perfusion	2% ²	2% ²	---
Pulse Rate Accuracy, No Motion	3 bpm ³	3 bpm ³	3 bpm ⁴
Pulse Rate Accuracy, Low Perfusion	3 bpm ²	3 bpm ²	---

NOTE: A_{RMS} accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within $\pm A_{RMS}$ of the reference measurements in a controlled study.

¹ The Masimo SET technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory CO-Oximeter.

² The Masimo SET technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

³ The Masimo SET technology has been validated for pulse rate accuracy for the range of 25–240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁴ Not validated for Motion and Low Perfusion accuracy using Nellcor Technology.

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use DBI sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

 For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE DBI SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Greater than		Storage temperature range
	Do not discard		Less than		Keep dry
	Non-Sterile		Storage humidity limitation		Do not use if package is damaged
	Not made with natural rubber latex		Atmospheric pressure limitation		Fragile, handle with care
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.				

Patents: <http://www.masimo.com/patents.htm>

rainbow SET is a trademark of Masimo Corporation.

Masimo, SET, , X-Cal and LNCS are federally registered trademarks of Masimo Corporation.

Nellcor and Oximax are federally registered trademarks of Nellcor Puritan Bennet, LLC.

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© 2018 Masimo Corporation
Philips is a trademark of Koninklijke Philips Electronics, N.V.
EU Authorized Representative for Masimo Corporation:
MDSS GmbH