M-LNCS® / LNCS® Specialty Sensors

Adult, Pediatric, Infant, and Neonatal SpO₂ Disposable Sensors







Images	2-3
en English	4-8
fr French	9-13
de German	14-18
it Italian	19-23
es Spanish	24-28
sv Swedish	29-33
nl Dutch	34-38
da Danish	39-43
pt Portuguese	44-48
zh Chinese	49-52
ja Japanese	53-57
fi Finnish	58-62
no Norwegian	63-67
cs Czech	68-72
hu Hungarian	73-77
pl Polish	78-82
ro Romanian	83-87
sk Slovak	88-92
tr Turkish	93-97
el Greek	98-102
ru Russian	103-108
ko Korean	109-112
ar Arabic	116-113
Performance Specifications	117









M-LNCS[®] / LNCS[®] Specialty Sensors

Adult, Pediatric, Infant, and Neonatal SpO₂ Disposable Sensors

M-LNCS® / LNCS® Specialty Sensors

Adult, Pediatric, Infant, and Neonatal SpO2 Disposable Sensors

M-LNCS® / LNCS® Specialty Sensors

en

Adult, Pediatric, Infant, and Neonatal SpO₂ Disposable Sensors

DIRECTIONS FOR USE

 Single patient use only

 Not made with natural rubber latex

 Non sterile

Prior to using this sensor, the user should read and understand the Operator's Manual for the Device and this Directions for Use.

INDICATIONS - When Used With Masimo Set® and Masimo compatible Pulse Oximeters:

The M-LNCS®, LNCS® Adult, Pediatric, Infant and Neonatal Specialty sensors are indicated for continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

INDICATIONS- When used with Nellcor® and Nellcor Compatible Pulse Oximeters:

The M-LNCS, LNCS Adult, Pediatric, Infant, and Neonatal Specialty sensors are indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

M-LNCS, LNCS Specialty sensors are contraindicated for patients who exhibit allergic reactions to foam rubber products and/or adhesive tape.

CAUTION: WITH SPECIALTY SENSORS, THE SENSOR OFF DETECTION PERFORMANCE MAY BE COMPROMISED.

DESCRIPTION

These Specialty sensors enable the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting. This results in a faster and more sensitive response to physiological changes.

The M-LNCS, LNCS Specialty sensors are for use with devices containing Masimo SET® oximetry (Version 4.1 or higher) or licensed to use M-LNCS, LNCS Specialty sensors and also with Nellcor and Nellcor compatible pulse oximeters, except Nellcor OxiMax® enabled devices. Consult individual device manufacturer for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether its devices are compatible with each sensor model.

When used with Nellcor and Nellcor compatible Pulse Oximeters, M-LNCS, LNCS Specialty sensors will not enable the algorithms to give new data higher priority than old data.

The M-LNCS, LNCS Specialty sensors have been verified on Masimo SET Oximetry Technology and on Nellcor's N-200 Pulse Oximeter.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- Circulation distal to the sensor site should be checked routinely.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Misapplications due to wrong sensor type can cause inaccurate or no readings.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid value regurgitation).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.

- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High levels of COHb or MethHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MethHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated Total Bilirubin levels may lead to inaccurate SpO₂ measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate SpO₂ measurements.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize the sensor.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reuse on multiple patients, reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- High oxygen concentrations may predispose a premature infant to retinopathy. Therefore, the upper alarm limit for the oxygen saturation must be carefully selected in accordance with accepted clinical standards.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 168 hours of patient monitoring time or up to 336 hours for sensors with a replaceable tape. After single-patient use, discard sensor.

INSTRUCTIONS:

A) Site Selection

- Always choose a site that is well perfused and will completely cover the sensor's detector window.
- Site should be cleaned of debris and dry prior to sensor placement.

Trauma Sensor (Adult)

> 30 kg: The preferred site is the middle or ring finger of the non-dominant hand.

Newborn Sensor (Infant/Pediatric)

3–10 kg: The preferred site is the great toe. Alternatively, the toe next to the great toe, or the thumb can be used.

10–30 kg: The preferred site is the middle or ring finger of the non-dominant hand.

Newborn Sensor (Neonatal)

< 3 kg: The preferred site is the foot. Alternatively, across the palm and back of the hand can be used.

B) Attaching the sensor to the patient

1. Open the pouch and remove the sensor. Remove the backing from the sensor, if present.

ADULTS (> 30 kg), PEDIATRIC (10–30 kg)

2. Refer to **Fig. 1a**. Direct the sensor cable so that it runs along the top of the hand or points away from the patient. Position the detector onto the fleshy part of the finger. Alternatively, the sensor may also be applied to the toe (not shown). Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to **Fig. 1b**. Wrap the foam wrap around the finger and ensure that the emitter window (red star) aligns on the top of the finger directly opposite of the detector.
4. Refer to **Fig. 1c**. Check the sensor to verify correct positioning and reposition if necessary.

INFANTS (3–10 kg)

2. Refer to **Fig. 2a**. Direct the sensor cable so that it runs along the bottom of the foot or points away from the patient. Position the detector onto the fleshy part of the great toe. Alternatively, the sensor may also be applied to the thumb (not shown). Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to **Fig. 2b**. Wrap the foam wrap around the toe and ensure that the emitter window (red star) aligns on the top of the toe directly opposite of the detector.
4. Refer to **Fig. 2c**. Check the sensor to verify correct positioning and reposition if necessary.

NEONATES (< 3 kg)

2. Refer to **Fig. 3a**. For fragile skin, the stickiness of the medical grade adhesive can be diminished or eliminated by daubing the adhesive areas with a cotton ball or with gauze.
3. Refer to **Fig. 3b**. Direct the sensor cable so that it either points away from the patient or runs along the bottom of the foot. Apply the detector onto the fleshy part of the lateral aspect sole of the foot aligned with the fourth toe. Alternatively, the detector may also be applied to the top of the foot (not shown). Complete coverage of the detector window is needed to ensure accurate data.

4. Refer to **Fig. 3c**. Wrap the foam wrap around the foot and ensure that the emitter window (red star) aligns directly opposite of the detector. Be careful to maintain proper alignment of the detector and emitter windows while attaching the foam wrap to secure the sensor.
5. Refer to **Fig. 3d**. Check the sensor to verify correct positioning and reposition if necessary.

C) Attaching the Sensor to the Patient Cable

M-LNCS

Refer to **Fig. 4**. Insert the sensor connector completely into the patient cable connector and lock into place.

LNCS

Refer to **Fig. 5**. Insert the sensor connector completely into the patient cable connector (1). Completely close the protective cover (2).

D) Reattachment

- The sensor may be reapplied to the same patient if the emitter and detector windows are clear and the adhesive still adheres to the skin.
 - If the adhesive no longer adheres to the skin, use a new sensor.
- NOTE:** When changing application sites, or reattaching the sensor, first disconnect sensor from the patient cable.

E) Disconnecting the Sensor from the Patient Cable

M-LNCS

Refer to **Fig. 7**. Pull firmly on the sensor connector to remove it from the patient cable.

LNCS

Refer to **Fig. 8**. Lift the protective cover to gain access to the sensor connector (1). Pull firmly on the sensor connector to remove from the patient cable (2).

CAUTION: TO PREVENT DAMAGE, DO NOT SOAK OR IMMERSE THE SENSOR IN ANY LIQUID SOLUTION. DO NOT STERILIZE BY IRRADIATION, STEAM, AUTOCLAVE OR ETHYLENE OXIDE.

SPECIFICATIONS

When Used With Masimo Set and Masimo compatible Pulse Oximeters:				
M-LNCS/LNCS Specialty Sensor	Trauma (Adult)	Newborn (Infant/Pediatric)		Newborn (Neonatal) ⁵
 Body Weight	> 30 kg	3–10 kg	10–30 kg	< 3 kg
Application Site	Finger or toe	Thumb or great toe	Finger or toe	Hand or foot
SpO ₂ Accuracy, No Motion ¹	2%	2% ⁵	2%	3%
SpO ₂ Accuracy, Motion ²	3%	3%	3%	3%
SpO ₂ Accuracy, Low Perfusion ³	2%	2%	2%	3%
Pulse Rate Accuracy, No Motion ⁴	3 bpm	3 bpm	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ⁴	5 bpm	5 bpm	5 bpm	5 bpm
Pulse Rate Accuracy, Low Perfusion ⁴	3 bpm	3 bpm	3 bpm	3 bpm
When used with Nellcor and Nellcor Compatible Pulse Oximeters: ⁶				
SpO ₂ Accuracy, No Motion	2%	2%	2%	3%
Pulse Rate Accuracy, No Motion	3 bpm	3 bpm	3 bpm	3 bpm

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within \pm ARMS of the reference measurements in a controlled study.

¹ The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory co-oximeter.

² The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory co-oximeter.

³ The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁴ The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25–240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁵ The saturation accuracy of the Neonate and Preterm sensors were validated on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

⁶ Not validated for Motion and Low Perfusion accuracy using Nellcor Technology.

COMPATIBILITY

This sensor is intended for use only with devices containing Masimo SET oximetry (Version 4.1 or higher) or pulse oximetry monitors licensed to use M-LNCS, LNCS sensors and also with Nellcor and Nellcor compatible pulse oximeters. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

THIS SINGLE-PATIENT SENSOR IS LICENSED TO YOU UNDER THE PATENTS OWNED BY MASIMO FOR SINGLE-PATIENT USE ONLY. BY ACCEPTANCE OR USE OF THIS PRODUCT, YOU ACKNOWLEDGE AND AGREE THAT NO LICENSE IS GRANTED FOR USE OF THIS PRODUCT WITH MORE THAN A SINGLE PATIENT.

AFTER SINGLE-PATIENT USE, DISCARD SENSOR. PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE M-LNCS, LNCS SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
 For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events. The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Greater than		Storage temperature range
	Do not re-use/Single patient use only		Less than		Keep dry
	Non-Sterile		Storage humidity limitation		Do not use if package is damaged
	Not made with natural rubber latex		Fragile, handle with care		Atmospheric pressure limitation
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available for CE mark countries.				

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , M-LNCS, X-Cal, and LNCS are federally registered trademarks of Masimo Corporation.

Nellcor and OxiMax are federally registered trademark of Nellcor Puritan Bennett Incorporated.

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© 2018 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH