# NomoLine®

## LH & HH Series

Single-patient use sampling lines





Images	2-3
<b>en</b> English	4-6
<b>fr</b> French	7-9
<b>de</b> German	10-12
<b>it</b> Italian	13-15
<b>es</b> Spanish	16-18
<b>sv</b> Swedish	19-21
<b>nl</b> Dutch	22-24
<b>da</b> Danish	25-27
<b>pt</b> Portuguese	28-30
<b>zh</b> Chinese	31-33
<b>ja</b> Japanese	34-36
<b>fi</b> Finnish	37-39
<b>no</b> Norwegian	40-42
<b>cs</b> Czech	43-45
<b>hu</b> Hungarian	46-48
<b>pl</b> Polish	49-51
<b>ro</b> Romanian	52-54
<b>sk</b> Slovak	55-57
<b>tr</b> Turkish	58-60
<b>el</b> Greek	61-63
<b>ru</b> Russian	64-66
<b>ko</b> Korean	67-69
<b>et</b> Estonian	70-72
<b>lt</b> Lithuanian	73-75
<b>lv</b> Latvian	76-78
<b>ar</b> Arabic	80-79
<b>fa</b> Farsi	83-81

# NomoLine®

## LH & HH Series

Single-patient use sampling lines





	-40°C – +70°C
	5% – 100%
	+200hPa – +1200hPa

# NomoLine®

## LH & HH Series

Single-patient use sampling lines

## Single-patient use sampling lines

### DIRECTIONS FOR USE

② Single Patient Use Only

☒ Not made with natural rubber latex

#### INDICATIONS

The NomoLine® Product Family is indicated for the measurement of respiratory rate and respiratory and anesthetic gases in adult, pediatric, infant and neonatal patients. The NomoLine Product Family includes single-patient use and multi-use devices for gas sampling and/or oxygen delivery. The NomoLine Product Family is indicated for use by clinical professionals in healthcare environments, including mobile environments.

#### WARNINGS

- Carefully route the NomoLine Sampling Line to reduce the risk of patient entanglement or strangulation.
- Do not apply negative pressure to remove condensed water.
- Dispose in accordance with local regulations for biohazardous waste.
- Do not submerge the NomoLine Sampling Line in any cleaning solution or attempt to sterilize by autoclave, irradiation, steam, gas, ethylene oxide or any other method as this may seriously damage the NomoLine Sampling Line.
- Do not re-use due to cross contamination risk.
- Always replace the NomoLine Sampling Line if an occlusion alarm occurs.
- Do not use NomoLine Cannulas on patients with obstructions in one or both nostrils.
- The CO<sub>2</sub> measurement can be diminished if NomoLine Cannulas are used during mouth breathing or on a patient that has a perforated septum.
- Do not use NomoLine Adult/Pediatric Airway Adapter Sets for infants/neonates as the adapter adds 6 ml dead space.
- Do not use NomoLine Infant/Neonatal Airway Adapter Sets for adults/pediatrics as this may cause excessive flow resistance (0.7 ml dead space).
- Do not use NomoLine Airway Adapter Sets or Cannulas with inhalers or nebulized medications as this may clog the bacteria filter.

#### DESCRIPTION

NomoLine Sampling Lines include nasal and nasal/oral cannulas for non-intubated patients and airway adapter sets for intubated patients. The NomoLine Sampling Lines are accessories to the Masimo ISA sidestream gas analyzers. Refer to the ISA User's Guide for technical specification and full system use instructions for the ISA product system.

#### INSTRUCTIONS

- Open the bag and remove the NomoLine Sampling Line. For convenient preparations, do not yet remove the tape that keeps the tubing coiled.

##### A. Fitting a NomoLine Cannula to the patient

- If fitting a NomoLine Nasal CO<sub>2</sub> Cannula refer to **Fig. A1**: Insert the prongs into the nostrils.
- If fitting a NomoLine Nasal/Oral CO<sub>2</sub> Cannula refer to **Fig. A2**: Insert the prongs into the nostrils and position the oral scoop in front of the mouth.
- If fitting a NomoLine Single Nasal Prong CO<sub>2</sub> Cannula refer to **Fig. A3**: Insert the prong in one of the nostrils.
- Refer to **Fig. A4**: Pass the Cannula lines over the ears (1) and adjust the slider for comfort and fit under the chin (2). (This step does not apply to NomoLine Single Nasal Prong CO<sub>2</sub> Cannula.)

##### B. Connecting a NomoLine Airway Adapter Set to the patient circuit

- Refer to **Fig. B1**: Connect the small diameter end of the Airway Adapter to the Y-piece of the ventilation source.
- Refer to **Fig. B2**: Connect the large diameter end of the Airway Adapter to the patient's endotracheal tube.

##### C. Connecting a NomoLine Sampling Line to the ISA gas analyzer

- Refer to **Fig. C**: Remove the tape that keeps the tubing coiled. Connect the NomoLine Sampling Line to the ISA gas inlet connector. A steady green light indicates that the NomoLine Sampling Line is correctly connected to ISA.
- Ensure that the tubing is properly connected and that it is not twisted or crimped.

##### D. Connecting a NomoLine Cannula to an oxygen source

- Refer to **Fig. D**: Remove the tape that keeps the tubing coiled. Attach the oxygen supply tubing to the oxygen source. Adjust the oxygen flow. (This step applies only to NomoLine Cannulas with O<sub>2</sub>.)
- Ensure that the tubing is properly connected and that it is not twisted or crimped.

#### REPLACEMENT

NomoLine Sampling Lines are single patient use products and should after use be disposed of in accordance with local regulations for biohazardous waste and replaced for each new patient with a maximum of 14 days of use\*. They should also be replaced if the sampling line becomes occluded. Occlusion is indicated by a red flashing light on the ISA gas inlet connector. Replace the NomoLine Sampling Line and wait until the gas inlet connector switches to a steady green light.

\* The lifetime of the product depends on the application and product type.

- Remove and discard the disposable, single-patient use NomoLine Sampling Line according to local regulations for biohazardous waste.

- Connect the new Sampling Line to the ISA gas inlet connector.

---

## **CLEANING**

NomoLine Sampling Lines are single-patient use products and should after use be disposed of in accordance with local regulations for biohazardous waste.

---

## **DISPOSAL**

The NomoLine Sampling Line may after use contain accumulated moisture and fluids. Dispose the used product in accordance with local regulations for biohazardous waste.

---

## **WARRANTY**

MASIMO WARRANTS TO THE INITIAL BUYER ONLY THAT THESE PRODUCTS, WHEN USED IN ACCORDANCE WITH THE DIRECTIONS PROVIDED WITH THE PRODUCTS BY MASIMO, WILL BE FREE OF DEFECTS IN MATERIALS AND WORKMANSHIP FOR A PERIOD OF SIX (6) MONTHS. SINGLE USE PRODUCTS ARE WARRANTED FOR SINGLE PATIENT USE ONLY. THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

---

## **WARRANTY EXCLUSIONS**

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to products that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

**The following symbols may appear on the product or product labeling:**

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Gas sample	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code	<b>EC REP</b>	Authorized representative in the European community
	Manufacturer		Catalogue number (model number)		Body weight
	Use-by YYYY-MM-DD		Masimo reference number		Keep dry
	Do not re-use/Single patient use only		Less than or equal to		Do not use if package is damaged and consult instructions for use
	Not made with natural rubber latex		Storage temperature range	<b>LH</b> (orange background)	Collects water. Short term applications with low humidity
	Caution		Storage humidity limitation	<b>HH</b> (yellow background)	Removes aspired and condensed water. Long term applications with high humidity
	European Union Conformity Mark		Atmospheric pressure limitation		Single patient - multiple use
	Medical device		Unique device identifier		Importer
	Distributor		Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eFU is not available in all countries.</b>		

Patents: <http://www.masimo.com/patents.htm>

Masimo, NomoLine and  are federally registered trademarks of Masimo Corporation.

ISA is a trademark of Masimo Corporation.

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© 2021 Masimo Corporation
for Masimo Corporation:
MDSS GmbH