RD SET™ DCI® & DCI-P Series

Adult & Pediatric Reusable Finger Clip Sensors








Images	2
en English	3-5
fr French	6-8
de German	9-11
it Italian	12-14
es Spanish	15-17
sv Swedish	18-20
nl Dutch	21-23
da Danish	24-26
pt Portuguese	27-29
zh Chinese	30-32
ja Japanese	33-35
fi Finnish	36-38
no Norwegian	39-41
cs Czech	42-44
hu Hungarian	45-47
pl Polish	48-50
ro Romanian	51-53
sk Slovak	54-56
tr Turkish	57-59
el Greek	60-62
ru Russian	63-65
ar Arabic	68-66
fa Farsi	71-69
Performance Specs	72

RD SET™ DCI® & DCI-P Series

Adult & Pediatric Reusable Finger Clip Sensors

RD SET™ DCI® & DCI-P Series

Adult & Pediatric Reusable Finger Clip Sensors

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

INDICATIONS

The RD SET DCI® and DCI-P reusable sensors are indicated for either “spot check” or continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

The RD SET DCI and DCI-P sensors are contraindicated for use on mobile patients or for prolonged periods of use. They are not intended for long-term monitoring. They must be removed and repositioned to a different monitoring site at least every four (4) hours. Because individual skin conditions and perfusion levels affect the ability of the site to tolerate sensor placement, it may be necessary to move the sensor more frequently. If extended monitoring is required, use of an RD SET disposable sensor is recommended.

DESCRIPTION

The RD SET DCI and DCI-P reusable sensors are for use only with devices containing Masimo SET® oximetry or licensed to use RD SET DCI and DCI-P sensors. Consult individual device manufacturer for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether its device is compatible with each sensor model.

The RD SET DCI and DCI-P sensors have been verified with Masimo SET Oximetry Technology.

WARNING: Masimo sensors and cables are designed for use with instruments containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be checked frequently or per clinical protocol to ensure adequate circulation, skin integrity and correct optical alignment.
- Exercise extreme caution; skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess site as frequently as every (1) hour and move the sensor if there are signs of loss of skin integrity and/or loss of circulation or perfusion.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Carefully route sensor and patient cable to reduce the possibility of patient entanglement or strangulation.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect readings.
- Misapplications due to wrong sensor types can cause inaccurate or no readings.
- Inaccurate SpO₂ readings may be caused by abnormal venous pulsation or venous congestion.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can affect the pulse rate displayed on the oximeter. Verify patient's pulse rate against the ECG heart rate.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment as it may result in physical harm.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- Inaccurate readings may be caused by EMI radiation interference.
- Abnormal fingers, intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate or no readings.
- High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ readings.
- Elevated levels of Methemoglobin (MetHb) may lead to inaccurate SpO₂ readings.
- Elevated levels of total bilirubin may lead to inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by severe anemia, very low arterial perfusion, or extreme motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by vasospastic disease such as Raynaud's, and peripheral vascular disease.
- Inaccurate SpO₂ readings may be caused by elevated levels of dysmethemoglobin, hypoxic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- SpO₂ readings may be affected under very low perfusion conditions at the monitored site.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the sensor in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide as it will damage the sensor.
- Do not attempt to reprocess, precondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Caution:** Replace the sensor when a replace sensor or a persistent low SiQ message is displayed when the patient monitoring time is exhausted.
- Note:** The sensor is provided with X-Cal™ technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A. Site Selection

- Choose a site that is well perfused and least restricts a conscious patient's movements. The ring finger of the non-dominant hand is preferred. Alternatively, the other digits on the non-dominant hand may be used.
- Always choose a site that will completely cover the sensor's detector window. The great toe or long toe (next to the great toe) may be used on restrained patients or patients whose hands are unavailable.
- Site should be cleaned of debris prior to sensor placement.

CAUTION: Before using the sensor, ensure that the sensor is physically intact, with no broken or frayed wires or damaged parts.

B. Attaching the sensor to the patient

1. Refer to Fig. 1a. Open the sensor by pressing on hinge tabs. Place the selected digit over the sensor window of the RD SET DCI or DCI-P. The fleshiest part of the digit should be covering the detector window in the lower half of the sensor. The top half of the sensor is identified by the cable. On a finger site, the tip of the finger should touch the raised digit stop inside the sensor. If the fingernail is long, it may extend over and pass the finger stop.
2. Refer to Fig. 1b. The hinged tabs of the sensor should open to evenly distribute the grip of the sensor along the length of the finger. Check position of sensor to verify correct positioning. Complete coverage of the detector window is needed to ensure accurate data.
3. Refer to Fig. 1c. Orient the sensor so that the cable will be running towards the top of the patient's hand.

NOTE: With smaller digits, in order to completely cover the detector window, it may be necessary to use a lower weight range sensor. The sensor is not intended for use on the thumb or across a child's hand or foot.

C. Attaching the Sensor to the Patient Cable

1. Refer to Fig. 2a. Properly orient the sensor connector and insert the sensor connector completely into the patient cable connector.
2. Refer to Fig. 2b. Close the protective latch completely.

D. Disconnecting the Sensor from the Patient Cable

1. Refer to Fig. 3a. Lift up the protective latch.
2. Refer to Fig. 3b. Pull firmly on the sensor connector to remove it from the patient cable.

NOTE: To avoid damage, pull on the sensor connector, not the cable.

CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the sensor by wiping with a 70% isopropyl alcohol pad.
3. Allow the sensor to dry thoroughly prior to placement on a patient.

or

1. If low-level disinfection is required, use a 1:10 bleach / water solution.
2. Saturate a cloth or gauze pad with the cleaning solution and wipe all surfaces of the sensor and cable.
3. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the sensor and cable.
4. Dry the sensor and cable with a clean cloth or dry gauze pad.

CAUTIONS:

- Do not use undiluted bleach (5% - 5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor may occur.
- Do not immerse the sensor or connector in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave or ethylene oxide.

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors or with licensed Masimo SET pulse oximetry modules RD SET DCI[®] and DCI-P reusable sensors have the following specifications:

Used with:	Masimo SET Technology	
	RD SET DCI	RD SET DCI-P
Sensor		
 Body Weight	> 30 kg	10 - 50 kg
Application site	Finger	Finger or toe
SpO ₂ Accuracy, No Motion, (70 - 100%) ¹	2%	2%
SpO ₂ Accuracy, Motion ²	3%	3%
SpO ₂ Accuracy, Low Perfusion ³	2%	2%
Pulse Rate ⁴ Accuracy, No Motion, (25 - 240 bpm)	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ⁴	5 bpm	5 bpm
Pulse Rate Accuracy, Low Perfusion ⁴	3 bpm	3 bpm

NOTE: The Arms Accuracy is calculated based upon measurement values that are statistically distributed; approximately 68% of the measured values fell within +/- the Arms value when compared to the reference device under a controlled study.

¹The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

²The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

³The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

⁴The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use RD SET sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

 For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.


NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RD SET SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult instructions for use		Authorized representative in the European community	>	Greater than
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	<	Less than
	Manufacturer		Light Emitting Diode (LED) LED emits light when current flows through		Storage humidity limitation
	Date of manufacture YYYY-MM-DD		Do not discard		Storage temperature range
	Non-sterile		Lot code		Do not use if package is damaged
	Not made with natural rubber latex		Catalogue number (model number)		Atmospheric pressure limitation
Rx ONLY	Caution: Federal law (U.S.A) restricts this device to sale by or on the order of a physician		Masimo reference number		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.
	Mark of conformity to European medical device directive 93/42/EEC		Body weight		

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, DCI and  are federally registered trademarks of Masimo Corporation.

RD SET and X-Cal are trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH