
RD SET™ TC-I






Reusable Tip Clip Sensor

Images	2
en English	3-6
fr French	7-11
de German	12-16
it Italian	17-20
es Spanish	21-24
sv Swedish	25-28
nl Dutch	29-33
da Danish	34-38
pt Portuguese	39-43
zh Chinese	44-47
ja Japanese	48-52
fi Finnish	53-56
no Norwegian	57-60
cs Czech	61-64
hu Hungarian	65-69
pl Polish	70-74
ro Romanian	75-79
sk Slovak	80-83
tr Turkish	84-88
el Greek	89-93
ru Russian	94-98
ko Korean	99-102
ar Arabic	106-103
Performance Specifications	107

RD SET™ TC-I

Reusable Tip Clip Sensor

Reusable Tip Clip Sensor

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

Prior to using this Sensor the user should read and understand the Operator's Manual for the Device or Monitor and this Directions for Use.

INDICATIONS

The RD SET™ TC-I reusable sensor is indicated for either "spotcheck" or continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult and pediatric patients, who are well or poorly perfused in hospitals, hospital- type facilities, mobile, and home environments.

CONTRAINDICATIONS

The RD SET TC-I reusable sensor is contraindicated for patients with pierced ears at the measuring site.

The RD SET TC-I reusable sensor must be removed and repositioned to a different monitoring site at least every four (4) hours. Because individual skin conditions and perfusion levels affect the ability of the site to tolerate sensor placement, it may be necessary to move the sensor more frequently.

DESCRIPTION

The RD SET TC-I reusable sensor is for use only with devices containing Masimo® SET® oximetry or licensed to use RD SET TC-I sensors and compatible pulse oximeters. Consult individual device manufacturer for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether its device is compatible with each sensor model.

The RD TC-I sensors have been verified using Masimo SET Oximetry Technology.

WARNING: Masimo sensors and cables are designed for use with instruments containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- Do not use the RD SET TC-I on any site other than the ear lobe or pinna. This may result in inaccurate readings due to tissue thickness.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise extreme caution; skin erosion, tissue ischemia, and/or pressure necrosis can be caused when the sensor is not frequently moved, applied too tightly, or becomes too tight due to edema. Assess site as frequently as every (1) hour and move the sensor if there are signs of loss of skin integrity and/or loss of circulation or perfusion.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Carefully route sensor and patient cable to reduce the possibility of patient entanglement or strangulation.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Misapplications due to wrong sensor types can cause inaccurate or no readings.
- Inaccurate SpO₂ readings may be caused by abnormal venous pulsation or venous congestion.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor, Trendelenburg position).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment as it may result in physical harm.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.

- Inaccurate readings may be caused by EMI radiation interference.
- Intravascular dyes such as indocyanine green or methylene blue may lead to inaccurate or no readings.
- High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated levels of total bilirubin may lead to inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by severe anemia, very low arterial perfusion, or extreme motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by vasospastic disease such as Raynaud's, and peripheral vascular disease.
- Inaccurate SpO₂ readings may be caused by elevated levels of dyshemoglobin, hypocapnic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- SpO₂ readings may be affected under very low perfusion conditions at the monitored site.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the connector in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide as it will damage the sensor.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor or a persistent low SIQ message is displayed when the patient monitoring time is exhausted.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A. Site Selection

- The preferred measuring sites are either the ear lobe or pinna.
- The site should be free of debris prior to sensor placement.

CAUTION: If extended monitoring is required, use of a RD SET disposable sensor is recommended.

CAUTION: Before using the sensor, ensure that the sensor is physically intact, with no broken or frayed wires or damaged parts.

B. Attaching the sensor to the patient

1. To improve perfusion to the ear, rub the earlobe vigorously for 25-30 seconds. The ear lobe can also be rubbed with rubefacient cream (10-30% methylsalicylate and 2-10% menthol).

CAUTION: Do not use strong vasodilator creams such as nitroglycerine paste.

2. Refer to Fig. 1. Clip the sensor onto the ear lobe or pinna. Orient the cable so that it runs down the neck toward the body. If the RD SET TC-I sensor does not fit properly on the ear, consider using an RD SET disposable sensor or RD SET reusable finger clip on another measuring site.

CAUTION: Discontinue use of the ear sensor if there is loss of spring tension sufficient enough to allow slippage or movement of the sensor from its proper position on the ear lobe or pinna.

CAUTION: Avoid positioning the patient so that external pressure is applied to the sensor at the measuring site.

C. Attaching the Sensor to the Patient Cable

1. Refer to Fig. 2. Properly orient the sensor connector and insert the sensor connector completely into the patient cable connector.
2. Refer to Fig. 3. Close the protective latch cover completely.

D. Stress and Exercise Testing

1. Attach the RD SET TC-I sensor to the patient following the steps in B (Attaching the Sensor to the Patient).
2. Minimize unnecessary sensor motion during excessive patient movement. Using a headband, loop the ear sensor cable under the chin, and secure the cable under the headband on the side of the head opposite the ear sensor.
3. Clip the clothing clip to the patient gown.

E. Disconnecting the Sensor

From the patient cable

1. Refer to Fig. 4. Lift up the protective latch.
2. Refer to Fig. 5. Pull firmly on the sensor connector to remove it from the patient cable.

NOTE: To avoid damage, pull on the sensor connector, not the cable.

From the patient

1. To remove the RD SET TC-I from the patient's ear, gently press open the ear clip and remove.

CAUTION:

To avoid discomfort to the patient or damage to the sensor, do not remove the sensor from the patient's ear by pulling on the cable.

CLEANING

To surface clean the sensor

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the RD SET TC-I sensor by wiping it with a 70% isopropyl alcohol pad or mild detergent.
3. Allow the sensor to dry prior to placement on a patient.

or

1. If low-level disinfection is required, wipe all surfaces of the RD SET TC-I sensor and cable with a cloth or gauze pad saturated with a 1:10 bleach/water solution.
2. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the RD SET TC-I sensor and cable.
3. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

To clean or disinfect the sensor using a soaking method

1. Place the sensor in the cleaning solution (1:10 bleach/water solution), so that the sensor and desired length of cable are completely immersed.

WARNING: Do not immerse the connector end of the sensor cable as this may damage the sensor.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours. Do not immerse the connector.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable with a clean cloth or dry gauze pad.

CAUTION:

- Do not use undiluted bleach (5%-5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.
- Do not immerse the connector on the RD SET TC-I cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, Gluteraldehyde (Cidex) or ethylene oxide.

ENVIRONMENTAL SPECIFICATIONS

Storage/Transport Temperature	-40°C to +70°C, ambient humidity
Operating Temperature	+5°C to +40°C, ambient humidity
Storage/Transport Humidity	10% to 95%, non-condensing
Operating Humidity	10% to 95%, non-condensing

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, or with licensed Masimo SET pulse oximetry modules the RD SET TC-I reusable sensor has the following specifications:

When used with compatible pulse oximeters the RD SET TC-I reusable sensor has the following specifications:

Sensor	RD SET TC-I
Used with:	Masimo SET Technology
 Body Weight	> 30 kg
Application Site	Ear Lobe or Pinna
SpO ₂ Accuracy, No Motion (70-100%) ¹	± 3.5%
SpO ₂ Accuracy, Low Perfusion ²	± 3.5%
Pulse Rate ³ Accuracy, No Motion (25-240 bpm)	± 3 bpm
Pulse Rate Accuracy, Low Perfusion	± 3 bpm

NOTE: The ARMS Accuracy is calculated based upon measurement values that are statistically distributed; approximately 68% of the measured values fell within ± the ARMS value when compared to the reference device under a controlled study.

¹The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory co-oximeter.

²The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

³The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

COMPATIBILITY

This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use RD SET TC-I sensors and compatible pulse oximeters. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.


NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RD SET TC-I SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult Instructions for Use		Mark of Conformity to European Medical Device Directive 93/42/EEC		Do not discard
	Follow Instructions for Use		Authorized representative in the European community		Lot code
	Manufacturer		Separate collection for electrical and electronic equipment (WEEE).		Catalogue number (model number)
	Date of Manufacture YYYY-MM-DD		Masimo reference number		Storage Humidity Limitation
	Non-sterile		Body weight		Storage Temperature Range
	Not made with natural rubber latex		Do not use if package is damaged		Instructions/Directions for Use/ Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available in all countries.
Rx ONLY	Federal law (U.S.A) restricts this device to sale by or on the order of a physician		Atmospheric Pressure Limitation		

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, X-Cal, and  are federally registered trademarks of Masimo Corporation.

RD SET is a trademark of Masimo Corporation.

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© 2019 Masimo Corporation
MDSS GmbH