# RD SET™ TF-I®

## SpO<sub>2</sub> Reusable Transflectance Forehead Sensor










Images	2
<b>en</b> English	3-6
<b>fr</b> French	7-11
<b>de</b> German	12-16
<b>it</b> Italian	17-21
<b>es</b> Spanish	22-26
<b>sv</b> Swedish	27-30
<b>nl</b> Dutch	31-35
<b>da</b> Danish	36-39
<b>pt</b> Portuguese	40-44
<b>zh</b> Chinese	45-48
<b>ja</b> Japanese	49-52
<b>fi</b> Finnish	53-56
<b>no</b> Norwegian	57-60
<b>cs</b> Czech	61-64
<b>hu</b> Hungarian	65-69
<b>pl</b> Polish	70-74
<b>ro</b> Romanian	75-79
<b>sk</b> Slovak	80-84
<b>tr</b> Turkish	85-88
<b>el</b> Greek	89-93
<b>ru</b> Russian	94-98
<b>ko</b> Korean	99-102
<b>ar</b> Arabic	106-103
Performance Specifications	107

# RD SET™ TF-I®

## SpO2 Reusable Transflectance Forehead Sensor

### DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

**Prior to using this sensor, the user should read and understand the Operator's Manual for the monitor and this Directions for Use.**

#### INDICATIONS

The RD SET™ TF-I® Reusable Transflectance forehead sensor is indicated for the continuous noninvasive monitoring of arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate for adult and pediatric patients weighing > 30 kg.

#### CONTRAINDICATIONS

The RD SET TF-I is contraindicated for certain patient positions that affect the monitoring site - see the Warnings Section below. The RD SET TF-I must be removed and repositioned to a different monitoring site at least every two (2) hours. If extended monitoring is required, use of an RD SET single patient adhesive digit sensor is recommended.

#### DESCRIPTION

The RD SET TF-I sensor is applied to the sensor site using a headband.

The RD SET TF-I is for use only with devices containing Masimo SET® oximetry or licensed to use TF-I sensors. Consult individual oximetry system manufacturers for compatibility of particular devices and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor model.

The RD-SET TF-I series has been verified using Masimo SET® Oximetry Technology.

**WARNING:** Masimo® sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

#### WARNINGS, CAUTIONS, AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use immediately. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be observed and the sensor repositioned at least every two (2) hours to ensure proper skin integrity and proper alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- SpO<sub>2</sub> readings may be affected under very low perfusion conditions at the monitored site.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Sensors applied too tightly may cause erroneously low readings. Applying the headband too tightly can lead to inaccurate saturation measurements, or possibly lead to temporary pressure marks from the sensor.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Do not use this sensor under conditions where the patient is in the Trendelenburg position (head lower than the heart).
- Do not use the sensor during MRI scanning or in a MRI environment.
- Inaccurate readings may be caused by EMI radiation interference.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.
- Venous pulsations may cause inaccurate readings (e.g. tricuspid valve regurgitation, Trendelenburg position). Mechanical ventilation can cause increased intra-thoracic pressure and increase venous pooling in the head. This may cause inaccurate readings.
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. Trendelenburg position).

- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO<sub>2</sub> measurements.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO<sub>2</sub> measurements.
- High levels of COHb or MetHb may occur with a seemingly normal SpO<sub>2</sub>. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated Total Bilirubin levels may lead to inaccurate SpO<sub>2</sub> measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring may lead to inaccurate readings.
- Inaccurate SpO<sub>2</sub> readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO<sub>2</sub> readings.
- Inaccurate SpO<sub>2</sub> readings may be caused by elevated levels of dyshemoglobin, hypocapnic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- To prevent interference from ambient light, ensure that the sensor is properly applied. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the connector in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide.
- Do not attempt to remanufacture, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

---

## INSTRUCTIONS

### A. Site Selection

- The preferred measuring site is the forehead just above either eye. Do not place the RD SET TF-I on sites with a palpable pulse.
- Site should be cleaned of debris and dry prior to sensor placement.

**CAUTION:** Before using the sensor, ensure the sensor is physically intact, with no broken or frayed wires or damaged parts.

### B. Attaching the sensor to the patient

1. Ensure sensor surface contacting patient skin is dry, clean, free of debris and oil.
2. OPTIONAL Refer to Fig. 1. For added sensor stability, place one adhesive pad on sensor surface to contact patient skin, being careful to align the pad such that it does not cover sensor lenses. Once pad is securely placed, remove top release liner by pulling from release liner tab.

**WARNING:** Do not use adhesive pads if the patient exhibits allergic reactions to adhesive tape.

3. Refer to Fig. 2. Orient the sensor with the emitter and detector 'bumps' toward the patient. Orient the sensor so the cable runs upwards and away from the patient. Gently press on the sensor so that adhesive pad forms a good contact with patient skin.

### C. Applying the Headband

**CAUTION:** Headband must be used to avoid inaccurate readings.

1. Refer to Fig. 3. Wrap the Masimo Headband around the patient's head loosely enough to avoid restricting circulation around the site and tightly enough to maintain proper alignment of the sensor.
2. Refer to Fig. 4. Secure the Headband using the "hook and loop" tab.
3. Refer to Fig. 5. Verify that the sensor is covered and secure behind the Masimo Headband.
4. Use the clothing clip to secure the cable to the patient's gown (not shown).
5. Carefully route the patient cable to reduce the possibility of patient entanglement or strangulation.

### D. Attaching the sensor to the patient cable

1. Refer to Fig. 6. Properly orient the sensor connector and insert the sensor connector completely into the patient cable connector.
2. Refer to Fig. 7. Close the protective latch cover completely.

## E. Disconnecting the sensor from the patient cable

1. Refer to Fig. 8. Lift up the protective cover.
2. Refer to Fig. 9. Pull firmly on the sensor connector to remove it from the patient cable.

**NOTE:** To avoid damage, pull on the sensor connector, not the cable.

## CLEANING

### CAUTION:

- Do not use undiluted bleach (5–5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.
- To prevent damage, do not soak or immerse the connector in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, Gluteraldehyde (Cidex) or ethylene oxide.

### To surface clean the sensor

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Remove the adhesive pad by pulling on white tab and discard. Any residue left by adhesive can easily be removed and cleaned by wiping with 70% isopropyl alcohol.
3. Clean the RD SET TF-I sensor by wiping it with a 70% isopropyl alcohol pad or mild detergent.
4. Allow the sensor to dry prior to placement on a patient.

or

1. If low-level disinfection is required, wipe all surfaces of the RD SET TF-I sensor and cable with a cloth or gauze pad saturated with a 1:10 bleach/water solution.
2. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the RD SET TF-I sensor and cable.
3. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

### To clean or disinfect the sensor using a soaking method

1. Place the sensor in the cleaning solution (1:10 bleach/water solution), so that the sensor and desired length of cable are completely immersed.

**WARNING:** To prevent damage, do not soak or immerse the connector in any liquid solution.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours. Do not immerse the connector.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable with a clean cloth or dry gauze pad.

## SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, or with licensed Masimo SET pulse oximetry modules and patient cables, the RD SET TF-I sensors have the following specifications:

Sensor:	RD SET TF-I
 Body Weight	> 30 kg
Application Site	Forehead
SpO <sub>2</sub> Accuracy, No Motion <sup>1</sup>	2%
Pulse Rate Accuracy, No Motion <sup>2</sup>	3 bpm

**NOTE:** ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within  $\pm$  ARMS of the reference measurements in a controlled study.

<sup>1</sup> The Masimo SET<sup>®</sup> technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO<sub>2</sub> against a laboratory CO-Oximeter.

<sup>2</sup> The Masimo SET<sup>®</sup> technology has been validated for pulse rate accuracy for the range of 25–240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

## COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET<sup>®</sup> oximetry or pulse oximetry monitors licensed to use RD SET TF-I sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: [www.Masimo.com](http://www.Masimo.com)

## WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

## NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RD SET YI SENSORS.

## CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult instructions for use		Authorized representative in the European community	>	Greater than
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	<	Less than
	Manufacturer		Light Emitting Diode (LED). LED emits light when current flows through.		Storage humidity limitation
	Date of manufacture YYYY-MM-DD		Do not discard		Storage temperature range
	Non-sterile		Lot code		Do not use if package is damaged
	Not made with natural rubber latex		Catalogue number (model number)		Atmospheric pressure limitation
<b>Rx ONLY</b>	<b>Caution:</b> Federal law (U.S.A) restricts this device to sale by or on the order of a physician		Masimo reference number		Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eFU is not available for CE mark countries.</b>
	Mark of conformity to European medical device directive 42/93/EEC		Body weight		

Patents: <http://www.masimo.com/patents.htm>

RD SET is a trademark of Masimo Corporation.

Masimo, SET, Masimo SET, TF-I, X-Cal and  are federally registered trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH