
# RD SET™ E1® Ear Sensor

## Single Patient Use Ear Sensor







Images	2
<b>en</b> English	3-7
<b>fr</b> French	8-12
<b>de</b> German	13-17
<b>it</b> Italian	18-22
<b>es</b> Spanish	23-27
<b>sv</b> Swedish	28-32
<b>nl</b> Dutch	33-37
<b>da</b> Danish	38-42
<b>pt</b> Portuguese	43-47
<b>zh</b> Chinese	48-51
<b>ja</b> Japanese	52-56
<b>fi</b> Finnish	57-61
<b>no</b> Norwegian	62-66
<b>cs</b> Czech	67-71
<b>hu</b> Hungarian	72-76
<b>pl</b> Polish	77-81
<b>ro</b> Romanian	82-86
<b>sk</b> Slovak	87-91
<b>tr</b> Turkish	92-96
<b>el</b> Greek	97-102
<b>ru</b> Russian	103-107
<b>ko</b> Korean	108-111
<b>ar</b> Arabic	115-112
Performance Specifications	116

# RD SET™ E1® Ear Sensor

## Single Patient Use Ear Sensor

# RD SET™ E1® Ear Sensor

en

## Single Patient Use Ear Sensor

### DIRECTIONS FOR USE

Single patient use only

Not made with natural rubber latex

Non-sterile

**Prior to using this sensor, the user should read and understand the Operator's Manual for the monitor and this Directions for Use.**

#### INDICATIONS

The RD SET™ E1® Ear Sensor is indicated for single patient use for continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult and pediatric patients, (weighing >30 kg), who are well or poorly perfused, in hospitals, hospital-type facilities, mobile, and home environments.

#### CONTRAINDICATIONS

The RD SET E1 Ear Sensor is contraindicated for patients who exhibit allergic reactions to foam rubber products and/or adhesive tape or for patients with signs of redness, swelling, infection or skin breakdown at the sensor application site, including the inner aspect of the ear and behind the ear. The sensor cannot be placed over areas of the ear that have piercings.

Refer to **Fig. 1**. The sensor should not be used on patients if the inner aspect of the ear (cavum conchae [A]) (reference is also made to the tragus [B], ear lobe [C] and crus of helix [D].) is not large enough to accommodate the emitter limb of the sensor without touching the tragus and/or crus of helix. The sensor is not intended for use on a monitoring site for prolonged periods. If extended monitoring is required, the sensor must be moved to the opposite ear. If the sensor cannot be moved, use of a Masimo adhesive sensor on a peripheral site, such as the digit, is recommended.

Because individual skin conditions and perfusion levels affect the ability of the site to tolerate sensor placement, it may be necessary to move the sensor frequently.

#### DESCRIPTION

The RD SET E1 Ear Sensor are for use only with devices containing Masimo SET® MS-2000 (Version 4.8 or higher) technology, Masimo rainbow® SET® MX technology.

The RD SET E1 Ear Sensor has been verified using Masimo SET Oximetry technology.

**WARNING:** Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

#### WARNINGS

- Do not use the E1 Ear Sensor on any site other than the inner aspect of the ear (cavum conchae). This may result in inaccurate readings due to tissue thickness.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor, Trendelenburg position).
- Venous pulsations may cause erroneous low SpO<sub>2</sub> readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.

- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High levels of COHb or MethHb may occur with a seemingly normal SpO<sub>2</sub>. When elevated levels of COHb or MethHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO<sub>2</sub> measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO<sub>2</sub> measurements.
- Elevated Total Bilirubin levels may lead to inaccurate SpO<sub>2</sub> measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring (such as indelible ink) may lead to inaccurate SpO<sub>2</sub> measurements.
- Inaccurate SpO<sub>2</sub> readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize the sensor.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reuse on multiple patients, reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 168 hours of patient monitoring. After single-patient use, discard sensor.

## INSTRUCTIONS

### A) Site Selection

**CAUTION:** Before using the sensor, ensure that the sensor is physically intact, with no broken or frayed wires or damaged parts.

- Refer to **Fig. 1**. The preferred measuring site is the hollow of the inner aspect of the ear (cavum conchae [A]). (Reference is also made to the tragus [B], ear lobe [C] and crus of helix [D].)
- The site should be free of debris. Prior to sensor placement, the site should be checked to ensure that it is clean, without signs of redness, swelling, infection or breakdown.

### B) Attaching the sensor to the patient

1. Open the pouch and remove the sensor.
2. Refer to **Fig. 2**. Orient the sensor to ensure that the star (1) on the emitter limb is in front of the patient's ear. Gently pull the emitter and detector limb apart to separate, so they will clear the ear lobe. Place the emitter limb inside the ear so that it rests on the hollow of the inner aspect of the ear (cavum conchae).
3. Place the detector limb on the back of the ear (2). The detector limb pad should lay flat on the back of the ear, without any part of the pad folding under.
4. Position the bottom of the sensor notch (3) so that it fits under the smallest portion of the ear lobe. For the patient's comfort, the E1 Ear Sensor can be adjusted after it is placed. Ensure that the sensor is not compressing any skin. It is recommended that the emitter not touch the tragus.
5. If the E1 Ear Sensor does not fit properly on the ear, consider using a Masimo reusable or adhesive sensor on another measuring site.
6. Refer to **Fig. 3**. The patient should turn their head to the opposite side from the ear where the sensor is attached. Attach the anchor tab on the cable to the anchor pad. Remove the release liner from the anchor pad. Secure the pad to the skin on the patient's shoulder area. Do not secure to the patient's clothing.

C) Attaching the Sensor to the Patient Cable

- 1. Refer to **Fig. 4**. Orient the sensor's connector tab so that the side with the "shiny" contacts is facing up. Orient the patient cable with the color bar and finger grips facing up.
- 2. Refer to **Fig. 5**. Insert the sensor tab into the patient cable until there is a tactile or audible click of connection. Gently tug on the connectors to ensure a positive contact. Tape may be used to secure the cable to the patient for ease of movement.

**Note:** When changing application sites, or reattaching sensor, first attach the sensor to the application site, then connect the patient cable to the sensor.

D) Stress and Exercise Testing

- 1. Attach the E1 Ear Sensor to the patient following the steps in B (Attaching the Sensor to the Patient).  
To minimize unnecessary sensor motion during excessive patient movement, use a headband (not shown). Loop the ear sensor cable under the chin, and secure the cable under the headband on the side of the head opposite the ear sensor.

E) Disconnecting the Sensor

From the Patient Cable:

- 1. Refer to **Fig. 6**. Pull firmly on the sensor connector to remove it from the patient cable.  
**Note:** To avoid damage, pull on the sensor connector, not the cable.

From the Patient:

- 1. Pull gently on the anchor pad and remove from the patient's shoulder.
- 2. To remove the E1 Ear Sensor from the patient's ear, gently pull apart the emitter and detector and remove.  
**CAUTION:** To avoid discomfort to the patient or damage to the sensor, do not remove the sensor from the patient's ear by pulling on the cable.  
**CAUTION:** Avoid positioning the patient so that external pressure is applied to the sensor at the measuring site.  
**CAUTION:** To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not sterilize by irradiation, steam, autoclave or ethylene oxide.

SPECIFICATIONS

When used with Masimo RD SET pulse oximetry monitors, Masimo RD rainbow SET pulse CO-Oximetry monitors, or with licensed Masimo SET pulse oximetry or Masimo rainbow SET pulse CO-Oximetry modules, and patient cables, during no motion, the RD E1 Ear sensor has the following specifications:

Sensor	RD SET E1 Ear Sensor
 Body Weight	> 30 kg
Application Site	Ear
SpO2 Accuracy, No Motion <sup>1</sup>	2.5%
SpO2 Accuracy, Low Perfusion <sup>2</sup>	2.5%
Pulse Rate Accuracy, No Motion <sup>3</sup>	3 bpm
Pulse Rate Accuracy, Low Perfusion <sup>3</sup>	3 bpm

**NOTE:** The ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within +/- ARMS of the reference measurements in a controlled study.

<sup>1</sup> The Masimo SET technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70%-100% SpO2 against a laboratory CO-Oximeter.

<sup>2</sup> The Masimo SET technology has been validated for low perfusion accuracy in bench top testing against a BioteK Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

<sup>3</sup> The Masimo SET technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a BioteK Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

---

## COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use RD SET E1 sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance. For Compatibility Information Reference: [www.Masimo.com](http://www.Masimo.com)

---

## ENVIRONMENTAL

Operating Temperature	+5°C to +40°C
Storage Temperature	-40°C to +70°C
Relative Humidity	10% to 95% noncondensing

---

## WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

---

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

---

## NO IMPLIED LICENSE

THIS SINGLE-PATIENT SENSOR IS LICENSED TO YOU UNDER THE PATENTS OWNED BY MASIMO FOR SINGLE-PATIENT USE ONLY. BY ACCEPTANCE OR USE OF THIS PRODUCT, YOU ACKNOWLEDGE AND AGREE THAT NO LICENSE IS GRANTED FOR USE OF THIS PRODUCT WITH MORE THAN A SINGLE PATIENT.

AFTER SINGLE-PATIENT USE, DISCARD SENSOR. PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE RD SENSORS.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**  
For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

**The following symbols may appear on the product or product labeling:**

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE)	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician.
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use by YYYY-MM-DD		Storage humidity limitation		Storage temperature range
	Do not re-use/ Single patient use only		Do not use if package is damaged		Keep dry
	Non-Sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Greater than		Less than		Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eIFU is not available in all countries.</b>

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , E1, X-Cal, Pulse CO-Oximeter and rainbow are federally registered trademarks of Masimo Corporation.

RD SET and rainbow SET are trademarks of Masimo Corporation.

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© 2019 Masimo Corporation
MDSS GmbH