
M-LNCS™/LNCS®/LNOP® YI









SpO₂ Multisite Reusable Sensor and Single Use Attachment Wraps

Images	2
en English	5-10
fr French	11-16
de German	17-22
it Italian	23-28
es Spanish	29-34
sv Swedish	35-40
nl Dutch	41-46
da Danish	47-52
pt Portuguese	53-58
zh Chinese	59-64
ja Japanese	65-70
fi Finnish	71-76
no Norwegian	77-82
cs Czech	83-88
hu Hungarian	89-94
pl Polish	95-100
ro Romanian	101-106
sk Slovak	107-112
tr Turkish	113-118
el Greek	119-124
ru Russian	125-130
ko Korean	131-135
ar Arabic	141-136
Performance Specifications	142














M-LNCS™/LNCS®/LNOP® YI

SpO₂ Multisite Reusable Sensor and Single Use Attachment Wraps

M-LNCS™/LNCS®/LNOP® YI

SpO₂ Multisite Reusable Sensor and Single Use Attachment Wraps

M-LNCS™/LNCS®/LNOP® YI

SpO₂ Multisite Reusable Sensor and Single Use Attachment Wraps


Fig. 10


Fig. 11a


Fig. 11b


Fig. 12a


Fig. 12b


Fig. 13

	-40°C – +70°C
	5% – 95%

M-LNCS™/LNCS®/LNOP® YI

en

SpO₂ Multisite Reusable Sensor and Single Use Attachment Wraps

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

INDICATIONS

The M-LNCS™/LNCS®/LNOP® YI Reusable Sensor is indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

CONTRAINDICATIONS

The M-LNCS/LNCS/LNOP YI is contraindicated for patients who exhibit allergic reactions to foam urethane products and/or adhesive tape.

DESCRIPTION

The YI sensor is applied to the sensor site using Masimo attachment wraps. The attachment wraps are for single-patient use only.

The M-LNCS/LNCS/LNOP YI sensor is for use only with devices containing Masimo SET® oximetry or licensed to use M-LNCS/LNCS/LNOP sensors. The Masimo Attachment Wraps are for use only with the M-LNCS, LNCS or LNOP YI Reusable Multisite Sensors. Consult individual device manufacturer for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether its devices are compatible with each sensor model. The YI series has been verified with Masimo SET Oximetry Technology.

The sensor must be removed and the site inspected at least every four (4) hours or sooner, and, if indicated by circulatory condition or skin integrity, reapplied to a different monitoring site.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- Circulation distal to the sensor site should be checked routinely.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor, Trendelenburg position).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid valve regurgitation, Trendelenburg position).
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.

- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or the unit might read zero for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High levels of COHb or MethHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MethHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated Total Bilirubin levels may lead to inaccurate SpO₂ measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate SpO₂ measurements.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- High oxygen concentrations may predispose a premature infant to retinopathy. Therefore, the upper alarm limit for the oxygen saturation must be carefully selected in accordance with accepted clinical standards.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal™ technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A) Site Selection

- Always choose a site that will completely cover the sensor's detector window.
- The site should be free of debris prior to sensor placement.
- Choose a site that is well perfused and least restricts a conscious patient's movements.
- The sensor is not intended for placement on the ear, if the ear is the desired monitoring site the Masimo M-LNCS/LNCS/LNOP TC-I reusable sensor is recommended.

NEONATES (1-3 kg), Standard Wrap / CleanShield Wrap / Foam Wrap:

The preferred sites are the outer aspect of the foot, under the fourth toe or the outer aspect of the palm of the hand.

INFANTS (3-10 kg), Standard Petite Wrap / Foam Wrap:

The preferred sites are the outer aspect of the foot, under the fifth toe or the outer aspect of the palm of the hand, under the fifth finger. For infants with fat or edematous feet, the great toe or thumb is recommended

ADULTS (> 30 kg) and PEDIATRICS (10-50 kg), Standard Wrap / CleanShield Wrap / Petite Wrap / Foam Wrap:

The preferred sites are the middle or ring finger of the non-dominant hand. Alternatively, the the great or second toe may be used.

B) Attachment Wraps

Body Weight	Clean Shield® Multisite Wrap	Standard Petite Wrap	Standard Wrap	Foam Wrap
1 kg ~ 3 kg	■		■	■
3 kg ~ 10 kg	■		■	■
10 kg ~ 30 kg		■		■
10 kg ~ 50 kg		■		■
> 30 kg	■	■	■	■

C) Attaching the adhesive squares to the sensor (not required for CleanShield Multisite Wrap)

NOTE: For improved adherence of the adhesive squares to the sensor wipe the sensor pads with 70% isopropyl alcohol and allow to dry prior to attaching the adhesive squares.

1. Refer to **Fig. 1a**. Remove the adhesive squares from the backing
2. Refer to **Fig. 1b**. Attach one square to each window of the sensor pads (emitter and detector). Avoid touching the sticky side prior to applying to the sensor pads.
3. Do not remove the release liner until ready to apply the sensor to the site.

CAUTION: Do not use adhesive squares on fragile skin.

D) Attaching the Sensor to the Patient

Refer to **Fig. 2**. The foam wrap can be trimmed to the shorter length for smaller site applications such as a child's finger or toe or a pre-term infant's foot or hand.

1. Remove the backing from the adhesive wrap, leaving the backing on the tab end (Standard Wrap, CleanShield Wrap and Standard Petite Wrap only).
2. Refer to **Fig. 3a**. Push the "button" on the emitter sensor pad (cable is marked with a red indicator) through the hole on the left **Fig. 3b**, and the remaining button on the detector sensor pad through the hole on the right **Fig. 3c**.

NOTE: Use the center hole of the Petite Wrap for patients 10 – 30 kg. Refer to **Fig. 3d**.

3. Refer to **Fig. 4**. Fold the upper portion of the attachment tape down over the sensor pads and remove the release liner from the folded over section of the adhesive wrap prior to applying the sensor to the site (CleanShield Wrap only).

NEONATES (1-3 kg) and INFANTS (3-10 kg), Standard Wrap / CleanShield Wrap / Foam Wrap

1. Direct the sensor cable toward the patient. Orient the YI on the outer aspect of the foot **Fig. 5a**, or hand **Fig. 5b**, with the center of the wrap **Fig. 5c**, under the 5th digit and the detector window on the fleshy portion and the emitter window (cable is marked with a red indicator) directly opposite.
2. Secure the small tape end and remove the remaining backing (Standard Wrap, CleanShield Wrap and Standard Petite Wrap only).
3. Wrap the tape loosely enough to avoid restricting circulation around the site and to maintain proper alignment of the detector and emitter windows; foot **Fig. 6a**, hand **Fig. 6b**.
4. If using the Foam Wrap, secure the Foam Wrap **Fig. 6c**, using the hook and loop tab.

NOTE: When placing the YI on the hand, it is recommended that the EMITTER (cable is marked with a red indicator) pad be placed on the palm of the hand, under the 4th or 5th finger.

ADULTS (> 30 kg) and PEDIATRICS (10-50 kg), Standard Wrap / CleanShield Wrap / Petite Wrap / Foam Wrap

1. Refer to **Fig. 7a**. Direct the sensor cable toward the patient. Orient the YI on the selected site so the detector window is on the fleshy portion and the emitter window (cable is marked with a red indicator) directly opposite on the nail bed.
2. Refer to **Fig. 7b**. Secure the small tape end and remove the remaining backing (Standard Wrap, CleanShield Wrap and Standard Petite Wrap only).
3. Refer to **Fig. 7c**. Wrap the adhesive or foam wrap around the finger loosely enough to avoid restricting circulation around the site and to maintain proper alignment of the detector and emitter windows.

NOTE: When placing the YI on a finger or toe, the EMITTER (cable is marked with a red indicator) pad must be placed on the NAIL side.

CAUTION: Replace the attachment wrap whenever the adhesive does not adequately adhere the sensor to the patient or when the attachment wrap becomes soiled from excessive use.

NOTE: Clip the clothing clip to the patient gown to reduce cable pulling on the measuring site.

E) Attaching the Sensor to the Patient Cable

M-LNCS

1. Refer to Fig. 8a. Orient sensor connector to patient cable connector as shown.
2. Refer to Fig. 8b. Insert sensor connector completely into patient cable connector. Close sensor connector cover over patient cable connector until it locks in place.

LNCS

1. Refer to Fig. 9a. Insert the sensor connector completely into the patient cable connector.
2. Refer to Fig. 9b. Completely close the protective cover.

LNOP

1. Refer to Fig. 10. Orient the connecting tab so that the "shiny" contacts are facing up and mate the logo to the logo on the patient cable.
2. Insert the LNOP YI connector into the patient cable connector until there is a tactile or audible click of connection. Gently tug on the connectors to ensure a positive connection. Tape may be used to secure the cable to the patient for ease of movement.

F) Disconnecting the sensor from the patient cable

M-LNCS

1. Refer to Fig. 11a. While holding the sides of patient cable connector, lift sensor connector cover from patient cable connector.
2. Refer to Fig. 11b. Pull firmly on the sensor connector to remove from the patient cable.

LNCS

1. Refer to Fig. 12a. Lift the protective cover to gain access to the sensor connector.
2. Refer to Fig. 12b. Pull firmly on the sensor connector to remove from the patient cable.

LNOP

1. Place thumb and index finger on buttons on either side of the patient cable connector.
2. Refer to Fig. 13. Press firmly on the buttons and pull to remove the sensor.

NOTE: When changing the sensor site disconnect the sensor from the patient cable and reconnect once the sensor is positioned on the new site.

CLEANING

CAUTION:

- Do not use undiluted bleach (5%-5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.
- Do not immerse the connector of the YI cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, Gluteraldehyde (Cidex) or ethylene oxide.
- Using excessive force when removing the attachment wrap may damage the sensor.

To surface clean the sensor

1. Remove the sensor from the patient and disconnect it from the attachment wrap and patient cable. Remove the adhesive squares.
2. Clean the YI sensor by wiping it with a 70% isopropyl alcohol pad or mild detergent.
3. Allow the sensor to dry prior to placement on a patient.

or

1. If low-level disinfection is required, wipe all surfaces of the YI sensor and cable with a cloth or gauze pad saturated with a 1:10 bleach/water solution.
2. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the YI sensor and cable.
3. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

To clean or disinfect the sensor using a soaking method

1. Place the sensor in the cleaning solution (1:10 bleach/water solution), so that the sensor and desired length of cable are completely immersed.

WARNING: Do not immerse the connector end of the sensor cable as this may damage the sensor.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours. Do not immerse the connector.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable with a clean cloth or dry gauze pad.

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, or with licensed Masimo SET pulse oximetry modules and patient cables the YI Sensor has the following specifications:

YI Sensor	Neonates	Adult / Pediatric / Infants
Weight Range	1 - 3 kg	>3 kg
SpO ₂ Accuracy, No Motion, (70 - 100%) ¹	3%	2%
SpO ₂ Accuracy, Motion ²	3%	3%
SpO ₂ Rate Accuracy, No Motion, (25 - 240 bpm) ³	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ³	5 bpm	5 bpm
SpO ₂ Accuracy, Low Perfusion ⁴	3%	2%
Pulse Rate Accuracy, Low Perfusion ⁴	3 bpm	3 bpm

NOTE: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within +/- ARMS of the reference measurements in a controlled study.

SpO₂ accuracy was determined by testing on healthy adult volunteers with light to dark skin in the range of 70% -100% SpO₂ against a laboratory co-oximeter. Accuracy specifications are statistically distributed, and only about two-thirds of the measurements fall within the 1 Std. Dev. specification.

¹ The Masimo rainbow SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

² The Masimo rainbow SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70%-100% SpO₂ against a laboratory co-oximeter.

³ The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁴ The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use YI sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

 For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.



NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE YI SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Consult instructions for use		Lot code
	Follow instructions for use		Catalogue number (model number)
	Manufacturer		Masimo reference number
	Date of Manufacture YYYY-MM-DD		Do not discard
	Use by YYYY-MM-DD		Pediatric patient
	Non-sterile		Body weight
	Not made with natural rubber latex		Greater than
Rx ONLY	Federal law (USA) restricts this device to sale by or on the order of a physician		Less than
	Mark of Conformity to European Medical Device Directive 93/42/EEC		Storage humidity limitation
	Authorized representative in the European community		Storage temperature range
	Separate collection for electrical and electronic equipment (WEEE).		Keep dry
	Do not use if package is damaged		Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eFU is not available for CE mark countries.

Patents: referenced at www.masimo.com/patents.htm. Other patents pending.

M-LNCS and X-Cal are trademarks of Masimo Corporation.

Masimo, SET, , CleanShield, LNCS and LNOP are federally registered trademarks of Masimo Corporation.

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© 2019 Masimo Corporation
MDSS GmbH