# RD SET®

## Adapter Cable Series






Images	2-3
<b>en</b> English	4-6
<b>fr</b> French	7-9
<b>de</b> German	10-12
<b>it</b> Italian	13-15
<b>es</b> Spanish	16-18
<b>sv</b> Swedish	19-21
<b>nl</b> Dutch	22-24
<b>da</b> Danish	25-27
<b>pt</b> Portuguese	28-30
<b>zh</b> Chinese	31-33
<b>ja</b> Japanese	34-36
<b>fi</b> Finnish	37-39
<b>no</b> Norwegian	40-42
<b>cs</b> Czech	43-45
<b>hu</b> Hungarian	46-48
<b>pl</b> Polish	49-51
<b>ro</b> Romanian	52-54
<b>sk</b> Slovak	55-57
<b>tr</b> Turkish	58-60
<b>el</b> Greek	61-63
<b>ru</b> Russian	64-66
<b>ar</b> Arabic	69-67
<b>fa</b> Farsi	72-70






# RD SET®

## Adapter Cable Series

# RD SET®

## Adapter Cable Series

# RD SET<sup>®</sup>

## Adapter Cable Series

en

### DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

#### INDICATIONS

The RD SET<sup>®</sup> Adapter Cable Series have the same indications for use as the applicable sensors. Please see the compatible sensor's directions for use for indications for use and prescribed information.

#### DESCRIPTION

The RD SET Adapter cables allow Masimo RD SET sensors to be used on PC or M-LNC and RC patient cables on Masimo<sup>®</sup> SET<sup>®</sup> or Masimo rainbow<sup>®</sup> SET devices. The RD SET Adapter cables allow Masimo<sup>®</sup> SET<sup>®</sup> LNOP<sup>®</sup>, LNCS<sup>®</sup> or M-LNCS<sup>™</sup> sensors to be used with RD SET series patient cables. RD EXT cables extend the length of RD SET and RD rainbow SET series patient cables.

**RD SET Adapter Cable Series consist of the following:**

- RD to PC Adapter Cable
- LNCS to RD Adapter Cables
- RD SET EXT-4
- RD to M-LNC Adapter Cable
- M-LNCS to RD Adapter Cables
- LNCS to RD Adapter Cables

The RD SET Adapter Cable Series have been verified using Masimo SET or Masimo rainbow SET technology.

**WARNING:** Masimo sensors and cables are designed for use with devices containing Masimo SET or Masimo rainbow SET oximetry or licensed to use Masimo sensors.

**NOTE:** When using the RD SET EXT cable with a RD rainbow SET series patient cable, the maximum length for the sensor, cable and extension cable combination is 15 feet.

#### WARNINGS

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully route patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the sensor or the oximeter module to the cable will result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- To prevent damage, do not soak or immerse the sensor or cable in any liquid solution. Do not attempt to sterilize the cable or sensor.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.

#### INSTRUCTIONS

##### A) Attaching the RD SET Adapter Cable to a Patient Cable

###### RD to PC Adapter Cable

1. Refer to **Fig. 1**. Orient the connecting tab so that the "shiny" contacts are facing up and mate to the patient cable.
2. Insert the connector into the patient cable connector until there is a tactile or audible click of connection.
3. Gently tug on the connectors to ensure a positive contact.

###### RD to M-LNC Adapter Cable

1. Refer to **Fig 2**. Properly orient and insert the connector completely into the patient cable connector.
2. Refer to **Fig. 3**. Completely close the protective latch.

###### PC to RD Adapter Cable / LNCS to RD Adapter Cable / M-LNCS to RD Adapter Cable / RD SET EXT-4:

1. Refer to **Fig. 4**. Properly orient the connector and insert the connector completely into the patient cable connector.
2. Refer to **Fig. 5**. Close the protective latch completely.

## B) Disconnecting the RD SET Adapter Cable from a Patient Cable

### RD to PC Adapter Cable

1. Refer to **Fig. 6**. Place thumb and index finger on the buttons on either side of the connector.
2. Press firmly on the buttons and pull to remove.

### RD to a M-LNC Adapter Cable

1. Refer to **Fig. 7**. Lift the protective latch to gain access to the connector.
2. Refer to **Fig. 8**. Pull firmly on the connector to remove from the patient cable. To avoid damage, pull on the connector, not the cable.

### PC to RD Adapter Cable / LNCS to RD Adapter Cable / M-LNCS to RD Adapter Cable / RD SET EXT-4:

1. Refer to **Fig. 9**. Lift up the protective latch.
2. Refer to **Fig. 10**. Pull firmly on the connector to remove it from the patient cable. To avoid damage, pull on the connector, not the cable.

---

## CLEANING

1. Remove the sensor from the patient and disconnect it from the adapter cable.
2. Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad.
3. Wipe all surfaces of the cable.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.
5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

## CAUTION:

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

---

## ENVIRONMENTAL

Operating Temperature	41°F to 104°F (5°C to 40°C)
Storage Temperature	-40°F to 158°F (-40°C to 70°C)
Relative Humidity	10% to 95% non-condensing

---

## WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

---

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

---

## NO IMPLIED LICENSE

Purchase or possession of this RD SET Adapter Cable does not carry any express or implied license to use this RD SET Adapter Cable with any device that is not an authorized device or separately authorized to use RD SET Adapter Cables.

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**  
 For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	<b>Rx ONLY</b>	<b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		European Union Conformity Mark
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Caution		Storage humidity limitation		Storage temperature range
	Do not discard		Do not use if package is damaged and consult instructions for use		Keep dry
	Fragile, handle with care		Greater than		Less than
	Non-sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Medical device		Unique device identifier		Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a> <b>Note: eIFU is not available in all countries.</b>
	Importer		Distributor		

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , LNCS, LNOP, RD SET, M-LNCS, rainbow SET and rainbow are federally registered trademarks of Masimo Corporation.

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© 2021 Masimo Corporation
MDSS GmbH