Masimo® NK Series











Patient Cables

| | |
|----------------------|-------|
| Images | 2-3 |
| en English | 4-6 |
| fr French | 7-9 |
| de German | 10-12 |
| it Italian | 13-15 |
| es Spanish | 16-18 |
| sv Swedish | 19-21 |
| nl Dutch | 22-24 |
| da Danish | 25-27 |
| pt Portuguese | 28-30 |
| zh Chinese | 31-33 |
| ja Japanese | 34-36 |
| fi Finnish | 37-39 |
| no Norwegian | 40-42 |
| cs Czech | 43-45 |
| hu Hungarian | 46-48 |
| pl Polish | 49-51 |
| ro Romanian | 52-54 |
| sk Slovak | 55-57 |
| tr Turkish | 58-60 |
| el Greek | 61-63 |
| ru Russian | 64-66 |
| ar Arabic | 69-67 |
| fa Farsi | 72-70 |

Masimo® NK Series Patient Cables


Fig. 11


Fig. 12


Fig. 13


Fig. 14


Fig. 15


Fig. 16


Fig. 17

Masimo® NK Series

Patient Cables

en

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-Sterile

INDICATIONS

The Masimo® NK Series patient cables have the same indications for use as the applicable sensors. Please see the compatible sensor's directions for use for indications for use and prescribed information.

DESCRIPTION

The Masimo NK Series (LNC, M-LNC, and PC) patient cables are used with LNCS®, M-LNCS® and LNOP® sensors.

The Masimo NK Series patient cables and Masimo sensors have been functionally verified on Nihon Kohden devices with Masimo SET® technology.

Masimo NK Series patient cables have been verified using Masimo SET technology.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully route patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the sensor or the oximeter module to the cable may result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- To prevent damage, do not soak or immerse the sensor or cable in any liquid solution. Do not attempt to sterilize the cable or sensor.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- **Caution:** Replace the cable when a replace cable message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The cable is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The cable will provide up to 17,520 hours of patient monitoring time. Replace the cable when the patient monitoring time is exhausted. Some legacy monitoring devices do not distinguish between the cable and the sensor. In these devices, a replace sensor or similar message could pertain either to the sensor or the cable.

INSTRUCTIONS

A) Attach the Masimo NK Series Patient Cable to the Device

1. Orient the connector of the cable to mate with the patient cable connector on the device and insert.

B) Attach the Masimo NK Series Patient Cable connector to the sensor connector

Attaching to a LNCS sensor

1. Refer to **Fig. 1**. Insert the sensor connector completely into the patient cable connector.
2. Refer to **Fig. 2**. Completely close the protective cover.

Attaching to a M-LNCS disposable sensor

1. Refer to **Fig. 3**. Orient sensor connector to the patient cable connector as shown.
2. Refer to **Fig. 4**. Insert the sensor connector completely into the patient cable connector until it locks in place as shown.

Attaching to a rainbow® disposable sensor

1. Refer to **Fig. 5**. Orient sensor connector (1) to patient cable connector (2) as shown.
2. Refer to **Fig. 6**. Insert the sensor connector completely into the patient cable connector until it locks in place as shown.

Attaching to a rainbow or M-LNCS reusable sensor

NOTE: rainbow spot check reusable sensors are not for use with the Masimo NK Series patient cable.

1. Refer to **Fig. 7**. Orient sensor connector (1) to patient cable connector (2) as shown.
2. Refer to **Fig. 8**. Insert sensor connector (1) completely into patient cable connector (2). Close sensor connector cover (3) over patient cable connector (2) until it locks in place.

Attaching to a LNOP reusable or disposable sensor

1. Refer to **Fig. 9**. Orient the connecting tab of the sensor so that the “shiny” contacts are pointed up. Mate the logo on the sensor to the logo on the patient cable.
2. Refer to **Fig. 10**. Insert the tab of the sensor into the patient cable connector until there is a tactile or audible click of connection. Gently tug on the patient cable connector to ensure a positive contact. Tape may be used to secure the cable to the patient for ease of movement.

C) Disconnecting the Masimo NK Series Patient Cable from the sensor connector

Disconnecting from a LNCS sensor

1. Refer to **Fig. 11**. Lift the protective cover to gain access to the sensor connector.
2. Refer to **Fig. 12**. Pull firmly on the sensor connector to remove it from the patient cable.

Disconnecting from a M-LNCS disposable sensor

1. Refer to **Fig. 13**. Pull firmly on the sensor connector to remove it from the patient cable.

Disconnecting from a rainbow® disposable sensor

1. Refer to **Fig. 14**. While holding the sides of the cable connector, place thumb on latch (1) and press. Hold the latch down while pulling the connectors apart.

Disconnecting from a rainbow or M-LNCS reusable sensor

1. Refer to **Fig. 15**. While holding the sides of patient cable connector, lift sensor connector cover from patient cable connector.
2. Refer to **Fig. 16**. Pull firmly on the sensor connector to remove from the patient cable.

Disconnecting from a LNOP reusable or disposable sensor

1. Refer to **Fig. 17**. Place thumb and index finger on the buttons on either side of the patient cable connector. Press firmly on the buttons and pull to remove it from the patient cable.

CLEANING

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad.
3. Wipe all surfaces of the cable.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.
5. Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

CAUTION

- Do not immerse the connector on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

ENVIRONMENTAL

| | |
|-----------------------|--------------------------------|
| Operating Temperature | 41°F to 104°F (5°C to 40°C) |
| Storage Temperature | -40°F to 158°F (-40°C to 70°C) |
| Relative Humidity | 10% to 95% noncondensing |

COMPATIBILITY

 Masimo NK Series patient cables are intended for use only with devices containing Masimo SET technology and with compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular devices and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor and cable model. Use of this cable with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

Purchase or possession of this NK patient cable does not carry any express or implied license to use this cable with any device that is not an authorized device or separately authorized to use NK patient cables.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

| SYMBOL | DEFINITION | SYMBOL | DEFINITION | SYMBOL | DEFINITION |
|---|------------------------------------|---|---|---|---|
|  (blue background) | Follow instructions for use |  | Separate collection for electrical and electronic equipment (WEEE). | Rx ONLY | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician |
|  | Consult instructions for use |  | Lot code |  | Mark of conformity to European Medical Device Directive 93/42/EEC |
|  | Manufacturer |  | Catalogue number (model number) |  | Authorized representative in the European community |
|  | Date of manufacture YYYY-MM-DD |  | Masimo reference number |  | Body weight |
|  | Do not discard |  | Caution |  | Storage temperature range |
|  | Non-sterile |  | Greater than |  | Keep dry |
|  | Not made with natural rubber latex |  | Less than |  | Do not use if package is damaged and consult instructions for use |
|  | Fragile, handle with care |  | Storage humidity limitation |  | Atmospheric pressure limitation |
|  | Medical device |  | Unique device identifier |  Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available in all countries. | |
|  | Importer |  | Distributor | | |

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, LNCS, LNOP, X-Cal, M-LNCS, rainbow and  are federally registered trademarks of Masimo Corporation.

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© 2021 Masimo Corporation
MDSS GmbH