# NomoLine-O™ Sampling Lines








## Single-patient-use Airway Adapter Sets and Cannulas

|                      |       |
|----------------------|-------|
| Images               | 2     |
| <b>en</b> English    | 3-4   |
| <b>fr</b> French     | 5-7   |
| <b>de</b> German     | 8-10  |
| <b>it</b> Italian    | 11-13 |
| <b>es</b> Spanish    | 14-16 |
| <b>sv</b> Swedish    | 17-19 |
| <b>nl</b> Dutch      | 20-22 |
| <b>da</b> Danish     | 23-25 |
| <b>pt</b> Portuguese | 26-28 |
| <b>zh</b> Chinese    | 29-30 |
| <b>ja</b> Japanese   | 31-33 |
| <b>fi</b> Finnish    | 34-35 |
| <b>no</b> Norwegian  | 36-38 |
| <b>cs</b> Czech      | 39-41 |
| <b>hu</b> Hungarian  | 42-44 |
| <b>pl</b> Polish     | 45-47 |
| <b>ro</b> Romanian   | 48-50 |
| <b>sk</b> Slovak     | 51-53 |
| <b>tr</b> Turkish    | 54-56 |
| <b>el</b> Greek      | 57-59 |
| <b>ru</b> Russian    | 60-62 |
| <b>ko</b> Korean     | 63-65 |
| <b>et</b> Estonian   | 66-68 |
| <b>lt</b> Lithuanian | 69-71 |
| <b>lv</b> Latvian    | 72-74 |
| <b>ar</b> Arabic     | 76-75 |
| <b>fa</b> Farsi      | 78-77 |

# NomoLine-O™ Sampling Lines

## Single-patient-use Airway Adapter Sets and Cannulas

# NomoLine-O™ Sampling Lines

en

## Single-patient-use Airway Adapter Sets and Cannulas

### DIRECTIONS FOR USE

② Single-Patient Use Only

⊗ Not made with natural rubber latex

#### INDICATIONS

The NomoLine® Product Family is indicated for the measurement of respiratory rate and respiratory and anesthetic gases in adult, pediatric, infant and neonatal patients. The NomoLine Product Family includes single-patient use sampling lines for gas sampling and/or oxygen delivery. The NomoLine Product Family is indicated for use by clinical professionals in healthcare environments, including mobile environments.

#### DESCRIPTION

The NomoLine Product Family includes nasal and nasal/oral cannulas for non-intubated patients and airway adapter sets for intubated patients. NomoLine-O sampling lines are provided with a connector that is compatible with side-stream gas analyzers that use a single twist luer connector with an optical identifying reflector. See Masimo website for a list of compatible devices.

#### WARNINGS

- Carefully route the NomoLine-O sampling line to reduce the risk of patient entanglement or strangulation.
- Do not apply negative pressure to remove condensed water.
- Dispose in accordance with local regulations for biohazardous waste.
- Do not submerge the NomoLine-O sampling line in any cleaning solution or attempt to sterilize by autoclave, irradiation, steam, gas, ethylene oxide or any other method as this may seriously damage the NomoLine-O sampling line.
- Do not re-use due to cross contamination risk.
- Always replace the NomoLine-O sampling line if an occlusion alarm occurs.
- Do not use NomoLine-O Cannulas on patients with obstructions in one or both nostrils.
- The CO<sub>2</sub> measurement can be diminished if NomoLine-O Cannulas are used during mouth breathing or on a patient that has a perforated septum.
- Do not use NomoLine-O Adult/Pediatric Airway Adapter Sets for infants/neonates as the adapter adds 6 ml dead space.
- Do not use NomoLine-O Infant/Neonatal Airway Adapter Sets for adults/pediatrics as this may cause excessive flow resistance (0.7 ml dead space).
- Do not use NomoLine-O Airway Adapter Sets with inhalers or nebulized medications as this may clog the bacteria filter.

#### INSTRUCTIONS

Open the bag and remove the NomoLine-O Sampling Line. For convenient preparations, do not yet remove the tape that keeps the tubing coiled.

##### A. Fitting a NomoLine-O Cannula to the patient

- If fitting a NomoLine-O Nasal CO<sub>2</sub> Cannula refer to **Fig. 1**. Insert the prongs into the nostrils.
- If fitting a NomoLine-O Nasal/Oral CO<sub>2</sub> Cannula refer to **Fig. 2**. Insert the prongs into the nostrils and position the oral scoop in front of the mouth.
- Refer to **Fig. 3**. Pass the cannula lines over the ears (1) and adjust the slider for comfort and fit under the chin (2).

##### B. Connecting a NomoLine-O Airway Adapter Set to the patient circuit

- Refer to **Fig. 4**. Connect the small diameter end of the airway adapter to the Y-piece of the ventilation source.
- Refer to **Fig. 5**. Connect the large diameter end of the airway adapter to the patient's endotracheal tube.

##### C. Connecting a NomoLine-O sampling line to the compatible gas analyzer

- Refer to **Fig. 6**. Remove the tape that keeps the tubing coiled. Turn the NomoLine-O sampling line connector clockwise into the device's gas sampling inlet until it can't be turned anymore.
- Ensure that the tubing is properly connected and that it is not twisted or crimped.

##### D. Connecting a NomoLine-O Cannula to an oxygen source

- Refer to **Fig. 7**. Remove the tape that keeps the tubing coiled. Attach the oxygen supply tubing to the oxygen source. Adjust the oxygen flow. (This step applies only to NomoLine-O Cannulas with O<sub>2</sub>).
- Ensure that the tubing is properly connected and that it is not twisted or crimped.

#### ENVIRONMENTAL

| Storage                 |                                                               |
|-------------------------|---------------------------------------------------------------|
| Temperature             | -40 – +70°C                                                   |
| Humidity                | 10–95% RH<br>(95% RH at 40°C)                                 |
| Atmospheric Pressure    | 20–120 kPa                                                    |
| Operation               |                                                               |
| Temperature             | 0 – +50°C                                                     |
| Humidity                | < 5 kPa H <sub>2</sub> O (non-condensing)<br>(41% RH at 50°C) |
| Ambient CO <sub>2</sub> | ≤ 800 ppm                                                     |
| Atmospheric Pressure    | 52.5–120 kPa (< 4572 m)                                       |

## REPLACEMENT

The NomoLine-O sampling lines are single-patient use products and should after use be disposed of in accordance with local biohazardous waste and replaced for each new patient with a maximum of 14 days of use\*. They should also be replaced if the sampling line becomes occluded. Occlusion is indicated by a message on the monitor. Replace the NomoLine-O sampling line and ensure that the tubing is properly connected and that it is not twisted or crimped.

\*The lifetime of the product depends on the application and product type.

## DISPOSAL

The NomoLine-O sampling line may after use contain accumulated moisture and fluids. Dispose the used product in accordance with local regulations for biohazardous waste.

## WARRANTY

MASIMO WARRANTS TO THE INITIAL BUYER ONLY THAT THESE PRODUCTS, WHEN USED IN ACCORDANCE WITH THE DIRECTIONS PROVIDED WITH THE PRODUCTS BY MASIMO, WILL BE FREE OF DEFECTS IN MATERIALS AND WORKMANSHIP FOR A PERIOD OF SIX (6) MONTHS. SINGLE USE PRODUCTS ARE WARRANTED FOR SINGLE-PATIENT USE ONLY. THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

## WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to products that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

**CAUTION:** FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

| SYMBOL         | DEFINITION                                                                                        | SYMBOL | DEFINITION                                                                                                                                                                                                       |
|----------------|---------------------------------------------------------------------------------------------------|--------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                | Follow Instructions for Use                                                                       |        | Masimo reference number                                                                                                                                                                                          |
|                | Manufacturer                                                                                      |        | Body weight                                                                                                                                                                                                      |
|                | Use By YYYY-MM-DD                                                                                 |        | Less than or equal to                                                                                                                                                                                            |
|                | Do not re-use/Single-Patient Use only                                                             |        | Storage Temperature Range                                                                                                                                                                                        |
|                | Not made with natural rubber latex                                                                |        | Storage Humidity Range                                                                                                                                                                                           |
| <b>Rx ONLY</b> | <b>Caution:</b> Federal law (USA) restricts this device to sale by or on the order of a physician |        | Storage Atmospheric Pressure Range                                                                                                                                                                               |
|                | Mark of Conformity to European Medical Device Directive 93/42/EEC                                 |        | Gas sample                                                                                                                                                                                                       |
|                | Keep Dry                                                                                          |        | Collects water. Short term applications with low humidity                                                                                                                                                        |
|                | Lot code                                                                                          |        | Removes aspirated and condensed water. Long term applications with high humidity                                                                                                                                 |
|                | Catalogue number (model number)                                                                   |        | Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.Masimo.com/TechDocs">http://www.Masimo.com/TechDocs</a><br><b>Note: eIFU is not available in all countries.</b> |
|                | Do not use if package is damaged                                                                  |        |                                                                                                                                                                                                                  |

Patents: <http://www.masimo.com/patents.htm>

Masimo, NomoLine, and are federally registered trademarks of Masimo Corporation.

NomoLine-O is a trademark of Masimo Corporation.

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© 2019 Masimo Corporation
Masimo Sweden AB