




M-LNC™ MP Patient Cable

Images	2
en English	3-7
fr French	8-12
de German	13-17
it Italian	18-22
es Spanish	23-27
sv Swedish	28-32
nl Dutch	33-37
da Danish	38-42
pt Portuguese	43-47
zh Chinese	48-51
ja Japanese	52-56
fi Finnish	57-61
no Norwegian	62-66
cs Czech	67-71
hu Hungarian	72-76
pl Polish	77-81
ro Romanian	82-86
sk Slovak	87-91
tr Turkish	92-96
el Greek	97-101
ru Russian	102-106
ko Korean	107-110
ar Arabic	114-111

M-LNC™ MP Patient Cable

M-LNC MP Patient Cable

en

M-LNC MP Series

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this Cable, the user should read and understand the Operator's Manual for the Device, this Directions for Use and the Directions for Use for the Sensor.

INDICATIONS

The M-LNC™ MP Series cables and Masimo sensors are indicated for the continuous, noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients in hospitals, hospital-type facilities, mobile, and home environments.

DESCRIPTION

The M-LNC MP Series cables are for use on Philips modules with Masimo SET® Technology and devices with Philips FAST-SpO₂ Technology.

On Philips modules with Masimo SET Technology, M-LNC Patient Cables are for use with M-LNCS™ adhesive and reusable sensors, rainbow® adhesive, disposable and rainbow ReSpasable™ sensors. If rainbow reusable sensors are connected, an Incompatible/Invalid message will display on the monitor.

On devices with Philips FAST-SpO₂ Technology, M-LNC MP Patient Cables are for use with M-LNCS adhesive and reusable sensors. rainbow adhesive, ReSpasable and reusable sensors are not compatible with Philips FAST-SpO₂ Technology. Masimo SET Technology is not available on devices with Philips FAST-SpO₂ Technology.

M-LNC MP Series Patient Cables and the above Masimo Sensors have been validated with Masimo SET technology and Philips FAST-SpO₂ technology. The saturation accuracy of the neonatal sensors were validated on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® oximetry or licensed to use Masimo sensors.

WARNINGS

- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discoloration are found.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Refer to the specific sensor's Directions for Use for detailed sensor instructions and warnings.
- Probes and cables are designed for use with specific monitors. The user and/or operator needs to verify the compatibility of the monitor, probe, and cable before use, otherwise patient injury can result.
- Do not attempt to reprocess, recondition or recycle Masimo patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Failure to properly connect the sensor or the oximeter module to the cable will result in intermittent readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either end.
- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.

INSTRUCTIONS

A) Attaching the M-LNC MP Patient Cable directly to the Instrument

Refer to **Fig. 1**. Orient the connector of the cable to mate with the patient cable connector on the Module.

Caution: Ensure you are connecting a M-LNC MP Patient Cable to either a Masimo SET or Philips FAST-SpO₂ cable connector on the Module. There are a number of similar connectors with different colors and different mechanical keying. Never force a patient cable connector into the module. Failure to use a M-LNC MP patient cable may result in damage to the module, inaccurate readings, or no readings.

Note: There is a difference in color shade between the M-LNC MP Patient Cable and the Philips Module patient cable connector. However, this is an acceptable configuration. Masimo SET performance is not available when connected to an instrument with Philips FAST-SpO₂ Technology.

B) Attaching the M-LNC MP Patient Cable to the Sensor

M-LNCS Reusable Sensor

1. Refer to **Fig. 2**. Orient the reusable sensor connector (1) to the M-LNC MP patient cable (2) as shown.
2. Insert the reusable sensor connector (1) completely into the M-LNC MP patient cable (2). Close the reusable sensor connector cover (3) over the M-LNC patient cable connector (2).

M-LNCS or rainbow Adhesive Sensor

1. Refer to **Fig. 3**. Orient the adhesive sensor connector (1) to the M-LNC MP patient cable connector (2) as shown.
2. Insert the adhesive sensor connector (1) completely into the M-LNC MP patient cable connector (2) until it locks in place.

Disconnecting the M-LNC MP Patient Cable from the Sensor

M-LNCS Reusable Sensor

1. Refer to **Fig. 4**. While holding the sides of the M-LNC MP patient cable (2), lift the sensor cover (3) from the reusable sensor.
2. Disconnect the reusable sensor from the M-LNC MP patient cable.

M-LNCS or rainbow Adhesive Sensor

1. Refer to **Fig. 5**. While holding the sides of the connectors, place thumb on the latch of the adhesive sensor connector and press.
2. Hold the latch down while pulling the connectors apart.

SPECIFICATIONS

M-LNC MP Patient Cable Compatibility:

When used with Masimo SET Technology:

Sensors	 Body Weight	Saturation Accuracy (70–100% SpO ₂)		Pulse Rate Accuracy (25–240 bpm)		Low Perfusion Accuracy (70–100% SpO ₂)	
		No Motion	Motion	No Motion	Motion	No Motion	Motion
M-LNCS Aadx/Aadx-3	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS Pdx/Pdx-3	10–50 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS Inf/Inf-3/Inf-L	3–20 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS Neo/Neo-3/ Neo L	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
	> 40 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS NeoPt/ NeoPt-3/NeoPt-L	< 1 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
M-LNCS Trauma ¹	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS Newborn Neonatal ¹	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
M-LNCS Newborn Infant/Pediatric ¹	3–10 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
	10–30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS Blue	2.5–3 kg	3% (for 80-100)	-	3 bpm	-	3%	3 bpm
		4% (for 60-80)					
		3.3% (for 70-100)					
M-LNCS E1	> 30 kg	2.5%	-	3 bpm	-	2.5%	3 bpm
M-LNCS DCI	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS DCIP	10–50 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
M-LNCS TC-I	> 30 kg	3.5%	N/A	3 bpm	N/A	3.5%	3 bpm
M-LNCS TF-I [®]	> 30 kg	2%	N/A	3 bpm	N/A	2%	3 bpm
M-LNCS DBI	> 30 kg	2%	N/A	3 bpm	N/A	2%	3 bpm
M-LNCS Y-I	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
	> 3 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
rainbow R1 25, R1 25L, R25, R25L	> 30 kg	2% (60–80% 3%)	3%	3 bpm	5 bpm	2%	3 bpm
rainbow R25L	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
rainbow R1 20, R20	10-50 kg	2% (60–80% 3%)	3%	3 bpm	5 bpm	2%	3 bpm
rainbow R1 20L, R20L	3-10 kg	2% (60–80% 3%)	3%	3 bpm	5 bpm	2%	3 bpm
	10-30 kg	2% (60–80% 3%)	3%	3 bpm	5 bpm	2%	3 bpm
rainbow ReSposable R2 20a, R2 20r	10–50 kg	2%	N/A	3 bpm	N/A	2%	3 bpm
rainbow ReSposable R2 25a, R2 25r	> 30 kg	2%	N/A	3 bpm	N/A	2%	3 bpm
M-LNCS S-ROS 3U	universal	See applicable S-DOS sensor below for accuracy specifications when used in combination with S-ROS optical sensor.					
Resposable S-DOS 25	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
Resposable S-DOS 20	10-50 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
Resposable S-DOS 25L	< 3 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
	> 30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
Resposable S-DOS 20L	3–10 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
	10–30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
Resposable S-DOS 30L	3–10 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
	> 40 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
Resposable S-DOS 25Pt-L	< 1 kg	3%	3%	3 bpm	5 bpm	3%	3 bpm
Resposable S-DOS 20L-500	3–10 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm
	10–30 kg	2%	3%	3 bpm	5 bpm	2%	3 bpm

Note: ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within $\pm$ ARMS of the reference measurements in a controlled study.

The M-LNCS TF-I and TC-I sensors were not validated under motion conditions.

When used with Philips FAST-SpO2 Technology:

Sensors	 Body Weight	Saturation Accuracy (70 - 100% SpO2)	Pulse Rate Accuracy (25 - 240 bpm)
		No Motion	No Motion
M-LNCS Adtx/Adtx-3	> 30 kg	2%	3 bpm
M-LNCS Pdtx/Pdtx-3	10 - 50 kg	2%	3 bpm
M-LNCS Inf/Inf-3/Inf-L	3 - 20 kg	2%	3 bpm
M-LNCS Neo/Neo-3/Neo L	< 3 kg	3%	3 bpm
	> 40 kg	2%	3 bpm
M-LNCS NeoPt/NeoPt-3/NeoPt-L	< 1 kg	3%	3 bpm
M-LNCS Trauma ¹	> 30 kg	2%	3 bpm
M-LNCS Newborn Neonatal ¹	< 3 kg	3%	3 bpm
M-LNCS Newborn Infant/Pediatric ¹	3 - 10 kg	2%	3 bpm
	10 - 30 kg	2%	3 bpm
M-LNCS DCI	> 30 kg	2%	3 bpm
M-LNCS DCIP	10 - 50 kg	2%	3 bpm
M-LNCS TC-I	> 30 kg	3.5%	3 bpm
M-LNCS Y-I	< 3 kg	3%	3 bpm
	> 3 kg	2%	3 bpm
M-LNCS DBI	> 30 kg	2%	3 bpm

¹Provides saturation and pulse rate values when used with Philips technology but does not operate in the specialty mode unless used with Masimo technology.

CLEANING

Clean the M-LNC MP Patient Cable by wiping it with a 70% isopropyl alcohol pad and allow it to dry.

CAUTION

Do not soak or immerse the cable in any liquid solution. Do not attempt to sterilize.

ENVIRONMENTAL

Operating Temperature	41°F to 104°F (5°C to 40°C)
Storage Temperature	-40°F to 158°F (-40°C to 70°C)
Relative Humidity	5% to 95% noncondensing

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended instrument or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

Purchase or possession of these patient cables does not carry any express or implied license to use these cables with any device that is not an authorized device or separately authorized to use M-LNC MP Patient Cable.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Caution		Storage humidity limitation		Storage temperature range
	Do not discard		Do not use if package is damaged		Keep dry
	Fragile, handle with care		Greater than		Less than
	Non-sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available for CE mark countries.				

Patents: <http://www.masimo.com/patents.htm>

Masimo, , SET, TF-I and rainbow, are federally registered trademarks of Masimo Corporation.

M-LNCS and rainbow ReSposable are trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
PHILIPS is a trademark of Koninklijke Philips Electronics, N.V.
EU Authorized Representative for Masimo Corporation:
MDSS GmbH