


Quick Reference Guide

iOS Masimo Doctella™ App



> Download and install the Masimo Doctella App from the App Store.



> Launch the app and follow the step-by-step instructions to create an account, engage in a shared CareProgram™, and begin collecting and sharing data.

For further assistance or questions visit:
www.Masimo.com/doctella





Not intended to be used in the prevention, diagnosis, or treatment of any disease.

Log In, Sign Up, Notifications

> Ensure notifications are allowed for proper app function.


> Sign Up

- If you received a CareProgram from a health provider, you must sign in with the email or phone number used to send your invitation.
- Sign up for an account using email, phone, or with Google single-sign on as your username.



ALLOW NOTIFICATIONS


Health Permissions



> **Select Next to review the Apple Health Kit menu and allow Doctella to gather data as follows:**

- 1 Select Turn All Categories On to ensure proper app function.
- 2 Select Allow for Health App permission.

Activity Center



> **Menu**

- View the App Menu.

> **Dashboard**

- View the Dashboard screen, which shows the user's health data.
- View health data shared by other users.


> **Sharing**

- Displays the healthcare providers, caregivers, and family members with whom the user is currently sharing data.
- Share data with another healthcare provider, caregiver, or family member.

> **CareProgram**

- Allows user to view or delete a CareProgram.

Menu



> Feedback

- Allows the user to send an email to Doctella support.

> Rate on App Store

- Allows the user to rate the app on the app store.

> Log Out

- Logs the user out of the app.

> Activity Center

- Displays the Activity Center screen.

> Dashboard

- Displays the Dashboard screen.

> Profile

- Displays the user profile information.
- Allows the user to modify the profile information.

> Change Password

- Allows the user to change the password for the app.

> Change Security Questions

- Allows the user to change the security questions in case of lost password.

> Notification Settings

- Allows the user to change the app notification and reminder settings.


> Settings

- Displays the user identification (ID).
- Allows the user to enable or disable Touch ID or Face ID to log into the app.
- Displays app information.

> Restore Purchases

- Users are allowed to view shared data from one user at no charge. Click this to reset the one free allocation, if previously used.

Dashboard



> Invite

- Allows the user to invite others to share information through the app.

> Sharing

- Displays the healthcare providers, caregivers, and family members with whom the user is currently sharing data.


> Edit

- Allows the user to rearrange data on the Dashboard.

> View data for the last 30 days or a weekly summary.


- Select a parameter card to view trends, date ranges, and additional data.

CareProgram




- > **Care Card**
 - View the Care Card data.
- > **Home**
 - Return to the Activity Center.
- > **Learn**
 - View additional information provided by a caregiver.
- > **Care Team**
 - View the list of shared healthcare providers, caregivers, and family members.

Sharing



- > **Displays the list of shared healthcare providers, caregivers, and family members.**
- > **Add**
 - Share with additional healthcare providers.
 - Share with additional caregivers or family members.
- > **Swipe to edit privileges.**



ADD